



Clinical Corner

Column Selection for Toxicological Analyses

Testing for the presence of drugs in biological fluids is a technique that is commonly used in a variety of applications. These applications are broken down into two main categories. The first includes drug overdose screening for clinical and forensic purposes. The second covers drug abuse screening for drug addiction clinics, employee screening and athletic testing. In all cases an analyst screens for as many substances as possible within a limited time frame. Capillary gas chromatography offers the greatest flexibility among the various methods for drug screening. A large number of compounds can be screened in a single chromatographic analysis within a rela-

Capillary gas chromatography offers the greatest flexibility among the various methods for drug screening.

tively short time period. Confirmation of a positive result is usually achieved by re-analyzing the specimen by a second, more specific technique, usually GC-Mass Spectrometry. GC-MS instrument time can be a valuable commodity in laboratories with limited budgets. Therefore, confirmational analyses of suspected positive results are performed only on those samples that have a clear indication of the presence of a positive result. Reliability of identification of an initial positive result can be increased without adding additional testing time through dual column analysis. In MS mode, one injection can be split between two columns that are attached to separate detectors. By using columns with dissimilar stationary phases, retention time shifts and elution order changes can be demonstrated for most analytes.

Several factors should be considered in the selection of appropriate stationary phases for screening and confirming positive results. Overall resolution and peak shape, analysis time, inertness, and thermal stability are all factors that are influenced by the composition of the stationary phase. Comparative values obtained for retention time and elution order also depend on phase composition. Previous studies employing dual column systems have compared methyl silicone and phenylmethyl silicone stationary phases. In this article we have compared these phases along with another intermediate polarity phase, the trifluoropropyl methyl silicone phase (Rtx-200).

Table I (page 8) shows the retention times for a group of drugs on a methyl silicone, several phenylmethyl silicones and the trifluoropropyl stationary phase. The polarity of the stationary phase can be changed by adding increasing percentages of phenyl groups to the silicone polymer. The resulting increase in polarity yields longer retention times for many analytes. This shifting of retention times when columns are run parallel can be used to tentatively identify an unknown substance. Substantial increases in polarity can also cause changes in the elution order for a given set of substances. However, increasing the polarity of columns to change separations can bring diminished returns in extended analysis times and poor peak shape.

Retention time shifts and elution order changes are also accomplished by incorporating different functional groups into the stationary phase. The Rtx-200 exhibits unique selectivity for specific drugs based upon its affinity for compounds that are nucleophilic in nature. Figures 1-3 (page 9) show the analysis of a select group of drugs on the Rtx-5, Rtx-35 and the Rtx-200 columns. Large changes in elution order can be attributed to the presence of specific functional groups in the analyte. Peaks 2 and 4, cotinine and caffeine, are retained significantly longer than the other early eluting compounds. This is due to the presence of a carbonyl group on the ring of both compounds. Morphine and diazepam, peaks 19 and 20, which had been difficult to separate on the phenyl columns are well resolved on the Rtx-200 column. Diazepam and the other benzodiazepines exhibit a strong affinity for the trifluoropropyl stationary phase because of a carbonyl group and an azo group in the molecule. The triazolo benzodiazepines alprazolam and triazolam, peaks 28 and 29, are particularly affected because of the large number of azo groups in each compound. Flunitrazepam, peak 24, also shows a large elution order change relative to the other benzodiazepines due to the nitro group in the molecule. Compounds with electron donating groups (carbonyl, azo and nitro) will be preferentially retained when compared to compounds with similar base structure but do not contain these groups.

The Rtx-200 column can be used in combination with any of the methyl silicone or phenylmethyl silicone stationary phases for a dual column screening and confirmation system for drug analysis. All of the columns have comparable temperature ranges and exhibit excellent inertness. The unique selectivity can be used to improve specific separations of hard to resolve compounds or to offer an alternative phase for confirmation of results by altering retention times and elution order.

Table I - Retention times for common drugs on various polarity stationary phases.

| Compound | Retention Times | | | | | Compound | Retention Times | | | | |
|------------------|-----------------|------------|-------------|------------|------------|-------------------|-----------------|------------|------------|------------|------------|
| | Rtx-1 | Rtx-5 | Rtx-35 | Rtx-50 | Rtx-200 | | Rtx-1 | Rtx-5 | Rtx-35 | Rtx-50 | Rtx-200 |
| | Cat.#10123 | Cat.#10223 | Cat.#10042: | Cat.#10523 | Cat.#15023 | | Cat.#10123 | Cat.#10223 | Cat.#10423 | Cat.#10523 | Cat.#15023 |
| ethosuximide | 6.06 | 7.49 | 10.02 | 11.92 | 9.96 | primidone | 31.76 | 35.2 | 42.83 | 45.88 | 35.15 |
| barbital | 12.96 | 14.87 | 18.27 | 20.31 | 16.29 | promethazine | 32.46 | 35.43 | 39.62 | 41.12 | 30.28 |
| benzocaine | 14.32 | 16.77 | 21.46 | 23.55 | 17.08 | bupivacaine | 32.7 | 35.57 | 38.57 | 39.86 | 34.92 |
| aprobarbital | 16.4 | 18.45 | 21.74 | 23.75 | 18.86 | pentazocine | 32.91 | 35.57 | 38.74 | 40.13 | 30 |
| butabarbital | 17.55 | 19.67 | 22.92 | 24.81 | 20.35 | promazine | 33.46 | 36.6 | 40.9 | 42.46 | 31.64 |
| butalbital | 17.79 | 19.84 | 22.73 | 24.58 | 19.98 | scopolamine | 33.4 | 36.6 | 41.62 | 43.48 | 35.49 |
| cotinine | 17.53 | 20.55 | 26.35 | 28.62 | 25.02 | carbamazepine | 33.23 | 36.64 | 43.76 | 46.4 | 36.79 |
| amobarbital | 19.13 | 21.2 | 23.84 | 25.53 | 21.6 | maprotylme | 33.89 | 36.88 | 40.91 | 42.03 | 32.12 |
| methyl phenidate | 19.56 | 21.76 | 25.4 | 26.72 | 18.3 | diphenylhydantoin | 34.05 | 37.27 | 43.76 | 46.5 | 36.39 |
| pentobarbital | 19.7 | 21.86 | 24.84 | 26.59 | 22.12 | codeine | 34.68 | 37.96 | 43.32 | 45.4 | 34.18 |
| mepredine | 19.58 | 21.95 | 24.71 | 25.86 | 18.49 | clomipramine | 35.38 | 38.41 | 41.83 | 43.04 | 33 |
| secobarbital | 21.08 | 23.25 | 26.09 | 27.86 | 22.9 | lorazepam | 35.26 | 38.65 | 44.07 | 46.3 | 37.34 |
| caffeine | 21.27 | 24.33 | 30.18 | 33.01 | 25.05 | dipam | 35.87 | 39.26 | 44.8 | 47 | 38.63 |
| ketamine | 22.12 | 24.93 | 29.4 | 31.09 | 23.71 | morphine | 36.17 | 39.26 | 44.75 | 47.05 | 36.09 |
| diphenhydramine | 22.76 | 25.31 | 28.4 | 29.41 | 20.04 | hydrocodone | 35.95 | 39.47 | 45.48 | 47.7 | 37.89 |
| lidocaine | 22.97 | 25.58 | 28.68 | 30.06 | 27.71 | hydromorphone | 36.51 | 39.9 | 45.93 | 48.29 | 38.95 |
| phencyclidine | 23.38 | 25.8 | 28.08 | 28.82 | 19.05 | chlorpromazine | 37.1 | 40.31 | 44.29 | 45.86 | 34.99 |
| mephobarbital | 23.66 | 26.29 | 30.7 | 32.66 | 25.58 | chlorprothixene | 37.1 | 40.31 | 44.29 | 45.86 | 34.4 |
| doxylamine | 24.01 | 26.64 | 30.05 | 31.2 | 20.76 | nordiazepam | 37.27 | 40.81 | 47.05 | 49.53 | 40.13 |
| phenyltoloxamine | 24.7 | 27.34 | 30.47 | 31.61 | 22.32 | oxycodone | 37.55 | 41.13 | 47.2 | — | 40.11 |
| phenobarbital | 25.36 | 28 | 32.97 | 35.43 | 27.64 | clobazam | 38.29 | 41.84 | 47.97 | 50.3 | 44.24 |
| tripelennamine | 25.67 | 28.39 | 31.98 | 33.21 | 23.07 | nalorphine | 38.68 | 41.96 | 47.2 | — | 37.99 |
| methapyriline | 25.8 | 28.56 | 32.4 | 33.73 | 23.44 | temazepam | 39.02 | 42.6 | 48.51 | 50.9 | 42.74 |
| chlorpheniramine | 26.27 | 28.96 | 32.09 | 33.3 | 23.9 | flutitrazepam | 39.36 | 43.02 | 48.88 | 51.26 | 45.33 |
| procaine | 26.77 | 29.61 | 33.84 | 35.65 | 27.91 | bromazepam | 39.48 | 43.37 | 50.77 | 53.56 | 41.68 |
| brompheniramine | 28.58 | 31.4 | 34.86 | 36.18 | 26.03 | pt=epam | 40.11 | 43.56 | 48.62 | 50.62 | 41.68 |
| dextromethorphan | 29.34 | 32.08 | 35.52 | 36.84 | 26.37 | trifluoperazine | 40.51 | 43.73 | 46.68 | 47.79 | 38.86 |
| methadone | 29.59 | 32.32 | 35.18 | 36.2 | 26.95 | dibucaine | 41.03 | 44.28 | 47.43 | 48.9 | 41.85 |
| propoxyphene | 30.59 | 33.31 | 35.93 | 36.8 | 28.09 | acetopromazine | 41.18 | 44.64 | 49.32 | 51.23 | 41.32 |
| amitriptyline | 30.71 | 33.55 | 36.93 | 38.06 | 27.6 | flurazepam | 42.44 | 45.81 | 49.96 | 51.61 | 43.45 |
| atropine | 30.86 | 33.74 | 37.83 | 39.41 | 30.61 | papaverine | 43.43 | 47.06 | 53.22 | 55.76 | 43.51 |
| nortriptyline | 31.09 | 33.99 | 38.01 | 39.11 | 28.5 | clonazepam | 43.6 | 47.61 | 54.68 | 57.58 | 49.6 |
| trimipramine | 31.27 | 34.1 | 37.21 | 38.18 | 28.3 | haloperidol | 45.31 | 49.07 | 53.72 | 55.58 | 46.61 |
| imipramine | 31.27 | 34.25 | 37.79 | 38.93 | 28.66 | alprazolam | 45.46 | 49.34 | 56.24 | 59.6 | 52.67 |
| tetracaine | 31.43 | 34.33 | 37.55 | 39.07 | 32.39 | prochlorperazine | 45.79 | 49.41 | 53.98 | 55.72 | 43.21 |
| doxepine | 31.46 | 34.8 | 38.17 | 39.5 | 29.46 | tiazolixl1 | 47.12 | 51.04 | 58.39 | 62.37 | 54.15 |
| pyrilamine | 31.59 | 34.56 | 38.54 | 40.04 | 29.37 | thioridazine | 48.45 | 52.16 | 57.55 | 60.09 | 45.81 |
| medazepam | 31.54 | 34.71 | 39.71 | 41.6 | 29.95 | strychnine | 48.84 | 52.95 | 61.53 | — | 49.69 |
| desipramine | 31.83 | 34.82 | 38.97 | 40.13 | 29.79 | verapamil | 49.17 | 52.89 | 57.49 | 59.8 | 49.95 |
| procainamide | 31.87 | 34.99 | 40.1 | 42.32 | 35.93 | | | | | | |

Run Conditions and Peak List for Figures 1 - 3.

30m. 0.25mm ID, 0.25um
 1.0ul split injection of 100ug/ml (2ng on-column)
 Oven temp.: 100°C to 320°C @ 4C/min. (hold 10 min.)
 Inj. temp.: 225°C
 Det. temp.: 300°C
 Carrier gas: helium
 Linear velocity: 30cm/sec. @ 100°C
 TSD sensitivity: 4 x 10-10 AFS
 Split ratio: 50:1

| COMPOUNDS | | | |
|-----------|------------------|----|----------------|
| 1 | benzocaine | 12 | amitriptyline |
| 2 | cotinine | 13 | trimipramine |
| 3 | mepredine | 14 | imipramine |
| 4 | caffeine | 15 | medazepam |
| 5 | lidocaine | 16 | pentazocine |
| 6 | phencyclidine | 17 | promazine |
| 7 | doxylamine | 18 | codeine |
| 8 | phenyltoloxamine | 19 | morphine |
| 9 | chlorphenamine | 20 | diazepam |
| 10 | daxtromethorphan | 21 | chlorpromazine |
| 11 | methadone | 22 | clorprothixene |
| 23 | clobazam | | |
| 24 | flunitrazepam | | |
| 25 | prazepan | | |
| 26 | flurazepm | | |
| 27 | haloperidol | | |
| 28 | alprazolam | | |
| 29 | triazolam | | |
| 30 | thioridazine | | |
| 31 | verapamil | | |
| 32 | strychnine | | |