



DEGASi – Stand Alone Degassers → DEGASi PLUS GPC

## DEGASi® PLUS GPC

## READ ABOUT OUR PRODUCTS

[IDEX Catalog](#)

[Column Hardware](#)

[Detectors](#)

[Degassing & Debubbling – No troubles with bubbles](#)

### **Degassing Applications**

[Degassing References](#)

[Degasi Finder](#)

### **DEGASi – Stand Alone Degassers**

[DEGASi PLUS Degassers](#)

[DEGASi PLUS CLASSIC](#)

[DEGASi PLUS GPC](#)

[DEGASi PLUS MICRO](#)

[DEGASi PLUS SEMI-PREP](#)

[DEGASi COMPACT Degassers](#)

[DEGASi High Flow 1000](#)

[DEGASi PREP Degassers](#)

[DEGASi PREP+ Degassers](#)

[OLD DEGASSERS](#)

### **Debubbler**

### **Degassing OEM Solutions**

### **Degassing Parts**



## THE RIGHT CHOICE

DEGASi® PLUS GPC is the right choice of degasser if you are working with 100 % organic solvents in your fluidic line. Example of application areas where this degasser is successfully used are GPC (Gel Permeation Chromatography) and Normal Phase Chromatography.

 [Visit our Shop](#)

***“WHEN YOU ARE USING 100 % ORGANIC SOLVENTS IN YOUR FLUIDIC LINE...”***

## THE DIFFERENCE BETWEEN THE NEW AND OLD DEGASi® GPC

- Much smaller footprint, height and weight
- New vacuum pump using a smart board
- Allmost silent
- Available with 1-6 channels



[Fittings](#)

[Tubing](#)

[HPLC Columns](#)

[HPLC Filters and Frits](#)

[Sensors](#)

[Pumps](#)

[Smart Saver – USB Solvent](#)

[Recycler](#)

[KOT – Tubular Coils](#)

[Technical Resources](#)

[Valves](#)

[Volumetric Precision Dispenser](#)

- Parallel connected vacuum chambers
- Easier accessibility and visibility
- Easier troubleshooting and replacement of parts
- Short leadtime

## WHAT IS UNIQUE

DEGASi® GPC uses the same 480 µl Systec AF™ degassing membrane as used in [DEGASi® Classic](#). The main difference compared with the [DEGASi® Classic](#) is that we use a stented version of the vacuum chamber in the DEGASi® GPC. The stent is a short piece of a stainless steel tube placed inside the ends of the degassing membrane, in order to make the internal ferrule in the bulkhead to get a better grip. This solution makes an even more secure internal connection when working with 100 % organic solvents in the fluidic line.

## FEATURES

- Ultra-high degassing efficiency
- Fits every system
- Long life expectancy 5+ years continuous (24/7) running capacity
- Easy to prime
- Almost silent
- Closed loop vacuum control means constant vacuum (variable RPM)
- ZHCR® patented control eliminates baseline fluctuations
- Excellent chemical compatibility flow path
- Continuous vacuum system monitoring to ensure optimum operational conditions are maintained
- Advanced error and leak checking functions
- CE & RoHS compliant



[Download our Brochure](#)