

CHAMOMILE – Convenzional GC vs FAST-GC


Courtesy of Prof. C. Bicchi, C. Brunelli – Università di Torino, Dipartimento Scienza e Tecnologia del Farmaco – Via P.Giuria, 9 – Torino

Column

Phase	MEGA-1701
I.D.	0.25 mm
Film Thickness	0.3 µm
Length	25 m

Chromatographic Conditions

Inlet	Split	230°C	Oven	T start	50°C (0.1 min)
Injected Volume	1.0 µL			Rate	3°C/min
Sample Dilution	1:200	in Cyclohexane		T end	250°C (5 min)
Carrier Gas	Hydrogen	1.5 mL/min	Detector	FID	250°C



What GC ???

CHAMOMILE – Convenzional vs FAST-GC


Courtesy of Prof. C. Bicchi, C. Brunelli – Università di Torino, Dipartimento Scienza e Tecnologia del Farmaco – Via P.Giuria, 9 – Torino

Column

Phase	MEGA-1701 FAST
I.D.	0.1 mm
Film Thickness	0.1 µm
Length	5 m

Condizioni

Inlet	Split	230°C	Oven	T start	50°C (0.1 min)
Injected Volume	1.0 µL			Rate	50°C/min
Sample Dilution	1:200	in Cyclohexane		T end	250°C (5 min)
Carrier Gas	Hydrogen	0.5 mL/min	Detector	FID	250°C



CHAMOMILE

Courtesy of Prof. C. Bicchi, C. Brunelli – Università di Torino, Dipartimento Scienza e Tecnologia del Farmaco – Via P.Giuria, 9 – Torino

Column

Phase	MEGA-WAX FAST
I.D.	0.1 mm
Film Thickness	0.1 µm
Length	5 m

Chromatographic Condition

Inlet	Split	230°C	Oven	T start	50°C (0.1 min)
Injected Volume	1.0 µL			Rate	3°C/min
Sample Dilution	1:200	in Cyclohexane		T end	250°C
Carrier Gas	Hydrogen	0.5 mL/min	Detector	FID	250°C

