

« [Shortcut: QuEChERS EN and Original Methods Protocol](#)
[Early Eluting Terpenes – GC – Medical Cannabis](#) »

Faster GC Analysis of Medical Cannabis Terpenes with Same 624Sil MS Selectivity

March 25th, 2014 by Jack Cochran

The chromatograms below show what happens when you translate a GC method (previously used for medical cannabis terpenes [here](#) and [here](#)) from a [30m x 0.25mm x 1.40µm Rxi-624Sil MS](#) GC column to a 30m x 0.25mm x 1.00µm Rxi-1301Sil MS column. Both of these columns have arylene-modified cyanopropylphenyl dimethyl polysiloxane-type stationary phases. As should be expected, the separation is approximately the same, but faster, due to the thinner film on the 1301Sil MS.

Another thing that the thinner film gives users is the ability to elute less volatile compounds without long isothermal hold times. The example here is for [Phytol](#), an acyclic diterpene alcohol that may be a therapeutically-active compound in cannabis. Phytol wasn't in the terpene mix I analyzed on the 624Sil MS, and in fact, I added several “new” terpenes to my qualitative medical cannabis terpenes standard for the 1301Sil MS work, including Sabinene, alpha-Phellandrene, Ocimene, p-Cymene, alpha-Humulene, alpha-Bisabolol, and Phytol. I've labelled some of those additional terpenes with small letters in the first 1301Sil MS chromatogram and the analyte list, and I show the later eluting additions by name in the second 1301Sil MS chromatogram.

I am continuing this work along several lines, including running on even thinner-film cyanopropylphenyl dimethyl-type columns. Why do that? Because we can pick up efficiency, i.e. more separation power through narrower peaks. Chromatographic efficiency is EXTREMELY important in this type of analysis because there are many, many terpenes in medical cannabis and *high efficiency GC is the only way to go* to achieve the best separations and avoid coelutions that will lead to inaccuracies in terpene quantification for medical marijuana.

I'd like to acknowledge Andrew Goldsmith of [SRI Instruments](#) for providing some of the terpene standards, and Don Rhoads of Restek for making the beta-version Rxi-1301Sil MS GC column.

Please, check out our [Medical Marijuana web page](#).

Analyte	Name
1	alpha-Pinene
2	Camphene
3	Myrcene
4	beta-Pinene
5	delta-3-Carene
6	alpha-Terpinene
7	Limonene
8	1,8-Cineole
9	gamma-Terpinene
10	Terpinolene
11	Linalool
12	Fenchone
13	1-Isopulegol
14	dl-Menthol
15	1-Borneol
16	alpha-Terpineol
17	Dihydrocarveol
18	Citronellol
19	Geraniol
20	2-Piperidone
21	Citral 1
22	Pulegone
23	Citral 2
24	Citral 3
25	Citral 4
26	beta-Caryophyllene
27	Nerolidol 1
28	Nerolidol 2
29	Caryophyllene oxide

a Sabinene

b alpha-Phellandrene


c Ocimene 1


d p-Cymene


e Ocimene 2

**Terpene
Mix New**

f alpha-Humulene







Carrier Gas		Original	Translation
		Helium	Helium
Column			
Length		30.00	30.00 m
Inner Diameter		0.25	0.25 mm
Film Thickness		1.40	1.00 µm
Phase Ratio		44.64	62.50
Control Parameters			
Outlet Flow	→	1.40	→ 1.40 mL/min
Average Velocity		33.32	33.23 cm/sec
Holdup Time		1.50	1.50 min
Inlet Pressure	PSI	16.96	16.80 PSI
Outlet Pressure (abs)		14.70	14.70 PSI
		Atm Vacuum	Atm Vacuum
Oven Program			
<input type="radio"/> Isothermal <input checked="" type="radio"/> Ramps			
	Ramp Rate (°C/min)	Temp (°C)	Hold Time (min)
Number of Ramps (1-4)		60	0.1
1	12.5	300	1.7
	15.9	300	1.35
Control Method			
		Constant Flow	
Results			
Solve for <input type="radio"/> Efficiency <input type="radio"/> Speed <input checked="" type="radio"/> Translate <input type="radio"/> Custom			
Run Time		21.00	16.89 min
Speed			1.24 x

This entry was posted on Tuesday, March 25th, 2014 at 2:38 pm and is filed under [Cannabis](#), [Faster Analyses](#), [Method Translator and Flow Calculator](#), [New GC Columns](#), [Optimizing Applications](#). You can follow any responses to this entry through the [RSS 2.0](#) feed. You can [leave a response](#), or [trackback](#) from your own site.

One Response to “Faster GC Analysis of Medical Cannabis Terpenes with Same 624Sil MS Selectivity”

1. [Restek's EZGC Online Suite that includes the Method Translator and Flow Calculator, and a Chromatogram Modeler, wins a TASIA Award – Chris Nelson, One of the Suite Builders](#) « [ChromaBLOGraphy: Restek's Chromatography Blog](#) says:
[January 26, 2015 at 7:58 pm](#)

[...] Faster GC Analysis of Medical Cannabis Terpenes with Same 624Sil MS Selectivity [...]

Leave a Reply

Name (required) Mail (will not be published) (required) Website

[Entries \(RSS\)](#) and [Comments \(RSS\)](#).

Check Out Our Industry Pages!

[Environmental](#)

[Clinical & Forensics](#)

[Pharmaceutical](#)

[Petrochemical & Chemical](#)

[Foods, Flavors & Fragrances](#)

[Food Safety](#)

[Medical Cannabis](#)

[ISO Quality Credentials](#)

[Employment Opportunities](#)

[Report a Problem / Help](#)

[Patents & Trademarks](#)

[Privacy Policy](#)

[Terms & Conditions](#)

[Subscribe](#)

[Unsubscribe](#)

[Contact Us](#)

[About Restek](#)

[Media Kit](#)

[Sitemap](#)

[News](#)



Restek Corporation, U.S., 110 Benner Circle, Bellefonte, PA 16823

Copyright © 2014 Restek Corporation. All Rights Reserved.



A Company of Owners



