


Products

Home (../index_en.html) / Products (../product_en.html) / What is CoreShell?

What is CoreShell?

A **core shell particle** consists of a non-porous core in the center and a porous layer outside of the core.

A core-shell structure makes theoretical plate increase 50%, so that a 2.6 μm core shell particle column shows the same theoretical plate as a 1.8 μm totally porous particle column.




Schematic view of a core-shell silica

Why core-shell particles, the better performance?

Van Deemter

$$H = Ad_p + B \frac{D_m}{u} + C \frac{d_p^2}{D_m} u$$

**A term**

Eddy diffusion (d_p is particle diameter)

B term

Longitudinal diffusion (D_m is diffusion coefficient)

C term

Mass transfer

H (height equivalent of one theoretical plate) represents the column length per one theoretical plate. The lower the H , the higher the performance.

Van Deemter plot is a curve using the value derived from the above equation as a function of the theoretical plate height and u (mobile phase linear velocity). Van Deemter Equation contributes A, B and C term individually.

A core-shell particles can be reduced all the values of these three terms.

- ✓ Reason why the A term reduces
- ✓ Reason why the B term reduces
- ✓ Reason why the C term reduces

SunShell series is here using the core-shell particles!

SunShell series ([./SunShell_en.html](#))

[about site \(./about_site_en.html\)](#) | [privacy policy \(./privacy_policy_en.html\)](#) | [old site \(./../old/englishindex.html\)](#)

✉ info@chromanik.co.jp

☎ **TEL:+81-6-6581-0885**

☎ **FAX:+81-6-6581-0890** (any time reception) [FAX sheet \(./../pdf/FAX_sheet_en.pdf\)](#)

