




**Natural Gas**

**Rt®-Alumina BOND/KCl**

(PLOT)



1. Methane
2. Ethane
3. Propane
4. Isobutane
5. n-Butane
6. neo-Pentane
7. Isopentane
8. n-Pentane
9. Hexane Isomers
10. Heptane Isomers
11. Benzene
12. C7+

Column Sample Rt®-Alumina BOND/KCl, 50 m, 0.53 mm ID, 10 µm (cat.# 19760)  
natural gas

**Injection**

Inj. Vol.: 500 µL split  
Liner: Gooseneck Splitless (2mm) (cat.# 20795)  
Inj. Temp.: 200 °C  
Split Vent  
Flow Rate: 50 mL/min.

**Oven**

Oven Temp: 45 °C (hold 1 min.) to 200 °C at 10 °C/min. (hold 8.5 min.)

**Carrier Gas**

H<sub>2</sub>, constant pressure (8.0 psi, 55.2kPa)  
Linear Velocity: 45 cm/sec. @ 45 °C

**Detector**

Make-up  
Gas Type: N<sub>2</sub>  
Data Rate: 20 Hz

**Instrument**


HP5890 GC



**Natural Gas**

**Rt®-Alumina BOND/Na<sub>2</sub>SO<sub>4</sub>**

(PLOT)



1. Methane
2. Ethane
3. Propane
4. Isobutane
5. n-Butane
6. neo-Pentane
7. Isopentane
8. n-Pentane
9. Hexane Isomers
10. Heptane Isomers
11. Benzene
12. C7+

Column Sample Rt®-Alumina BOND/Na<sub>2</sub>SO<sub>4</sub>, 50 m, 0.53 mm ID, 10 µm (cat.# 19756)  
natural gas

**Injection**

Inj. Vol.: 500 µL split  
Liner: Gooseneck Splitless (2mm) (cat.# 20795)  
Inj. Temp.: 200 °C  
Split Vent  
Flow Rate: 50 mL/min.

**Oven**

Oven Temp: 45 °C (hold 1 min.) to 200 °C at 10 °C/min. (hold 8.5 min.)

**Carrier Gas**

H<sub>2</sub>, constant pressure (8.0 psi, 55.2kPa)  
Linear Velocity: 45 cm/sec. @ 45 °C

**Detector**

Make-up  
Gas Type: N<sub>2</sub>  
Data Rate: 20 Hz

**Instrument**

HP5890 GC

