

# Passivation of Components used for Sample Transfer and Holding

Gary A. Barone, David Smith, David  
Shelow

Restek Corporation  
[www.restekcorp.com](http://www.restekcorp.com)



# Components That Benefit from Inert Coatings

- Sample-holding vessels
- Valves or regulators
- Transfer tubing
- Vaporizing or flash chambers
- Reactors
- Tubing connectors and fittings

# Restek Coating Solutions

- History of Restek passive coatings for sample pathways and sample chambers
- Silcosteel<sup>®</sup>, Sulfinert<sup>™</sup> and Siltek<sup>™</sup> Coatings

# Introduction

- Physical Properties of Sulfinert™, Siltek™ and Silcosteel® coatings
- Chemical Properties
- Siltek™-treated Liners for Gas Chromatographs
- Example of coated stainless steel used for sampling and storage of organosulfur compounds
- Corrosion Resistance
- Reduction of carbon buildup using coated stainless steel components

# Chemical Inertness Properties of Silcosteel<sup>®</sup> Coating

- Silicon based
- No reactivity to organosulfurs compounds
- No reactivity or adsorption of polar organic compounds such as alcohols, esters, ethers, etc.,.
- Some flaws such as pin-holes and porosity created during coating process
- Insoluble in acids such as HCl, HNO<sub>3</sub>, and H<sub>2</sub>SO<sub>4</sub>
- Highly soluble in caustic environments, pH>8

# Chemical Properties of Sulfinert™ Coating

- Silicon Based Coating; non-reactive
- Initial processing similar to Silcosteel® process
- Secondary and Tertiary coatings applied to reduce effects of pin-holes and porosity
- Improved performance for transfer & storage of organosulfur compounds
- Improved resistance to caustic environments  
pH=10



# Sulfinert™

Patented passivation treatment  
that provides a durable,  
thermally-stable and highly  
inert surface

# Siltek™ Deactivation

- Inertness for acids, bases, neutrals, pesticides (developed for pesticides analysis)
- Low bleed
- Thermal stability
- Durability (acids, bases, water)
- Regeneration. Easily cleaned with simple solvent sonication.
- Siltek™ for SemiVols+



# Physical Properties of Restek Coatings

- Durable and flexible
- Strong adhesion to substrate
- Non-permeable surface
- Stable to high temperatures

# Processing conditions for Restek Coatings

- Chemical vapor deposition
- Applied to stainless steels, carbon steels, high nickel containing alloys of steel, ceramics and glass
- Clean-room process, any oils or dust reduce coating quality
- Vacuum process. Uses either vacuum chambers or connections to items capable of holding vacuum
- Processing temperature of 400°C

# When to use Coatings

- Transfer and holding of adsorptive or reactive species such as:
  - Sulfurs
  - Alcohols
  - Pharmaceuticals
  - Explosives
  - chemical weapons
  - VOC's
  - OP Pesticides and Chlorinated Pesticides
  - Herbicides
- All Analytical Instrument components used for low-level analysis

# When to use Coatings

- Transfer of corrosive materials such as HCl, nitric acid, and sulfuric acid
- In components used to transfer or contain hydrocarbon streams that are prone to “coking”

# Examples

- Siltek™-treated glass inlet liners; reduce breakdown and improve durability
- Project review detailing the storage and transfer of organosulfur containing streams
- How far can we go?
- Other applications benefiting from Sulfinert™-treated components


# Siltek™-treated Inlet Liners



# Siltek™ Deactivation

- Surface modification, not deactivation “layer”
  - Does not require existing silanol groups (ie. can be applied to glass, metal, ceramic)
- Higher level of inertness for Endrin
- More resistant to acid or base attack
  - Stays inert for longer time
- Easily cleaned
  - Solvent rinsing usually acceptable

# Siltek™ Offers Extended Use Range




Restek

[www.restekcorp.com](http://www.restekcorp.com)



# Excellent Performance: 1.5% Endrin Breakdown




# Project Example: Sulfinert™

- To create passivation processes for stainless steel surfaces which will allow the analysis of low- ppbv sulfur gases
  - Chromatographic sampling system
  - Containment vessels (high pressure vessels and air sampling canisters)

# Organosulfur Compounds

- Certain species adsorb to steel surfaces (e.g., hydrogen sulfide)
- Reactions can occur on a non-coated stainless steel surface (e.g., methyl mercaptan)
- Importance of accurate quantitation (e.g., odorants, beverage grade CO<sub>2</sub>, impurities in ethylene and propylene)


# Block diagram of Analytical System



# Sulfur Compounds Investigated

Compound	Formula	Stock conc. (ppm v)	Standard conc. (ppbv)	Standard conc. as S (ppbv)
hydrogen sulfide	H <sub>2</sub> S	100	60	56
carbonyl sulfide	COS	100	60	30
methyl mercaptan	CH <sub>3</sub> SH	100	60	40
ethyl mercaptan	CH <sub>3</sub> CH <sub>2</sub> SH	100	60	30
dimethyl sulfide	CH <sub>3</sub> SCH <sub>3</sub>	100	60	30
dimethyl disulfide	CH <sub>3</sub> SSCH <sub>3</sub>	100	60	40


# Methyl Mercaptan Stability



# List of Sulfur Compounds

		Conc	Conc	Conc as S
Compound Name	Formula	(ppmv)	(ppbv)	(ppbv)
hydrogen sulfide	H <sub>2</sub> S	105	11.51	10.83
carbonyl sulfide	COS	98	10.74	5.73
methyl mercaptan	CH <sub>3</sub> SH	101	11.07	7.38
ethyl mercaptan	CH <sub>3</sub> CH <sub>2</sub> SH	101	11.07	5.71
dimethylsulfide	CH <sub>3</sub> SCH <sub>3</sub>	99	10.85	6.81
dimethyl disulfide	CH <sub>3</sub> SSCH <sub>3</sub>	100	10.96	7.46


# H<sub>2</sub>S at 11.51ppbv



- reference standard
- Silcosteel treated
- Sulfinert treated



# Methyl Mercaptan at 11.07ppbv




reference standard  
Silcosteel treated  
Sulfinert treated


# Extended Stability Study

- Is the Sulfinert™ surface capable of storing 11ppbv sulfurs longer than 6 days?

# H<sub>2</sub>S at 11ppbv for 14 days



# Methyl Mercaptan 11.07 ppbv for 14 Days




# Can the Sulfinert™ coated vessels go to lower levels?

- Is the Sulfinert™ surface treatment stable for sulfurs at 1.5ppbv?

# 1.5ppbv organosulfur Compounds in Sulfinert™-treated Canisters


## Sulfur Stability @ 1.5ppbv



# Example Applications Requiring Sulfinert™-Treatment for Optimum Performance


- Testing of Beverage Grade CO<sub>2</sub> for organosulfur impurities
- Quality testing of Beers/Wines/Distilled Spirits
- Natural Gas / Refinery Gas

# Beverage Grade CO<sub>2</sub> 20ppbv Sulfur standard added






# Typical Run to Determine Sulfur Content of Beer



# Sulfurs and Hydrocarbon analysis of Natural Gas / Refinery Gas



50ppb each

1. hydrogen sulfide
2. carbonyl sulfide
3. methyl mercaptan
4. ethyl mercaptan
5. dimethyl sulfide
6. dimethyl disulfide

- A. methane
- B. ethane
- C. propylene
- D. propane
- E. isobutane
- F. butane
- G. isopentane
- H. pentane
- I. hexane

# Conclusions on Stability of Organosulfur compounds

- Surface treatments for steel surface allow low-ppbv containment, transfer and subsequent analysis of sulfur gases
  - Silcosteel<sup>®</sup> coating good for ppm levels
  - Improved performance to single ppb levels using Sulfinert<sup>™</sup>
- Future directions include testing Sulfinert<sup>™</sup> in a broader array of applications

# Coating Resistance to Acids

- Good resistance to HCl using Silcosteel<sup>®</sup>-AC coating
- Initial data: 14 fold increase in resistance to 37% HCl of a coated vs. non-coated metal coupon
- Currently undertaking a program to quantify and graph resistance to acids

# Silcosteel<sup>®</sup> Coating as a Coking Inhibitor

- Studies being conducted by The Pennsylvania State University in conjunction with Wright Patterson Air Force Base
- Focus on jet fuels, JP-8
- Initial Data indicates tailored Restek coatings reduce deposits by 12-fold
- Silcosteel<sup>®</sup>-coated metal surfaces inhibit catalytic deposition from thermal decomposition of jet fuel

# Advantages of Inert, Coated Steel components

- Sulfur-containing compounds; polar organic compounds
- Coated storage vessels and any steel component in flowpath
- Improved durability compared to glass linings
- Improved lifetimes and reduced diffusion compared to Teflon linings

# Conclusion

- Consideration in sample flowpath construction
- Importance of shielding exposed metal surfaces
- Usefulness and advantages of using inert, coated components to reduce corrosion and carbon buildup