

Methods of Analysis of Polychlorinated Biphenyl Congeners Using an Application-Specific Capillary GC Column

Frank L. Dorman, Gary B. Stidsen, Chris M.
English and Lydia Nolan

Restek Corporation, Bellefonte, PA

Jack Cochran

LECO Corporation, Las Vegas, NV

209 Possible PCB Congeners

- 7 Indicator Congeners (European):
 - 28,52,101,118,138,153,180
 - Indicate presence/absence of PCB's
- 12 Most Toxic Congeners (WHO):
 - 77,81,126,169
 - 105, 114, 118, 123, 156, 157, 167, 189
 - Have Dioxin-like behavior
- 136 Congeners in Aroclor Products
 - Frame, et. al. (Fresenius J. Anal .Chem. 357, 714-22)
 - Add 126, 81, 169, and 209
- 140 Congeners should satisfy almost everyone

Analytical Techniques

● GC-MS

■ Pros:

- Allows for spectral separation of coeluting congeners of different chlorination level
- Single column
- Common instrumentation in many laboratories

■ Cons:

- Moderate sensitivity

Analytical Techniques

● GC-ECD

■ Pros:

- High Sensitivity
- Simplest instrumentation
- Common instrumentation in many laboratories

■ Cons:

- Dual column system usually necessary
- Must have chromatographic separation

Analytical Techniques

● GCxGC-TOFMS (GCxGC – ECD)

■ Pros:

- Allows for spectral separation of coeluting congeners of different chlorination level
- Sensitivity improvement –vs- GC-MS
- Increase in peak capacity for chromatographic separation

■ Cons:

- Instrumentally more complex
- Commercial instrumentation only recently available

Which GC Column?



- Identification of all 209 available on many phases
- **Rtx-PCB** column designed for this analysis
 - Greatest number of resolved congeners
 - High thermal stability
- Rtx-XLB column good for GC-MS
- Rtx-1, 5, CLPesticides, CLPesticides2, **440**, Dioxin2, 500, 35, 1301, etc...

Dioxin-like PCBs

PCB#	Cl#	Cl Pos.		PCB#	Cl#	Cl Pos.
77	4	34-34		156	6	2345-34
81	4	345-4		157	6	234-345
				167	6	245-345
105	5	234-34		169	6	345-345
114	5	2345-4				
118	5	245-34		189	7	2345-345
123	5	345-24				
126	5	345-34				


World Health Organization, non- and mono-ortho substituted

Toxicity of PCBs




- Planar shape
- Fit in a 3x10Å rectangle
- Chlorines on the corners
- Binds with Ah receptor, like dioxins
- Toxicity Equivalency Factors
 - 2378-TCDD = 1.0
 - PCB 77 = 0.0001
 - PCB 126 = 0.1
 - PCB 169 = 0.01

Toxic PCB Congeners (Carbon Cleanup)



GCxGC Schematic



Two independent separation mechanisms

GCxGC Setup

- Primary column (1st dimension)
 - Longer, wider bore, thicker film
 - Non-polar
- Modulator
 - Thermal in nature
 - Concentrates effluent from primary column
 - “Injects” this effluent onto secondary column
- Secondary column (2nd dimension)
 - Very short, narrow bore, thinner film
 - Polar or selective

GCxGC Column

- Integral column – No press fit
 - 50m x 0.18mm x 0.18 μ m Rtx-1
 - 5m x 0.18mm x 0.10 μ m Rtx-PCB
 - Installed so that all 5m of Rtx-PCB was past modulator and in a secondary oven
- Constant flow He at 1.5 mL/min
- Splitless injection
 - 1 μ L
 - 250°C
 - Purge time 60 sec


GCxGC Conditions

- Primary oven
 - 120°C (1 min), 5°/min to 330°
- Modulator (quad jet)*
 - Temperature offset: 40°C
 - Modulation time: 2 sec
- Secondary oven
 - 140°C (1 min), 5°/min to 350°

Run time = 43 min

*Built by LECO under license from Zoex Corporation


Modulator and Secondary Oven



TOFMS Conditions LECO Pegasus 4D


- Source temperature: 225°C
- Electron ionization: 70 eV
- Stored mass range: 120 to 520 u
- Acquisition rate: 50 spectra/sec

Dioxin-like PCBs




Dioxin-like PCBs vs. Other PCBs

Masses: 292 326 360 394




Dioxin-like PCBs in an Aroclor Mix


Masses: 292 326 360 394




Dioxin-like PCBs in an Aroclor Mix



Non-ortho PCBs in Siscowet



Mono-ortho PCBs in Siscowet



Non-ortho PCBs in Siscowet

PCB	HR	TOF	HR	TOF	HR	TOF
77	170	170	140	150	640	630
81	22	24	30	31	83	84
126	610	570	1200	1100	1700	1400
169	250	270	550	490	470	440

- Results in pg/g for three different Siscowet from Lake Superior
- Based on extracted amount and final extract volume:
 - 20 pg/g = ~2.5 pg/μL
 - 200 pg/g = ~25 pg/μL

Conclusions

- GCxGC-TOFMS is a viable way to determine dioxin-like PCBs in fish
 - Sensitive: modulation enhances detection
 - Selective: 2nd dimension separations
- Full mass spectra provide powerful confirmations
- Peak find and spectral deconvolution enable location, identification, and quantification (?) of other environmentally significant compounds

**GC-TOFMS Analysis
of Indicator Polychlorinated
Biphenyls**

Experiment

- Determine retention times for all **209 PCBs** on a new low-bleed capillary column
 - Analyze individual congener standards
 - Hexachlorobenzene for relative retention times
- Analyze mix of Aroclors 1221:1242:1254:1262
 - 1:2:2:2 ratio
 - Show chromatogram with congeners labeled
- GC-MS results obtained using a time-of-flight mass spectrometer to allow for possible spectral deconvolution

GC-TOFMS Conditions

● Gas Chromatography


- One microliter splitless 250°C, 60 sec valve
- 60 m x 0.18 mm x 0.18 μm **Restek Rtx-PCB** column
 - Constant flow helium, 1.5 mL/minute
 - 70°C (1 min), 50°/min to 120°, 3°/min to **360°** (1 min)


● Mass Spectrometry **LECO Pegasus III**


- Source temperature: 225°C
- Electron ionization: 70 eV
- Stored mass range: 120 to 520 u
- Acquisition rate: 5 spectra/sec


Aroclor Mix Chromatogram


- Overlaid ion chromatograms are plotted
 - For example, hexachlorobiphenyl is 360
 - Peaks are labeled with congener numbers
 - Text size and format is only for space and not meant to indicate anything about PCBs, except:
 - **European indicator congeners** are in red
- Chromatograms are not to one scale
 - Each panel is adjusted so most of the PCB peaks can be seen




















138







Important Congener Separations

- European indicator PCBs
- May coelute on 5% phenyl methylsilicone and other column phases

28

Time (seconds)

2285

2290

2295

2300


2305


2310

2315

2320

256





138

Time (seconds)

3185

3187.5

3190

3192.5


3195

3197.5

3200

3202.5

360



Conclusions


- Greater than 140 Aroclor PCBs can be determined with the new column and mass spectrometry
- The seven European indicator PCBs can be determined individually
- Spectral deconvolution available with TOFMS allows qualitative identification for coeluting PCBs

GC-ECD Analysis

- Columns should be “complementary”
 - As many congeners resolved as possible
 - Different coelutions (nothing will resolve all 140 environmentally-significant congeners)
 - Similar thermal properties
 - Operation under same temperature and pressure programs

Rtx-PCB


(40m x 0.18mm, 0.18um film)



300ppb Aroclor 1242, 1254, 1262
100°C(1 min) 20°C/min 200°C(0 min)
2°C/min 320°C(1 min)

Rtx-440

(40m x 0.18mm, 0.18um film)



300ppb Aroclor 1242, 1254, 1262
100°C(1 min) 20°C/min 200°C(0 min)
2°C/min 320°C(1 min)

Quantitation of Congeners

Rule#1: If compound elutes as single on both columns report lower value

	<u>Rtx-PCB</u>	<u>Rtx-440</u>
Peak	PCB 28	PCB 28

Quantitation of Congeners

Rule#2: If compound elutes as single on one column and multiple congeners on the second report the single column result if confirms on second column

	<u>Rtx-PCB</u>	<u>Rtx-440</u>
Peak	PCB 53 + 31	PCB 31

Quantitation of Congeners

Rule #3: If there are complicating coelutions on both columns calculate the results using normalized responses (ex. PCB 136)

	<u>Rtx-PCB</u>	<u>Rtx-440</u>
Peak	PCB 136 + 117	PCB 115 + 136
	PCB 115	

Quantitation of Congeners

Rule #3 continued:


PCB 136 Concentration =

$$\left([\text{PCB 115+136}] \times \frac{(\text{RF115} + \text{RF136})}{2 \times \text{RF115}} - [\text{PCB 115}] \right) \frac{\text{RF115}}{\text{RF136}}$$

RF = response factor of congener on ECD


Rtx-PCB

(40m x 0.18mm, 0.18um film)



Rtx-440

(40m x 0.18mm, 0.18um film)



PCB Congener Separation

141

Rtx-440		Rtx-PCB		50	72	19	
PCB#	s or m	PCB#	s or m	Case			
				1	2	3	
1	s	1	s	1			
3	s	3	s	1			
4	m	4	m			1	4+10 reported together
5	m	5	s		1		
6	s	6	s	1			
7	m	7	s		1		
8	m	8	s		1		
9	m	9	s		1		
10	m	10	m			1	4+10 reported together
12	m	12	m			1	ok
13	s	13	m		1		
15	s	15	s	1			
16	m	16	s		1		
17	m	17	s		1		
18	s	18	s	1			
19	s	19	s	1			
20	m	20	m			1	33+20 reported together
22	m	22	s		1		

GC-ECD Summary

- Rtx-PCB and Rtx-440 Combination
 - Quantification of the 141 PCB congeners possible using the 3-case quantification system
 - High thermal stability
 - Good choice for dual column operation

Overall Conclusions

- GC-MS analysis with Rtx-PCB column
 - European indicator unbiased
 - WHO-12 with carbon cleanup
- GCxGC-TOFMS
 - TOFMS has acquisition speed necessary to characterize the narrow peaks
 - WHO-12 resolved from non-toxic without carbon cleanup
- Dual-Column GC-ECD
 - Also successful, with high sensitivity
 - All congeners of significance can be quantified