

ArDSieve™

CHROMATOGRAPHIC COLUMN TO SEPARATE Ar AND O₂

Benefit from the new breakthrough in material science with the ArDSieve™ chromatographic column that separates argon and oxygen molecules.

Separating argon and oxygen in chromatography has always been difficult as both molecules coelute unless operated at cryogenic temperature using expensive setup. This is even more true when trace ppb argon need to be separated from pure oxygen.

This is a story of the past with the ArDSieve™ column.




UNIQUE FEATURES

- Resolve Ar and O₂ separation at 50 °C column temperature. No need for cryogenic temperature
- Improved separation compared to other columns due to proprietary plasma oxidation treatment
- Improved peak symmetry and reduced eddy diffusion with narrow mesh size range (60/65)
- Improved peak symmetry with ASDevices LipLOK™ end column fittings
- Increased durability with HydraGuard™ protection layer
- No need for consumable oxygen trap
- Achieve lower limit of detection compared to GC system using O₂ traps

APPLICATIONS

- ▶ Trace measurement argon in pure oxygen
- ▶ Trace measurement argon in electronics bulk gases
- ▶ Any other analysis that requires Argon and O₂ separation



PROPRIETARY ArDSieve™ MATERIAL

The material used in the ArDSieve™ column is the outcome of an intensive R&D program and decades of experience. The material is a combination of clinoptilite, an ion exchanged chabazite and proprietary treatments.



OXIDATION TREATMENT

Recent advances in material science have allowed ASDevices™ to better oxidise its solid phase material. This is made possible with a proprietary mixture containing O₂ which is introduced into a plasma chamber.


This treatment considerably improves the efficiency and provides a better argon and oxygen separation.

DEHYDRATION AND HYDRAGUARD™ LAYER

In order to obtain a good argon and oxygen separation, it is important to properly dry the column material. With our decades of experience, we have developed an enhanced dehydration process which further improves the column performance. We have also introduced a moisture protection layer which we have called HydraGuard™. The HydraGuard™ layer which is on both sides of the column avoids column contamination when manipulating the column or when a contaminated sample is injected.

ASDevices PACKED COLUMN DESIGN FOR PERFORMANCE

Doing pack column is easy. Doing great pack column is difficult, especially when trying to measure trace impurities. At ASDevices™, we have over 30 years of experience that we share with our customers.



MESH SIZE

Mesh size range is very often overlooked with packed column. Having a large mesh size range causes eddy diffusion which is well known in chromatography to increase peak broadening. The other issue is reproducibility from one column to another. We consequently manufacture all of our columns with a mesh size of 60/65 for optimum performance.




Compatible with the industrial standard double ferrule type

- ▶ Leak detection capability
- ▶ Better leak integrity
- ▶ No dead volume

TYPICAL CHROMATOGRAPHIC SETUP

The schematic below represents a typical chromatographic setup used with the ArDSieve™ column to measure trace argon in electronics bulk gases. Combining the Ar/O₂ separation properties of the ArDSieve™ column with the sensitivity of our Cubed™ detector based on our Epd technology* provides a powerful tool to achieve sub ppb limit of detection and this without the maintenance associated with traditional GC instruments based on O₂ traps.



MEASUREMENT	COLUMN 1	COLUMN 2	DETECTOR
Ar in O ₂ sample	ArDSieve™ 60/65	ArDSieve™ 60/65	Cubed™
Ar in H ₂ , N ₂ sample	Molecular Sieve 5A 60/65	ArDSieve™ 60/65	Cubed™
Ar in he	ArDSieve™ 60/65		Cubed™

All ASD products respect CE/RoHS (CSA optional) standards.

LEGAL DISCLAIMER

THE INFORMATION GIVEN IN THIS DOCUMENT (INCLUDING BUT NOT LIMITED TO CONTENTS OF REFERENCED WEBSITES) IS GIVEN AS A HINT FOR THE IMPLEMENTATION OUR PRODUCTS AND TECHNOLOGIES COMPONENT ONLY AND SHALL NOT BE REGARDED AS ANY DESCRIPTION OR WARRANTY OF A CERTAIN FUNCTIONALITY, CONDITION OR QUALITY OF THE ASDEVICES TECHNOLOGIES COMPONENT. THE RECIPIENT OF THIS DOCUMENT MUST VERIFY ANY FUNCTION DESCRIBED HEREIN IN THE REAL APPLICATION. ASDEVICES HEREBY DISCLAIMS ANY AND ALL WARRANTIES AND LIABILITIES OF ANY KIND (INCLUDING WITHOUT LIMITATION WARRANTIES OF NON-INFRINGEMENT OF INTELLECTUAL PROPERTY RIGHTS OF ANY THIRD PARTY) WITH RESPECT TO ANY AND ALL INFORMATION GIVEN IN THIS APPLICATION NOTE. THIS MATERIAL IS PROPRIETARY TO ANALYTICAL SENSING DEVICES AND CANNOT BE COPIED, REPRODUCED OR DISSEMINATED IN ANY WAY WITHOUT ITS PRIOR WRITTEN APPROVAL. THE TECHNOLOGY DESCRIBED HEREIN MAY BE SUBJECT TO PATENT PROTECTION OR OTHER FORMS OF INTELLECTUAL PROPERTY RIGHTS.

ASDevices is a fully owned Canadian company that is proud to design and manufacture its products in Canada

* = patent pending