



# CALIDUS

**ultraFAST** micro GAS CHROMATOGRAPH

from . . . **FALCON** Analytical

## Modular

### ultraFAST GC

**Detector Modules**  
microFID 100Hz  
microTCD 50Hz

portable  
Lab GC  
Process Control

less than 10Kg  
43x22x28cm

## Unique Features Dual Column System

**300watts max**

100/220VAC


50/60Hz

24VDC Ext

**Columns 2 metre**  
Restek MXT Capillary  
0.28, 0.53  
and new 0.18mmID  
many Liquid Phases  
PLOT  
MolSieve 5A  
Alumina  
Porous Polymer; Hayesep  
microPackedColumns

**Oven Modules**

- 2 direct heated  
600degC per min  
400degC max  
from 5degC  
above ambient




Sample Processing Unit  
Column Switching Module

**Operating System**  
ChromPerfect 100Hz 24bit  
/ Win PC(external)

### ultraFAST GC

. . . **10 to 50 times faster than conventional GC** **from gases**  
VOCs in 20 to 60seconds, semi-VOCs in 60-120seconds to **C60**



**LCD Touch Screen**  
for operating parameters

- temperatures
- pressures
- on/off functions
- analysis status

**Cooling Fans**



TCD

PTCM

FID



**AutoSampler Systems**

Gas  
Liquid  
Headspace

**CHROMalytic** +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)


**12/13**

Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-Mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel : 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ the Modular micro Gas Chromatograph

- Sample Processing Unit
- Plug & Play Programmed Temperature Column Modules
- Plug & Play Detector Modules
  - Flame Ionization
  - Thermal Conductivity



A CALIDUS micro gas chromatograph features fully-independent Programmed Temperature Column Modules (PTCM) that can be operated isothermally and interfaced to a Sample Processing Unit (SPU). Detector Modules (DM) are also fully-independent and can be mixed and matched within the CALIDUS Sample Processing Unit environment.

- **SPU** - standard with a split/splitless injection port (1:1 up to 1:200) suitable for gas or liquid samples via a syringe through the septum injection, optional automated sampling valves for gas or liquids or an optional auto-sampler capable of liquid or heated headspace gas samples. The inlet includes septum purge to prevent bleed components from entering the system.
- **PTCM** - resistively heated steel capillary chromatography column with necessary hardware, software and electronic control for temperature programming from 0.1°C to 10°C per second from 5°C above ambient to 400°C depending on the model and maximum temperature capability of the column material selected.
- **DM** - incorporates micro Flame Ionization Detector (FID) or micro Thermal Conductivity Detector (TCD) with the necessary hardware, software and electronic control to provide detector temperature control and digital output signal.
  - **FID** - fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. Control is provided for the fuel supply pressure and auto-ignition. The data rate is 100 Hz.
  - **TCD** - fully digital, universal detector consisting of a constant temperature filament sensing the change in power required to hold the filament temperature constant due to chromatographic component elution. The data rate is 50 Hz.

CALIDUS is controlled with ChromPerfect chromatography data system and fully integrated with LineUp running on a Windows PC.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518  
Ronceverte, WV 24971  
  
Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: [info@falconfast.net](mailto:info@falconfast.net)

### SIGNIFICANCE AND USE

The CALIDUS micro gas chromatograph provides a simple ultra-fast analysis (10-50 times faster than conventional lab or process GC's) of fixed gases and hydrocarbons up to C<sub>50</sub>. CALIDUS is available for laboratory, at-line, transportable or online use in the hydrocarbon processing industry, environmental labs, pharmaceuticals, food and beverage industry, military, medical industry, and educational markets.

The analyses are used for product specifications testing, product safety, environmental testing and measurements, process control, catalyst protection, educational tools, spot checks of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**



## CALIDUS™ micro Gas Chromatograph

- Model 101
- Model 101 HT
- Model 201
- Model 301
- Model CS



5 CALIDUS Models combine various standard modules to provide general or specific applications and expanded measurement capabilities. Each of these models can be installed as a plug and play module within the Calidus process analyzer enclosure.

**CALIDUS Model 101** - 3 modules, a Sample Processing Unit, a Programmed Temperature Column Module and a Detector Module interfaced with the ChromPerfect chromatography data system and fully integrated with LineUp running on a Windows PC. The user may select either a micro FID or TCD and one of the available different PTCMs to separate and measure fixed gases and hydrocarbons up to C<sub>44</sub>.

**CALIDUS Model 101 HT** - a 101 with a PTCM using High Temperature MXT-1HT Sim Dist column and a micro FID module specifically for determination of boiling range distribution of petroleum products and biodiesel formulations up to C<sub>50</sub> in boiling point. The analyzer is interfaced with the ChromPerfect chromatography data system, SimDis 2000 software, and fully integrated with LineUp running on a Windows PC. An ASTM method for Ultra Fast Micro GC 2887 is currently under development based on the CALIDUS Model 101 HT.

**CALIDUS Model 201** - 4 modules, 2 PTCMs in series with one Sample Processing Unit and one Detector Module interfaced with ChromPerfect chromatography data system and fully integrated with LineUp running on a Windows PC. There are two major advantages for having two PTCMs in series. First is the ability to leverage selectivity of different stationary phases. And second, is virtually doubling the column length for even greater separation power up to C<sub>44</sub>.

**CALIDUS Model 301** - 5 modules, a Sample Processing Unit with a single injector connected to a splitter dividing the sample between two PTCMs in parallel, each with a single micro FID or TCD Detector Module. The Model 301 handles hydrocarbon samples with a wide range of boiling points and a wide range of concentrations (% to ppm) with better separation and faster analysis all without complicated valve schemes and resultant additional hardware.

**CALIDUS Model CS (Column Switching)** - 5 modules, a Sample Processing Unit with a single injector connected to a 6-port diaphragm/plunger column valve, two PTCMs in parallel, and two detector Modules. This model can be plumbed to perform heartcutting from one PTCM with its own FID or TCD Detector Module to a second PTCM with its own FID or TCD Detector Module. Backflushing configurations are available too. This model is used for analysis where a specified discrete hydrocarbon (s) must be separated and measured from a defined stream or sample composition typically within a required time frame with optimum selectivity (up to C<sub>12</sub>).

**CALIDUS**  
micro GAS CHROMATOGRAPH



TCD



PTCM



FID

Faster, Smaller, Smarter, Easier, Greener



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ 101 micro Gas Chromatograph

GC analysis for virtually any fixed gas and hydrocarbons up to C<sub>44</sub> for laboratory, at-line, transportable or online use

- Upstream (E&P)
- Petrochemical/Chemical
- Food & Beverage
- Military
- Refining
- Pharmaceutical
- Medical
- Educational



**CALIDUS Model 101** - comprised of 3 modules

A single Sample Processing Unit with a standard split/splitless injection port (1:1 up to 1:200) suitable for gas and liquid samples via either syringe through the septum injections, optional gas, liquid or headspace auto-sampler, or automated sampling valves. The inlet includes septum purge to prevent bleed components from entering the system.

A single Programmed Temperature Column Module containing the resistively heated steel capillary chromatography column with necessary hardware, software and electronic control to enable temperature programming from 0.1°C to 10°C per second from 5°C above ambient to 350°C depending on the maximum temperature capability of the column material selected.

A single Detector Module incorporating either a micro Flame Ionization Detector (FID) or micro Thermal Conductivity Detector (TCD) with the necessary hardware, software and electronic control to provide detector temperature control, digital output signal and additionally for proper FID fuel supply pressure and auto-ignition.

The micro FID is a fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. The data rate is 100 Hz.

The micro TCD is a fully digital, universal detector. The TCD consists of a constant temperature filament that senses change in power required to hold the filament temperature constant when chromatographic components elute. The power measurement is used to determine the amount of the component eluting from the column. The data rate is 50 Hz.

CALIDUS is controlled with ChromPerfect chromatography data system fully integrated with LineUp running on a Windows PC.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518  
Ronceverte, WV 24971  
  
Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: [info@falconfast.net](mailto:info@falconfast.net)

### SIGNIFICANCE AND USE

The Model 101 is the simplest configuration of the CALIDUS micro GC systems. It provides ultra-fast analysis (10 to 50 times faster than conventional lab or process GC's) of fixed gases and hydrocarbons to C<sub>44</sub>. It is ideal for sample scouting, methods development and transportable uses. Simplicity doesn't negate the powerful capability of this model.

The analyses are used for product specifications testing, product safety, environmental measurements, process control, catalyst protection, educational tools, spot checks of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**



# CALIDUS™ 101 Specifications (global patents pending)

December 1, 2010

## Ambient Environment

Operating Temperature Range: 0°C to 35°C  
Storage Temperature Range: -20°C to 60°C  
Relative Humidity Range: 0 to 100% (non-condensing)

## Power Requirements

Less than 300 watts peak power at startup, practical use < 200 Watts for gas or liquid analyses  
24 VDC supplied from external power supply, 100-240VAC using 50/60Hz AC

## Safety

General purpose, light industrial (lab instrument environment)  
CE Mark and Nationally Recognized Testing Laboratory (NRTL) certified pending (TUV Rheinland)

## Gas Supplies

50 PSIG, 99.995% H<sub>2</sub> or He at up to 250 ml/min, 50 PSIG zero air for FID operation

## Sample Requirements (via split/splitless injector with septum purge)

Air or gaseous samples at 0 to 50 PSIG at ambient temperature  
Membrane, SPME and static and dynamic headspace extracts  
Direct liquid injections neat or dilute organic solvents (DCM, hexane, MEK, toluene, methanol, CS<sub>2</sub> etc.)

## Dimensions

17" wide by 8.5" deep by 11" high, ~ 20 lbs  
Uninterrupted power supply and data acquisition computer external to the base unit

## Controls/Outputs

All functions and parameters via ethernet or RS-232 using ChromPerfect™ software  
Start analysis from keyboard or GC  
Set method from external computer using ChromPerfect software  
50-100Hz digitization (detector dependent) on each column, 24 bit resolution, auto zero on each run  
Trigger in and ready out signals plus an array of others via ChromPerfect

## Front Panel Displays

Temperature and pressure readings, function on/off, other  
Power on/off  
Status of analysis columns (isothermal, programming, cool down, ready, cycles run, other)

## Standard Equipment

One capillary column, 2m long, 100µm to 320µm ID, temperature programmable from 0.1 to 10°C per second from 5°C above ambient to 350°C (maximum temperature software limited to be no greater than the limit for the columns installed, isothermal operation is available). Column modules are 2 meter columns in Mxt-1, Mxt-5, Mxt-1701, Mxt-Wax, Mxt-MoleSieve, Mxt-Alumina, < 320µm & various film thicknesses with others coming soon.  
Flame ionization or thermal conductivity (filament) detection  
Gas and liquid inlet for syringe injection of samples or automated gas and liquid sample valves available



## Performance (application dependent)

Repeatability of ± 1% RSD or better (area) and of ± 0.1% RSD or better (retention times)  
Analysis times for VOCs: can be <20 seconds and for SVOCs: can be <60 seconds  
Dynamic range: depends on detector used and application (FID typically 10<sup>5</sup>)

## Data Processing and Instrument Control

Note: computer system is integral and necessary component of the analysis system and includes the following requirements:

- RS-232 or USB to RS-232 adapter, ethernet
- Windows XP or newer operating environment
- ChromPerfect software for dual column data acquisition via RS-232 serial or ethernet ports



Faster, Smaller, Smarter, Easier, Greener



**HROM**alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ 101 HT micro Gas Chromatograph

Simulated Distillation GC analysis for virtually any hydrocarbon fuel or fuel blend component up to C<sub>50</sub> for laboratory, at-line, transportable or online use.

- **Upstream (E&P)**
- **Petrochemical/Chemical**
- **Educational**
- **Refining**
- **Military**



**CALIDUS Model 101 HT** - comprised of 3 modules fitted with high temperature inlet and column operational components

A single Sample Processing Unit with a standard split/splitless injection port (1:1 up to 1:200) suitable for gas and liquid samples via either syringe through the septum injections, optional gas, liquid or headspace auto-sampler, or automated sampling valves. The inlet includes septum purge to prevent bleed components from entering the system. The maximum operating temperature is 350°C.

A single Programmed Temperature Column Module containing a high temperature resistively heated steel capillary chromatography column with necessary hardware, software and electronic control to enable temperature programming from 0.1°C to 5°C per second from 5°C above ambient to 400°C. The column is Mxt 1-HT for high temperature simulated distillation gas chromatography.

A single Detector Module incorporating a micro Flame Ionization Detector (FID) with the necessary hardware, software and electronic control to provide detector temperature control (350°C maximum), digital output signal and additionally for proper FID fuel supply pressure and auto-ignition.

The micro FID is a fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. The data rate is 100 Hz.

CALIDUS 101-HT is controlled with ChromPerfect chromatography data system fully integrated with LineUp and SimDist 2000 running on a Windows PC. An ASTM method for Ultra Fast Micro GC D-2887 is currently under development based on CALIDUS 101-HT.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518  
Ronceverte, WV 24971  
  
Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: info@falconfast.net

### SIGNIFICANCE AND USE

The Model 101 HT provides a simple ultra-fast (10 to 50 times faster than conventional lab or process GC's) simulated distillation analysis for hydrocarbons to C<sub>50</sub>. This system is ideally configured for liquid fuels and fuel component characterization by boiling range distribution including gasoline range organics up through gas oil and even crude oil.

The analyses are used for exploration & production liquids characterization, fuels specification testing, regulatory evaluations, environmental measurements, process control, transportable spot check of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**



# CALIDUS™ 101 HT Specifications (global patents pending)

December 1, 2010

## Ambient Environment

Operating Temperature Range: 0°C to 35°C  
Storage Temperature Range: -20°C to 60°C  
Relative Humidity Range: 0 to 100% (non-condensing)

## Power Requirements

Less than 300 watts peak power at startup, practical use < 200 Watts for gas or liquid analyses  
24 VDC supplied from external power supply, 100-240VAC using 50/60Hz AC

## Safety

General purpose, light industrial (lab instrument environment)  
CE Mark and Nationally Recognized Testing Laboratory (NRTL) certified pending (TUV Rheinland)

## Gas Supplies

50 PSIG, 99.995% H<sub>2</sub> or He at up to 250 ml/min, 50 PSIG zero air for FID operation

## Sample Requirements (via split/splitless injector with septum purge)

Air or gaseous samples at 0 to 50 PSIG at ambient temperature  
Membrane, SPME and static and dynamic headspace extracts  
Direct liquid injections neat or dilute organic solvents (DCM, hexane, MEK, toluene, methanol, CS<sub>2</sub> etc.)

## Dimensions

17" wide by 8.5" deep by 11" high, ~ 20 lbs  
Uninterrupted power supply and data acquisition computer external to the base unit

## Controls/Outputs

All functions and parameters via ethernet or RS-232 using ChromPerfect™ software  
Start analysis from keyboard or GC  
Set method from external computer using ChromPerfect software  
50-100Hz digitization (detector dependent) on each column, 24 bit resolution, auto zero on each run  
Trigger in and ready out signals plus an array of others via ChromPerfect

## Front Panel Displays

Temperature and pressure readings, function on/off, other  
Power on/off  
Status of analysis columns (isothermal, programming, cool down, ready, cycles run, other)

## Standard Equipment

One capillary column, 2m long, 100µm, Mxt 1-HT, temperature programmable from 0.1 to 5°C per second from 5°C above ambient to 400°C (maximum temperature software limited to be no greater than the limit for the Mxt 1 HT column installed, isothermal operation is available).  
Flame ionization detection, data rate 100 Hz  
Gas and liquid inlet for syringe injection of samples or automated gas and liquid sample valves available




## Performance (application dependent)

Repeatability of ± 1% RSD or better (area) and of ± 0.1% RSD or better (retention times)  
Analysis times for fuel and fuel components: can be <60 seconds and up to 300 seconds  
Dynamic range: depends on detector used and application (FID typically 10<sup>5</sup>)

## Data Processing and Instrument Control

Note: computer system is integral and necessary component of the analysis system and includes the following requirements:

- RS-232 or USB to RS-232 adapter, ethernet
- Windows XP or newer operating environment
- ChromPerfect software for single column data acquisition via RS-232 serial or ethernet ports
- CALIDUS 101 HT comes with LineUp peak alignment and Simdist-2000 simulated distillation software fully integrated.



Faster, Smaller, Smarter, Easier, Greener



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ 201 micro Gas Chromatograph

GC analysis for virtually any fixed gas and hydrocarbons up to C<sub>44</sub> for laboratory, at-line, transportable or online use

- **Upstream (E&P)**
- **Petrochemical/Chemical**
- **Food & Beverage**
- **Military**
- **Refining**
- **Pharmaceutical**
- **Medical**
- **Educational**



**CALIDUS Model 201** - comprised of 4 modules

A single Sample Processing Unit with a standard split/splitless injection port (1:1 up to 1:200) suitable gas and liquid samples via either syringe through the septum injections, optional gas, liquid or headspace auto-sampler, or automated sampling valves. The inlet includes septum purge to prevent bleed components from entering the system.

Two Programmed Temperature Column Modules (PTCM) in series containing the resistively heated steel capillary chromatography column with necessary hardware, software and electronic control to enable temperature programming from 0.1°C to 5°C per second from 5°C above ambient to 350°C depending on the maximum temperature capability of the column material selected. Each column module is independently controlled by the method and can be any of the available column types.

A single Detector Module incorporating either a micro Flame Ionization Detector (FID) or micro Thermal Conductivity Detector (TCD) with the necessary hardware, software and electronic control to provide detector temperature control, digital output signal and additionally for proper FID fuel supply pressure and auto-ignition.

The micro FID is a fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. The data rate is 100 Hz.

The micro TCD is a fully digital, universal detector. The TCD consists of a constant temperature filament that senses change in power required to hold the filament temperature constant when chromatographic components elute. The power measurement is used to determine the amount of the component eluting from the column. The data rate is 50 Hz.

CALIDUS is controlled with ChromPerfect chromatography data system fully integrated with LineUp running on a Windows PC.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518  
Roncerverte, WV 24971

Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: info@falconfast.net

### SIGNIFICANCE AND USE

The Model 201 provides a simple ultra-fast analysis (10 to 50 times faster than conventional lab or process GC's) of fixed gases and hydrocarbons to C<sub>44</sub>. Use of two different PTCMs in series, for example one polar and one non-polar column material enables leveraging the selectivity differences for enhanced separations. Secondly, using two identical PTCMs virtually doubles the column length.

The analyses are used for product specifications testing, product safety, environmental measurements, process control, catalyst protection, educational tools, spot checks of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**





# CALIDUS™ 201 Specifications (global patents pending)

December 1, 2010

## Ambient Environment

Operating Temperature Range: 0°C to 35°C  
Storage Temperature Range: -20°C to 60°C  
Relative Humidity Range: 0 to 100% (non-condensing)

## Power Requirements

Less than 300 watts peak power at startup, practical use < 200 Watts for gas or liquid analyses  
24 VDC supplied from external power supply, 100-240VAC using 50/60Hz AC

## Safety

General purpose, light industrial (lab instrument environment)  
CE Mark and Nationally Recognized Testing Laboratory (NRTL) certified pending (TUV Rheinland)

## Gas Supplies

50 PSIG, 99.995% H<sub>2</sub> or He at up to 250 ml/min, 50 PSIG zero air for FID operation

## Sample Requirements (via split/splitless injector with septum purge)

Air or gaseous samples at 0 to 50 PSIG at ambient temperature  
Membrane, SPME and static and dynamic headspace extracts  
Direct liquid injections neat or dilute organic solvents (DCM, hexane, MEK, toluene, methanol, CS<sub>2</sub> etc.)

## Dimensions

17" wide by 8.5" deep by 11" high, ~ 20 lbs  
Uninterrupted power supply and data acquisition computer external to the base unit

## Controls/Outputs

All functions and parameters via ethernet or RS-232 using ChromPerfect™ software  
Start analysis from keyboard or GC  
Set method from external computer using ChromPerfect software  
50-100Hz digitization (detector dependent) on each column, 24 bit resolution, auto zero on each run  
Trigger in and ready out signals plus an array of others via ChromPerfect

## Front Panel Displays

Temperature and pressure readings, function on/off, other  
Power on/off  
Status of analysis columns (isothermal, programming, cool down, ready, cycles run, other)

## Standard Equipment


Two capillary columns, 2m long, 100µm to 320µm ID, temperature programmable from 0.1 to 5°C per second from 5°C above ambient to 350°C (maximum temperature software limited to be no greater than the limit for the columns installed, isothermal operation is available). Column modules are 2 meter columns in Mxt-1, Mxt-5, Mxt-1701, Mxt-Wax, Mxt-MoleSieve, Mxt-Alumina, < 320µm and various film thicknesses with others coming soon.  
Flame ionization or thermal conductivity (filament) detection  
Gas and liquid inlet for syringe injection of samples or automated gas and liquid sample valves available

## Performance (application dependent)

Repeatability of ± 1% RSD or better (area) and of ± 0.1% RSD or better (retention times)  
Analysis times for VOCs: can be <20 seconds and for SVOCs: can be <60 seconds  
Dynamic range: depends on detector used and application (FID typically 10<sup>5</sup>)

## Data Processing and Instrument Control

Note: computer system is integral and necessary component of the analysis system and includes the following requirements:  
RS-232 or USB to RS-232 adapter, ethernet  
Windows XP or newer operating environment  
ChromPerfect software for dual column data acquisition via RS-232 serial or ethernet ports



Faster, Smaller, Smarter, Easier, Greener



**HRM**alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ 301 micro Gas Chromatograph

GC analysis for virtually any fixed gas and hydrocarbons up to C<sub>44</sub> for laboratory, at-line, transportable or online use

- **Upstream (E&P)**
- **Petrochemical/Chemical**
- **Food & Beverage**
- **Military**
- **Refining**
- **Pharmaceutical**
- **Medical**
- **Educational**



**CALIDUS Model 301** - comprised of 5 modules

A single Sample Processing Unit with a standard split/splitless injection port (1:1 up to 1:200) suitable gas and liquid samples via either syringe through the septum injections, optional gas, liquid or headspace auto-sampler, or automated sampling valves. The inlet includes septum purge to prevent bleed components from entering the system. The sample is then delivered to a sample splitter for analysis on two independent column modules.

Two Programmed Temperature Column Modules (PTCM) in parallel containing the resistively heated steel capillary chromatography column with necessary hardware, software and electronic control to enable temperature programming from 0.1°C to 5°C per second from 5°C above ambient to 350°C depending on the maximum temperature capability of the column material selected. Each column module is independently controlled by the method and can be any of the available column types.

Two independent Detector Modules incorporating either a micro Flame Ionization Detector (FID) or micro Thermal Conductivity Detector (TCD) with the necessary hardware, software and electronic control to provide detector temperature control, digital output signal and additionally for proper FID fuel supply pressure and auto-ignition.

The micro FID is a fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. The data rate is 100 Hz.

The micro TCD is a fully digital, universal detector. The TCD consists of a constant temperature filament that senses change in power required to hold the filament temperature constant when chromatographic components elute. The power measurement is used to determine the amount of the component eluting from the column. The data rate is 50 Hz.

CALIDUS is controlled with ChromPerfect chromatography data system fully integrated with LineUp running on a Windows PC.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518

Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: [info@falconst.net](mailto:info@falconst.net)

### SIGNIFICANCE AND USE

The Model 301 provides a simple ultra-fast analysis (10 to 50 times faster than conventional lab or process GC's) of fixed gases and hydrocarbons to C<sub>44</sub>. Use of two different PTCMs in parallel with their own detectors for example one TCD and one FID with the appropriate column material enables leveraging the sensitivity and selectivity differences for enhanced separations. Proper choice of columns enables wide boiling range and concentration ranges with a single GC.

The analyses are used for product specifications testing, product safety, environmental measurements, process control, catalyst protection, educational tools, spot checks of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**



# CALIDUS™ 301 Specifications (global patents pending)

December 1, 2010

## Ambient Environment

Operating Temperature Range: 0°C to 35°C  
Storage Temperature Range: -20°C to 60°C  
Relative Humidity Range: 0 to 100% (non-condensing)

## Power Requirements

Less than 300 watts peak power at startup, practical use < 200 Watts for gas or liquid analyses  
24 VDC supplied from external power supply, 100-240VAC using 50/60Hz AC

## Safety

General purpose, light industrial (lab instrument environment)  
CE Mark and Nationally Recognized Testing Laboratory (NRTL) certified pending (TUV Rheinland)

## Gas Supplies

50 PSIG, 99.995% H<sub>2</sub> or He at up to 250 ml/min, 50 PSIG zero air for FID operation

## Sample Requirements (via split/splitless injector with septum purge)

Air or gaseous samples at 0 to 50 PSIG at ambient temperature  
Membrane, SPME and static and dynamic headspace extracts  
Direct liquid injections neat or dilute organic solvents (DCM, hexane, MEK, toluene, methanol, CS<sub>2</sub> etc.)

## Dimensions

17" wide by 8.5" deep by 11" high, ~ 20 lbs  
Uninterrupted power supply and data acquisition computer external to the base unit

## Controls/Outputs

All functions and parameters via ethernet or RS-232 using ChromPerfect™ software  
Start analysis from keyboard or GC  
Set method from external computer using ChromPerfect software  
50-100Hz digitization (detector dependent) on each column, 24 bit resolution, auto zero on each run  
Trigger in and ready out signals plus an array of others via ChromPerfect

## Front Panel Displays

Temperature and pressure readings, function on/off, other  
Power on/off  
Status of analysis columns (isothermal, programming, cool down, ready, cycles run, other)

## Standard Equipment


Two capillary columns, 2m long, 100µm to 320µm ID, temperature programmable from 0.1 to 5°C per second from 5° above ambient to 350°C (maximum temperature software limited to be no greater than the limit for the columns installed, isothermal operation is available). Column modules are 2 meter columns in Mxt-1, Mxt-5, Mxt-1701, Mxt-Wax, Mxt-MoleSieve, Mxt-Alumina, < 320µm and various film thicknesses with others coming soon.  
Flame ionization or thermal conductivity (filament) detection  
Gas and liquid inlet for syringe injection of samples or automated gas and liquid sample valves available

## Performance (application dependent)

Repeatability of ± 1% RSD or better (area) and of ± 0.1% RSD or better (retention times)  
Analysis times for VOCs: can be <20 seconds and for SVOCs: can be <60 seconds  
Dynamic range: depends on detector used and application (FID typically 10<sup>5</sup>)

## Data Processing and Instrument Control

Note: computer system is integral and necessary component of the analysis system and includes the following:  
RS-232 or USB to RS-232 adapter, ethernet  
Windows XP or newer operating environment  
ChromPerfect software for dual column data acquisition via RS-232 serial or ethernet ports



Faster, Smaller, Smarter, Easier, Greener



**HRM**alytic +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA



## CALIDUS™ CS micro Gas Chromatograph

GC analysis for virtually any fixed gas and hydrocarbons up to C<sub>12</sub> for laboratory, at-line, transportable or online use

- **Upstream (E&P)**
- **Petrochemical/Chemical**
- **Food & Beverage**
- **Military**
- **Refining**
- **Pharmaceutical**
- **Medical**
- **Educational**



**CALIDUS Model CS** - comprised of 5 modules

A single Sample Processing Unit with a standard split/splitless injection port (1:1 up to 1:200) suitable gas and liquid samples via either syringe through the septum injections, optional gas, liquid or headspace auto-sampler, or automated sampling valves. The inlet includes septum purge to prevent bleed components from entering the system. The sample is then delivered to a column switching valve for analysis on two independent column modules.

Two Programmed Temperature Column Modules (PTCM) separated by a column switching valve containing the resistively heated steel capillary chromatography column with necessary hardware, software and electronic control to enable temperature programming from 0.1°C to 5°C per second from 5°C above ambient to 180°C (the maximum temperature for the valve, higher available on request) depending on the maximum temperature capability of the column material selected. Each column module is independently controlled by the method and can be any of the available column types.

Two independent Detector Modules incorporating either a micro Flame Ionization Detector (FID) or micro Thermal Conductivity Detector (TCD) with the necessary hardware, software and electronic control to provide detector temperature control, digital output signal and additionally for proper FID fuel supply pressure and auto-ignition.

The micro FID is a fully digital carbon/hydrogen bond detector using the hydrogen flame to burn the sample components. It uses an electrometer to sense the current changes in the flame cell due to chromatographic component elution. The data rate is 100 Hz.

The micro TCD is a fully digital, universal detector. The TCD consists of a constant temperature filament that senses change in power required to hold the filament temperature constant when chromatographic components elute. The power measurement is used to determine the amount of the component eluting from the column. The data rate is 50 Hz.

CALIDUS is controlled with ChromPerfect chromatography data system fully integrated with LineUp running on a Windows PC.



100 AEI Drive  
Fairlea, WV 24901  
P.O. Box 518  
Ronceverte, WV 24971

Phone: 304-647-5860  
Fax: 304-645-4006  
E-mail: [info@falconfast.net](mailto:info@falconfast.net)

### SIGNIFICANCE AND USE

The Model CS provides a simple ultra-fast analysis (10 to 50 times faster than conventional lab or process GC's) of fixed gases and hydrocarbons to C<sub>12</sub>. Use of two different PTCMs with their own detectors separated by the column switching valve with the appropriate column material enables leveraging the sensitivity and selectivity differences for enhanced separations (for example heartcuts and backflushing). Proper choice of columns enables wide boiling range and concentration ranges with a single GC. Model CS is ideal for individual component speciation from other sample matrix components.


The analyses are used for product specifications testing, product safety, environmental measurements, process control, catalyst protection, educational tools, spot checks of fuels and many more.

**Faster, Smaller, Smarter, Easier, Greener**



# CALIDUS™ CS Specifications (global patents pending)

December 1, 2010



## Ambient Environment

Operating Temperature Range: 0°C to 35°C  
Storage Temperature Range: -20°C to 60°C  
Relative Humidity Range: 0 to 100% (non-condensing)

## Power Requirements

Less than 300 watts peak power at startup, practical use < 200 Watts for gas or liquid analyses  
24 VDC supplied from external power supply, 100-240VAC using 50/60Hz AC

## Safety

General purpose, light industrial (lab instrument environment)  
CE Mark and Nationally Recognized Testing Laboratory (NRTL) certified pending (TUV Rheinland)

## Gas Supplies

50 PSIG, 99.995% H<sub>2</sub> or He at up to 250 ml/min, 50 PSIG zero air for FID operation

## Sample Requirements (via split/splitless injector with septum purge)


Air or gaseous samples at 0 to 50 PSIG at ambient temperature  
Membrane, SPME and static and dynamic headspace extracts  
Direct liquid injections neat or dilute organic solvents (DCM, hexane, MEK, toluene, methanol, CS<sub>2</sub> etc.)

## Dimensions

17" wide by 8.5" deep by 11" high, ~ 20 lbs  
Uninterrupted power supply and data acquisition computer external to the base unit

## Controls/Outputs

All functions and parameters via ethernet or RS-232 using ChromPerfect™ software  
Start analysis from keyboard or GC  
Set method from external computer using ChromPerfect software  
50-100Hz digitization (detector dependent) on each column, 24 bit resolution, auto zero on each run  
Trigger in and ready out signals plus an array of others via ChromPerfect



## Front Panel Displays

Temperature and pressure readings, function on/off, other  
Power on/off  
Status of analysis columns (isothermal, programming, cool down, ready, cycles run, other)

## Standard Equipment

Two capillary columns, 2m long, 100µm to 320µm ID, temperature programmable from 0.1 to 5°C per second from 5° above ambient to 350°C (maximum temperature software limited to be no greater than the limit for the columns installed, isothermal operation is available). Column modules are 2 meter columns in Mxt-1, Mxt-5, Mxt-1701, Mxt-Wax, Mxt-MoleSieve, Mxt-Alumina, < 320µm and various film thicknesses with others coming soon. The column switching valve can be plumbed in several ways but always between the two column modules.  
Flame ionization or thermal conductivity (filament) detection  
Gas and liquid inlet for syringe injection of samples or automated gas and liquid sample valves available

## Performance (application dependent)

Repeatability of ± 1% RSD or better (area) and of ± 0.1% RSD or better (retention times)  
Analysis times for VOCs: can be <20 seconds and for SVOCs: can be <60 seconds  
Dynamic range: depends on detector used and application (FID typically 10<sup>5</sup>)

## Data Processing and Instrument Control

Note: computer system is integral and necessary component of the analysis system and includes the following:  
RS-232 or USB to RS-232 adapter, ethernet  
Windows XP or newer operating environment  
ChromPerfect software for dual column data acquisition via RS-232 serial or ethernet ports

Faster, Smaller, Smarter, Easier, Greener



**HROM**alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA