

# Improving Method Performance through Fast LC

C. Vernon Bartlett, B.A.; Terrence S. Reid,  
M.S.; Rebecca E. Wittreg, Ph.D

Restek Corporation  
[www.restekcorp.com](http://www.restekcorp.com)



**HR**OMalytic +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

12/13

Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

# Abstract

The analysis time needed for many separations can be drastically reduced by the use of fast HPLC techniques. Several separations were converted using fast HPLC techniques. Analysis times for some separations that previously took over 35 minutes were reduced to less than 12 minutes with improvement in selectivity between the components. Qualitative TLC techniques can be converted to truly quantitative HPLC. Since columns employed in the fast LC analysis were typically less than 100mm in length, the reduction in analysis time also resulted in increased sensitivity due to reduction of band spreading.

In addition to improving performance through reduction of particle size and column length, performance gains may also be realized by using columns equipped with an appropriately optimized and highly selective stationary phase. These phases allow improved separation of the components without the drastic increase in  $k$  that often results when reduction in mobile phase strength is used to improve selectivity.

# Fast LC Technique

- Highly selective stationary phase is desired to maximize alpha values.
- Elution of components is typically accomplished through the use of gradients to reduce retention of highly retained components.
- Simple resolution of methylene substitutions /additions may be accomplished isocratically.
- Good screening technique for unknowns.

# Fast LC Technique – Advantages




- Fast re-equilibration (when using gradients)
- Sensitivity improvements.
- Older qualitative techniques can be adapted to a highly automated quantitative technique.
- Allows potentially high increases in sample throughput.
- Great technique when performed by LC-MS
- Shorter analysis times reduce solvent consumption and waste.

# Fast LC Technique – Disadvantages

- Critical separations are more sensitive to extra-column volume (as post column reactors).
- Extremely selective stationary phase must be used to maximize selectivity – especially for structural isomers.
- May not be well suited to normal phase or ion pairing separations (with gradients).

# Fast LC improvement of USP TLC and HPLC Method


## Digitalis Extracts and Derivatives



Digitoxin:  $X_1=X_2=H$

Digoxin:  $X_1=H; X_2=OH$

Gitoxin:  $X_1=OH; X_2=H$



# Fast LC Separation of Digitalis Derivatives (3 minutes)

## Digitalis Extracts/Derivatives on Ultra Alkaloids Cartridge Column: Fast LC


**Column:**  
Catalog #: 30 x 4.0mm  
Dimensions:  
Particle Size: 3µm  
Pore Size: 100Å

**Conditions:**  
Mobile Phase: A: 100% Water  
B: Acetonitrile  
Time (min.) %B  
0.0 20  
1.5 20  
1.51 35  
3.0 35  
3.1 20

Flow: 2.0mL/min  
Temp.: 27°C  
Det.: UV @ 230nm

Peak List:	Conc.	Ret. Time (min.)
1. digoxigenin	100µg/mL	0.40
2. gitoxigenin	100µg/mL	0.80
3. digoxin	100µg/mL	1.10
4. gitoxin	~10µg/mL	2.20
5. digitoxin	100µg/mL	2.60

**Sample:**  
Inj.: 10µL  
Sample Diluent: water:acetonitrile (80:20 v/v)



# Advantages of Digitalis Fast LC over Current Methods

- Improved automation and analysis throughput.
- Reduction of analysis time – previously a 30cm length C18 column was required for resolution.
- Ability to analyze all materials by HPLC.
- Ability to precisely quantitate materials vs TLC.
- Highly selective stationary phase.
- Technique can also be applied to purification and analysis of digoxin labeled materials for biological activity.
- Perfect technique for use in high speed cleaning validations.



# Fast LC Analysis of Carbamates




Figure 3-Structures of Commonly Analyzed Carbamates


# Carbamate Analysis using Standard HPLC Methodology (About 40 minutes)

## Carbamate Pesticides on Pinnacle Carbamate

Applications Note: LC\_0192

- Peak List:**
1. aldicarb sulfoxide
  2. aldicarb sulfone
  3. oxamyl
  4. methomyl
  5. 3-hydroxycarbofuran
  6. aldicarb
  7. propoxur
  8. carbofuran
  9. carbaryl
  10. 1-naphthol
  11. methiocarb

**Column:** Pinnacle Carbamate  
**Catalog#:** 9173575  
**Dimensions:** 250 x4.6mm  
**Particle size:** 5µm  
**Pore size:** 120Å  
**Solvent A:** Water  
**Solvent B:** Methanol  
**Gradient Program:**  
20-68% B; 0-20 min.  
68-79% B; 0-34 min.  
**Flow Rate:** 1.0mL/min  
**Detection:** Fluorescence  
**Ex. 339**      **Em. 445nm**



*Courtesy of Analytical Laboratories, Inc.  
Boise, Idaho*

# Fast LC Separation of Carbamates (About 13 minutes)

## Fast LC Separation of 11 Carbamates on Ultra Carbamate

### Peak List:

1. aldicarb
2. aldicarb sulfoxide
3. oxamyl
4. methomyl
5. 3-hydroxycarbofuran
6. aldicarb
7. propoxur
8. carbofuran
9. carbaryl
10. methiocarb
11. 4-bromo-3,5-dimethylcarbamate

### Sample:

Inj.: 5 $\mu$ L  
Conc.: 50 $\mu$ g/mL  
Solvent: methanol

### Restek standards:

Catalog# 32274 and 32273 mixed 50:50


### Column:

Catalog #: 9177355  
Dimensions: 50 x 4.6mm  
Particle size: 3 $\mu$ m  
Pore size: 100Å

### Conditions:

Mobile Phase: A: 90:10 water:methanol  
B: 90:10 methanol:acetonitrile  
Time (min): %B  
0 10  
10 90

Flow: 1.5mL/min  
Temp.: 27°C  
Det.: UV @ 220nm



# Fast LC Analysis of Carbamates with MS Detection

Table II. Experimental conditions for the LC/MS analysis of carbamate compounds.


## HPLC Conditions

Column: Ultra Carbamate, 100 mm x 4.6 mm, 3  $\mu$ m  
 Mobile Phase A: 90% water:10% methanol with 10mM ammonium formate  
 Mobile Phase B: 10% acetonitrile:90% methanol with 10mM ammonium formate  
 Gradient: 90%A:10%B to 10%A:90%B from 0-15 minutes  
 Inj. Volume: 10  $\mu$ L  
 Flow Rate: 0.75 mL/min to UV detector, 0.25 mL/min to MSD

## MSD Conditions


	Compound	Ion	Cone V
Detector:	1	223.3	25V
Mode:	2	207.3	18V
Capillary V:	3	237.2*	10V
Extractor:	4	163.2	15V
Ion Energy:	5	238.3	15V
Multiplier:	6	191.2	8V
Source Temp:	7	210.2	18V
Desolv. Temp:	8	222.3	22V
Gas Flow:	9	202.2	18V
	10	226.3	19V

\*Ammonium adduct (all other are [M+H]<sup>+</sup> ions)



# Fast LC Analysis – Carbamate Separation Loss

- Post Column Volumes produced by external reactors can be detrimental to critical separations – especially to smaller bore columns.



# Conventional Vanilla Flavoring Analysis (About 20 minutes with gradient)

## Vanillin and Ethyl Vanillin on Ultra C8

Applications Note: LC\_0148


Peak List:	Conc. (mg/mL)
1. vanillin	0.12
2. ethyl vanillin	0.04

**Sample:**  
Inj.: 10mL  
Solvent: 40% ethanol

**Column:** Ultra C8  
Catalog#: 9103565  
Dimensions: 150x4.6mm  
Particle size: 5µm  
Pore size: 100Å

**Conditions:**  
Mobile phases: A: 1.2% acetic acid  
B: methanol  
Gradient: Minutes %B  
0.0 20.0  
5.0 20.0  
15.0 40.0  
20.0 40.0  
25.0 20.0

Flow: 1.0mL/min.  
Temp.: 28°C  
Det.: UV @ 254nm



## Vanillin on Ultra C8

Application Note: LC\_0149


**Column:** Ultra C8  
Catalog#: 9103565  
Dimensions: 150x4.6mm  
Particle size: 5µm  
Pore size: 100Å

**Conditions:**  
Mobile phases: A: 1.2% acetic acid  
B: methanol  
Gradient: Minutes %B  
0.0 20.0  
5.0 20.0  
15.0 40.0  
20.0 40.0  
25.0 20.0

Flow: 1.0mL/min.  
Temp.: 28°C  
Det.: UV @ 254nm

**Peak List:**  
1. vanillin

**Sample:**  
Inj.: 10mL  
Conc.: 5% solution of vanilla extract  
Solvent: 40% ethanol



# Highly Selective Fast LC Vanillin Analysis (Less than 5 minutes and isocratic)


## Vanillin and Ethylvanillin on Pinnacle™ DB C18: Fast LC

Peak List:	Ret. Time (min.)
1. vanillin	2.00
2. ethylvanillin	4.20

**Sample:**  
Inj.: 5µL  
Conc.: 200µg/mL each component  
Sample Diluent: 4% v/v ethanol in water

**Column:**  
Catalog #: 9414554  
Dimensions: 50 x 4.0mm  
Particle Size: 5µm  
Pore Size: 140Å

**Conditions:**  
Mobile Phase: A: 0.2% v/v conc. phosphoric acid in water  
B: methanol  
Isocratic: 80%A: 20%B  
Flow: 1.5mL/min  
Temp.: 27°C  
Det.: UV @ 254nm



LC\_0251


## Vanilla Bean Extract on Pinnacle™ DB C18: Fast LC

**Column:**  
Catalog #: 9414554  
Dimensions: 50 x 4.0mm  
Particle Size: 5µm  
Pore Size: 140Å

**Conditions:**  
Mobile Phase: A: 0.2% v/v conc. phosphoric acid in water  
B: methanol  
Isocratic: 80% A: 20% B  
Flow: 1.5mL/min  
Temp.: 27°C  
Det.: UV @ 254nm

Peak List:	Ret. Time (min.)
1. Vanillin	2.00

**Sample:**  
Inj.: 5µL  
Conc.: unknown  
Sample Diluent: 20% v/v commercial extract in water



LC\_0252



## Vanillin Fast LC Analysis Advantages

- Takes advantage of improved methylene selectivity by the use of C18 phase versus a C8
- Analysis can be conducted in 3 to 5 minutes versus older 15 minute methods.
- High but not excessive selectivity of C18 phases toward simple methylene substitutions allow the use of an isocratic mobile phase.
- Reduction of column length also reduces the time needed for more hydrophobic analytes to elute from system.



# Conclusion

- Fast LC techniques applied upon highly selective stationary phases create a viable, precise quantitative alternative for analyses previously performed by Thin Layer Chromatography. In addition, these techniques can be used to improve method sensitivity, reduce solvent waste, and enhance laboratory throughput. Simplification of methods from gradient elution to isocratic elution can also occur when the proper stationary phase is used with a drastic reduction in analysis time, however, extraneous column volume caused by items as post column reactors can have a greater adverse effect when using fast columns.

# Acknowledgements

- The authors would like to thank Randy Romesberg, Larry Peters, and Rahul Patil of Restek Corporation for their participation in creating columns and hardware.