

Develop Residual Solvent Methods Quickly By “Outsourcing” to Restek’s Technical Support Team

By Rick Lake, Pharmaceutical Market Development Manager and
Amanda Rigdon, Innovations Chemist

- **Methods are customized to your specific solvent list.**
- **Restek Technical Support provides the recommended column format and analytical conditions in minutes using our Pro ezGC software.**
- **Ensures fast, cost-effective method development with highly accurate results.**

Developing custom residual solvent methods can be extremely time-consuming, as it often requires screening multiple columns over several scouting runs. Rather than conducting time- and labor-intensive testing in-house, simplify and speed up residual solvent method development by “outsourcing” to Restek’s Technical Support group. Restek has benchmarked the ICH Class 1, 2, and 3 residual solvents on our most popular columns and, through the use of our Pro ezGC software, we can model your specific solvent separation quickly and provide you with an accurate set of analytical conditions within hours. Here we demonstrate the effectiveness of a Pro ezGC-modeled solution.

In this example, a customer requested method conditions to separate 29 specific residual solvents. Restek Technical Support used Pro ezGC software to produce conditions for optimal separation of these analytes (Table I). Using the recommended Rxi®-624Sil MS column, we first matched up dead times to confirm an accurate flow rate, then, using only the conditions provided by Pro ezGC software, we separated the target compounds in just 15 minutes (Figure 1). When comparing the actual and predicted retention times, the average difference was only 0.05 minutes (Table II), confirming that the software accurately predicted conditions that optimally separated the components of a specific, complex solvent mix.

Let Restek help you with your residual solvent method today and see why we are the company that chromatographers trust. Simply contact Restek’s Technical Support group at support@restek.com with your specific solvent set and instrument model and we can quickly provide you with the operational parameters needed for complete chromatographic resolution.

For more information on residual solvent testing, visit www.restek.com/usp467



We can help you develop custom residual solvent methods quickly & easily!

Table I Operating conditions provided by Restek’s Technical Support team using Pro ezGC chromatography modeling software.

Column:	Rxi®-624Sil MS, 20 m, 0.18 mm ID, 1.00 μm (cat.# 13865)
Sample:	Custom residual solvents mix
Injection:	
Inj. Vol.:	0.5 μL split (split ratio 50:1)
Liner:	5mm split (cat.# 22973-200.1)
Inj. Temp.:	250 °C
Oven:	
Oven Temp:	40 °C (hold 5 min.) to 260 °C at 50 °C/min.
Carrier Gas	He, constant flow
Flow Rate:	1.35 mL/min.
Linear Velocity:	45 cm/sec.
Dead Time:	0.74 min. @ 40 °C
Detector	FID @ 250 °C
Data Rate:	50 Hz
Instrument:	Agilent/HP6890 GC

Rxi®-624Sil MS Columns (fused silica)

ID	df (μm)	temp. limits	length	cat. #
0.18mm	1.00	-20 to 300/320°C	20-Meter	13865

Figure 1 Excellent separation of a customer's solvent list was obtained in less than 15 minutes using conditions set by Restek's Pro ezGC software.




Table II Highly accurate methods for customer-specified lists can be developed quickly and cost-effectively using Restek's Technical Support and Pro ezGC software.

Solvent	Predicted RT (minutes)	Actual RT (minutes)	Difference (minutes)	Solvent	Predicted RT (minutes)	Actual RT (minutes)	Difference (minutes)
1. acetonitrile	2.083	2.056	0.027	16. pyridine	7.500	7.439	0.061
2. dichloromethane	2.218	2.185	0.033	17. toluene	7.673	7.587	0.086
3. trans-1,2-dichloroethene	2.478	2.431	0.047	18. ethylene glycol	7.786	7.748	0.038
4. hexane	2.788	2.707	0.081	19. 2-hexanone	8.310	8.260	0.050
5. nitromethane	3.551	3.497	0.054	20. dimethylformamide	8.645	8.597	0.048
6. cis-1,2-dichloroethene	3.643	3.581	0.062	21. chlorobenzene	8.934	8.882	0.052
7. tetrahydrofuran	4.017	3.948	0.069	22. ethylbenzene	9.031	8.985	0.046
8. chloroform	4.219	4.125	0.094	23. p-xylene	9.142	9.096	0.046
9. 1,1,1-trichloroethane	4.413	4.329	0.084	24. m-xylene	9.142	9.096	0.046
10. cyclohexane	4.532	4.427	0.105	25. o-xylene	9.452	9.412	0.040
11. 2-methoxyethanol	4.972	5.079	-0.107	26. N,N-dimethylacetamide	9.678	9.654	0.024
12. trichloroethene	6.196	6.095	0.101	27. 1-methyl-2-pyrrolidinone	11.467	11.466	0.001
13. methyl cyclohexane	6.49	6.373	0.117	28. 1,2,3,4-tetrahydronaphthalene	12.12	12.116	0.004
14. 1,4-dioxane	6.656	6.554	0.102	29. sulfolane	13.234	13.266	-0.032
15. 2-ethoxyethanol	7.044	6.982	0.062				
					Average Difference =	0.050	

PATENTS & TRADEMARKS

Restek patents and trademarks are the property of Restek Corporation. Other trademarks appearing in Restek literature or on its website are the property of their respective owners.



Lit. Cat.# PHTS1262

© 2010 Restek Corporation.

Restek U.S. • 110 Benner Circle • Bellefonte, PA 16823 • 814-353-1300 • 800-356-1688 • fax: 814-353-1309 • www.restek.com

Restek France • phone: 33 (0)1 60 78 32 10 • fax: 33 (0)1 60 78 70 90 • e-mail: restek@restekfrance.fr

Restek Ireland • phone: 44 2890 814576 • fax: 44 2890 814576 • e-mail: resteurope@aol.com

Thames Restek U.K. LTD • phone: 44 1494 563377 • fax: 44 1494 564990 • e-mail: sales@thamesrestek.co.uk

Restek GmbH • phone: +49 (0) 6172 2797 0 • fax: +49 (0) 6172 2797 77 • e-mail: info@restekgmbh.de

