

## Restek's Exceptionally Inert (Rxi®) Capillary GC Columns

As GC detectors become more sensitive, accurately quantifying low concentrations of target compounds becomes much more challenging. We developed the Rxi® column line specifically to improve low-level quantification. Our goal was to develop a superior column that had the highest inertness, lowest bleed, and greatest reproducibility of any column available.

The exceptionally low bleed levels of Rxi® columns improve low-level analysis by reducing detector noise. This improves signal-to-noise ratios for low-level compounds leading to more accurate and reproducible results. A highly inert column improves results for active compounds by preventing adsorption of target analytes in the system. The inertness of Rxi® columns allows analysis of acidic and basic compounds on the same column—often under the same conditions—without the peak tailing, that can skew results for low-level analytes.

Finally, consistent column performance is critical to low-level analysis. In developing the Rxi® columns, we wanted to guarantee reproducibility so customers would always receive a column that worked as well as their previous column. To guarantee column-to-column reproducibility we redesigned the entire manufacturing process and used strict quality specifications. Every Rxi® column is individually tested for coating efficiency, selectivity, film thickness, inertness, and bleed level. As a result, Rxi® columns offer the most consistent retention times and highest level of inertness on the market. The data presented here demonstrate the unmatched performance of the Rxi® columns; we guarantee these columns, engineered to improve low-level analyses, are the most reliable columns available.

### Low Bleed


Our bleed test is performed using a flame ionization detector with a compound marker to ensure the accuracy of the comparison. Column bleed was evaluated at 330°C and also at 350°C. As shown, the Rxi®-5ms column exhibits the lowest bleed of any column at both 330°C and 350°C (Figure 1). Note that at 350°C the variation in the bleed levels of the columns tested increases significantly. This increase is due to the difference in how the stationary phases are cross-linked by different manufacturers. As shown, the Crossbond® technology used by Restek in the Rxi® columns, results in a very stable stationary phase that does not degrade or bleed, compared to other columns on the market.

### Highly Inert

We used pyridine (a basic compound) and 2,4-dinitrophenol (an acidic compound) to evaluate the activity level of our Rxi® columns. In this test, if the column was too acidic, the pyridine peak would tail; whereas if the column was too basic, the 2,4-dinitrophenol peak would tail and exhibit a low response factor. The excellent peak symmetry shown in Figure 2 demonstrates the neutrality of the Rxi®-5ms column for both acidic and basic compounds. Additionally, while many other commercially available columns are not able to detect 2,4-dinitrophenol at 0.5ng on-column, the Rxi®-5ms column produces a response factor of 0.14.

**Figure 1** Rxi®-5ms columns have the lowest bleed among all major column brands.

Comparison of 30m x 0.25mm ID, 0.25µm columns at 330°C through 350°C; hydrogen carrier gas; flame ionization detection.



**Figure 2** Peak symmetry for pyridine or 2,4-dinitrophenol is excellent from an Rxi®-5ms column, even with 0.5ng on-column!

