

Basic Compounds Analysis

restek
innovation!**Rtx®-35 Amine** (midpolarity phase; Crossbond® 35% diphenyl/65% dimethyl polysiloxane)


- Application-specific columns for amines and other basic compounds, including alkylamines, diamines, triamines, ethanolamines, and nitrogen-containing heterocyclics.
- Stable to 220°C.

Active basic compounds that otherwise require derivatization, or an alternative analytical technique, can be analyzed on an Rtx®-35 Amine column. The tubing surface is chemically altered to reduce tailing of basic compounds, eliminating the need for column priming. An Rtx®-35 Amine column is ideal for analyzing a wide variety of basic compounds, but breakthrough technology also allows the analysis of neutral compounds, adsorptive compounds with oxygen groups susceptible to hydrogen bonding. Every Rtx®-35 Amine column is tested to ensure that it meets the requirements for analyzing ppm levels of amines, without priming, and to ensure low bleed at maximum operating temperature.

Rtx®-35 Amine Columns (fused silica)

(Crossbond® 35% diphenyl/65% dimethyl polysiloxane)

ID	df (µm)	temp. limits	15-Meter	30-Meter
0.25mm	0.50	0 to 220°C	11335	11338
	1.00	0 to 220°C	11350	11353
0.32mm	1.00	0 to 220°C	11351	11354
	1.50	0 to 220°C	11366	11369
0.53mm	1.00	0 to 220°C	11352	11355
	3.00	0 to 220°C	11382	11385

Sharp ethanolamine peaks, low bleed: Rtx®-35 Amine column.

1. monoethanolamine
2. diethanolamine
3. triethyleneglycol monomethylether
4. triethanolamine

Column: Rtx®-35 Amine, 30m, 0.32mm ID, 1.0µm (cat.# 11354)
 Sample: 500µg/mL ethanolamine standard in water
 Inj.: 1.0µL split (split ratio 10:1),
 cup splitter inlet liner (cat.# 20709)
 Inj. temp.: 300°C
 Carrier gas: helium, constant pressure
 Linear velocity: 40cm/sec. @ 50°C
 Oven temp.: 50°C (hold 0.50 min.) to 280°C @15°C/min.
 Det.: FID @ 300°C

Table of Contents for
Applications

see pages 518-519