

Simulated Distillation (C5-C44) Analysis

Rtx®-2887 (nonpolar phase; Crossbond® 100% dimethyl polysiloxane)

- Application-specific columns for simulated distillation.
- Stable to 360°C.

Rtx®-2887 columns' stationary phase, column dimensions, and film thickness have been optimized to exceed the resolution and skewing factor requirements currently specified in ASTM method D2887. Each column is individually tested to guarantee a stable baseline with low bleed and reproducible retention times. The Crossbond® methyl silicone stationary phase has increased stability compared to packed columns, ensuring stable baselines and shorter conditioning times.

Rtx®-2887 Column (fused silica)

(Crossbond® 100% dimethyl polysiloxane—for simulated distillation)

ID	df (μm)	temp. limits	10-Meter
0.53mm	2.65	-60 to 360°C	10199

MXT®-2887 Column (Siltek® treated stainless steel)

(Crossbond® 100% dimethyl polysiloxane—for simulated distillation)

ID	df (μm)	temp. limits	10-Meter
0.53mm	2.65	-60 to 400°C	70199




Jason Fisher
GC Column
Manufacturing Supervisor

similar **phases**

DB-2887, Petrocol EX2887

Negligible baseline rise for C5 to C44 hydrocarbons on an Rtx®-2887 column.



1. C5
2. C6
3. C7
4. C8
5. C9
6. C10
7. C11
8. C12
9. C14
10. C16
11. C18
12. C20
13. C24
14. C28
15. C32
16. C36
17. C40
18. C44

Column: Rtx®-2887, 10m, 0.53mm ID, 2.65μm (cat.# 10199)
 Sample: 1μL direct injection of 0.01-0.1 wt. % C5 to C44 hydrocarbon standard in carbon disulfide
 Inj. temp.: 360°C
 Det. temp.: 360°C
 Carrier gas: helium (constant flow)
 Linear velocity: 15mL/min. (112cm/sec.)
 Oven temp.: 35°C to 360°C @ 15°C/min. (hold 5 min.)