

## Explosives Analysis

restek  
innovation!

**Rtx®-TNT/Rtx®-TNT2** (proprietary Crossbond® phase)

- Application-specific columns for explosives in US EPA Method 8095.
- Low bleed—ideal for ECD analysis.
- Complete analysis in less than 20 minutes.
- Rtx®-TNT2 confirmation column provides 8 elution order changes under same conditions.
- Economical 3-packs.
- Stable to 310°C.

## please note

Polymer specially designed for explosives analysis.

We designed Rtx®-TNT and Rtx®-TNT2 columns specifically for analyses of nitroaromatic compounds by GC/ECD, such as the 16 analytes listed in US EPA Method 8095. They provide better resolution and higher thermal stability than any other currently recommended columns. Operate the Rtx®-TNT primary column and Rtx®-TNT2 confirmation column under identical GC oven temperature programs.

**Rtx®-TNT Columns** (fused silica)

ID	df (µm)	temp. limits	6-Meter/3-pk.
0.53mm	1.50	-20 to 300/310°C	12998


**Rtx®-TNT2 Columns** (fused silica)

ID	df (µm)	temp. limits	6-Meter/3-pk.
0.53mm	1.50	-20 to 300/310°C	12999

## US EPA Method 8095 explosives on Rtx®-TNT and Rtx®-TNT2 columns.

**Rtx®-TNT**


6m, 0.53mm ID, 1.50µm (cat.# 12998)



1. nitrobenzene
2. 2-nitrotoluene
3. 3-nitrotoluene
4. 4-nitrotoluene
5. nitroglycerine
6. 1,3-dinitrobenzene
7. 2,6-dinitrotoluene
8. 2,4-dinitrotoluene
9. 3,4-dinitrotoluene (IS)
10. 1,3,5-trinitrobenzene
11. trinitrotoluene
12. PETN
13. RDX
14. 4-amino-2,6-dinitrotoluene
15. 3,5-dinitroaniline
16. 2-amino-4,6-dinitrotoluene
17. tetryl

**Rtx®-TNT2**

6m, 0.53mm ID, 1.50µm (cat.# 12999)



Inj.: Direct injection using a 1mm Siltek® Uniliner® (cat.# 21052-214.1)  
 On-column conc.: est. 200-1000pg/compound. 8095 Calibration Mix A (cat.# 31607),  
 8095 Calibration Mix B (cat.# 31608), and 3,4-dinitrotoluene (cat.# 31452)  
 Oven temp.: 80°C (hold 1 min.) to 180°C @ 10°C/min. to 300°C @ 30°C/min. (hold 3 min.)  
 Inj. temp.: 250°C  
 Det.: ECD @ 330°C with anode purge  
 Dead time: 4.4 sec.  
 Head pressure: helium @ 3psi (20.7 KPa)  
 Flow rate: helium @ 17mL/min. @ 80°C



**Jarl Snider**  
R&D Chemist  
12+ years of service!