

**pHidelity® C18 Columns**

- Stable under extreme pH conditions.
- Patented technology protects the silica particle from dissolution.
- True C18 selectivity.

Application	Page #
Basic Drugs .....	556
Serotonin Reuptake Inhibitors .....	552

Restek is pleased to offer pHidelity® C18 columns, designed for analyses that require extreme pH conditions. Using technology patented by Selerity Technologies (US Patent 2005/0191503A1), pHidelity® silica-based columns have exceptional stability under aggressive pH conditions.

In pHidelity® columns, a polycarbosilane barrier layer protects the silica particle from extremely basic conditions. This layer, with multiple points of attachment to the silica particle, yields a modified surface with enhanced stability. A second layer is then attached, providing the functional group (C18). Using this approach to shield the silica surface, a highly durable stationary phase is created.

Even after 50 hours of exposure to base (pH 10) at 60°C, a pHidelity® column maintains its original column efficiency. (A decrease in column efficiency indicates the dissolution of support silica and loss from the column.) This exceptional performance ensures pHidelity® columns offer greater resolution under harsh mobile phase conditions and will continue to resolve your target compounds over a long column lifetime. A pHidelity® column also maintains a very consistent capacity factor,  $k'$ , after exposure to pH 10 and 60°C. Under extremely basic conditions, a decrease in retention time indicates loss of both support silica and bonded phase. Backpressure in the pHidelity® C18 column is very stable, indicating little or no degradation of the packing material.

The patented process used to protect the support silica in pHidelity® columns ensures more consistent retention times over the lifetime of the column.

**for more info**

For more information about pH stable pHidelity® HPLC Columns, review the *Restek Advantage 2007.01* at [www.restek.com](http://www.restek.com).

lit. cat.# 580133

**pHidelity® C18 Columns****Physical Characteristics:**

particle size: 3µm or 5µm  
pore size: 140Å

pH limit: up to 12  
temperature limit: 80°C

**Chromatographic Properties:**

Excellent stability under extreme pH conditions. True C18 selectivity in a silica-based stationary phase.

Length	2.1mm ID	3.2mm ID	4.6mm ID
	cat.#	cat.#	cat.#
<b>3µm Columns</b>			
30mm	9579332	9579333	9579335
50mm	9579352	9579353	9579355
100mm	9579312	9579313	9579315
150mm	9579362	9579363	9579365
<b>5µm Columns</b>			
150mm			9579565
250mm			9579575

new!

**ordering note**

To order a 2.1mm, 3.2mm, or 4.6mm ID column with a Trident Integral Inlet Fitting, add "-700" to the catalog number for the column.

Nominal additional charge

Example: 100mm x 4.6mm ID Ultra C18 column with Trident Integral Inlet Fitting: 9174315-700

Also order an XG-XF fitting (cat.#25026 or 25062), see page 337.

For guard cartridges for these columns, see page 339.

## Restek Capillary HPLC Columns

### Physical Characteristics:

#### Allure® Basix:

particle size: 3µm, spherical  
pore size: 60Å  
carbon load: 12%  
endcap: fully endcapped  
pH range: 2.5 to 7.5  
temperature limit: 80°C

#### Allure® PFP Propyl:

particle size: 3µm, spherical  
pore size: 60Å  
carbon load: 17%  
endcap: fully endcapped  
pH range: 2.5 to 7.5  
temperature limit: 80°C

#### Ultra C18:

particle size: 3µm, spherical  
pore size: 100Å  
carbon load: 20%  
endcap: fully endcapped  
pH range: 2.5 to 7.5  
temperature limit: 80°C

#### Allure® Biphenyl:

particle size: 3µm, spherical  
pore size: 60Å  
carbon load: 23%  
endcap: yes  
pH range: 2.5 to 7.5  
temperature limit: 80°C

#### Ultra Aqueous C18:

particle size: 3µm, spherical  
pore size: 100Å  
carbon load: 15%  
endcap: no  
pH range: 2.5 to 7.5  
temperature limit: 80°C

#### Viva C18:

particle size: 3µm, spherical  
pore size: 300Å  
carbon load: 9%  
endcap: yes  
pH range: 2.5 to 10  
temperature limit: 80°C

### Chromatographic Properties:

Available in a wide variety of stationary phases for acidic to basic compounds. Columns are suitable for small molecules using Allure® and Ultra phases and large biomolecules using Viva C18.

Length	0.3mm ID cat.#
<b>Allure® Basix 3µm Capillary Columns</b>	
50mm	916135B
150mm	916136B
<b>Allure® Biphenyl 3µm Capillary Columns</b>	
50mm	916635B
150mm	916636B
<b>Allure® PFP Propyl 3µm Capillary Columns</b>	
50mm	916935B
150mm	916936B
<b>Ultra Aqueous C18 3µm Capillary Columns</b>	
50mm	917835B
150mm	917836B
<b>Ultra C18 3µm Capillary Columns</b>	
50mm	917435B
150mm	917436B
<b>Viva C18 3µm Capillary Columns</b>	
50mm	951435B
150mm	951436B

new!



- High quality, Restek manufactured packing materials.
- Superior packing technology ensures rugged, reproducible columns.
- Wide range of phases and dimensions available—please inquire.

### ordering note

#### Looking for another stationary phase or capillary dimension?

Please contact Restek Technical Service at 814-353-1300 or 800-356-1688 (ext. 4) or support@restek.com to inquire.



#### HPLC Group

Kevin Davey, Larry Peters, Terry Cressman, Randy Romesberg,  
Reck Wittrig, Bruce Albright, Vernon Bartlett, Frank Dorman

# Pinnacle™ DB PFP Propyl and Biphenyl Columns

## new! Pinnacle™ DB PFP Propyl Columns (USP L43)

### Physical Characteristics:

particle size: 1.9µm, 3µm or 5µm, spherical      endcap: yes  
pore size: 140Å      pH range: 2.5 to 7.5  
carbon load: 6%      temperature limit: 80°C

Restek  
Manufactured  
Silica

Application	Page #
Benzodiazepines . . . . .	527

### Chromatographic Properties:

Pinnacle™ DB PFP Propyl is a unique pentafluorophenyl phase with a propyl spacer, and uses a highly base-deactivated spherical silica manufactured by Restek. This highly base-deactivated packing exhibits excellent peak shapes for a wide range of compounds, including nucleosides, nucleotides, and halogenated compounds.

Length	1.0mm ID cat.#	2.1mm ID cat.#	3.2mm ID cat.#	4.6mm ID cat.#
<b>1.9µm Columns</b>				
30mm		9419232		
50mm		9419252		
100mm		9419212		
<b>3µm Columns</b>				
30mm	9419331	9419332	9419333	9419335
50mm	9419351	9419352	9419353	9419355
100mm	9419311	9419312	9419313	9419315
150mm	9419361	9419362	9419363	9419365
<b>5µm Columns</b>				
30mm	9419531	9419532	9419533	9419535
50mm	9419551	9419552	9419553	9419555
100mm	9419511	9419512	9419513	9419515
150mm	9419561	9419562	9419563	9419565
200mm	9419521	9419522	9419523	9419525
250mm	9419571	9419572	9419573	9419575

## for more info

See page 306 for more information on our new 1.9µm Pinnacle™ DB columns.

HPLC COLUMNS

## new! Pinnacle™ DB Biphenyl Columns (USP L11)

### Physical Characteristics:

particle size: 1.9µm, 3µm or 5µm, spherical      endcap: yes  
pore size: 140Å      pH range: 2.5 to 7.5  
carbon load: 8%      temperature limit: 80°C

Restek  
Manufactured  
Silica

Application	Page #
Steroids . . . . .	555

### Chromatographic Properties:

Pinnacle™ DB Biphenyl is a unique reversed phase material that displays both increased retention and selectivity for aromatic and/or unsaturated compounds when compared to conventional alkyl and phenyl phases. Highly base-deactivated spherical silica manufactured by Restek. An excellent choice for the analysis of steroids, tetracyclines, drug metabolites, and other compounds that contain some degree of unsaturation.

Length	1.0mm ID cat.#	2.1mm ID cat.#	3.2mm ID cat.#	4.6mm ID cat.#
<b>1.9µm Columns</b>				
30mm		9409232		
50mm		9409252		
100mm		9409212		
<b>3µm Columns</b>				
30mm	9409331	9409332	9409333	9409335
50mm	9409351	9409352	9409353	9409355
100mm	9409311	9409312	9409313	9409315
150mm	9409361	9409362	9409363	9409365
<b>5µm Columns</b>				
30mm	9409531	9409532	9409533	9409535
50mm	9409551	9409552	9409553	9409555
100mm	9409511	9409512	9409513	9409515
150mm	9409561	9409562	9409563	9409565
200mm	9409521	9409522	9409523	9409525
250mm	9409571	9409572	9409573	9409575



**Cathy Hetrick**  
Northwest States Sales  
Representative  
3+ years of service!

## Allure® Silica Columns (USP L3)

### Physical Characteristics:

particle size: 3µm or 5µm, spherical  
 pore size: 60Å  
 endcap: no  
 pH range: 2.5 to 7.5  
 temperature limit: 80°C

### Chromatographic Properties:

Highly retentive phase for normal phase separations. Very high surface area, Type B silica packing.

Length	1.0mm ID	2.1mm ID	3.2mm ID	4.6mm ID
	cat.#	cat.#	cat.#	cat.#
<b>3µm Columns</b>				
30mm	9160331	9160332	9160333	9160335
50mm	9160351	9160352	9160353	9160355
100mm	9160311	9160312	9160313	9160315
<b>5µm Columns</b>				
30mm	9160531	9160532	9160533	9160535
50mm	9160551	9160552	9160553	9160555
100mm	9160511	9160512	9160513	9160515
150mm	9160561	9160562	9160563	9160565
200mm	9160521	9160522	9160523	9160525
250mm	9160571	9160572	9160573	9160575

## Allure® AK Columns

### Physical Characteristics:


particle size: 5µm  
 pore size: 60Å  
 endcap: yes  
 pH range: 2.5 to 7.5  
 temperature limit: 80°C

### Chromatographic Properties:

This highly retentive, highly selective phase, unique to Restek, was developed specifically for the analysis of aldehydes and ketones as DNPH derivatives. Allure® AK is a reversed phase HPLC material that has the unique ability to separate all thirteen carbonyl compounds specified in California Air Resources Board (CARB) Method # 1004, using a simple acetonitrile/water gradient, in less than 15 minutes. Other columns require long analysis times or the use of tetrahydrofuran.

Application	Page #
Carbonyls . . . . .	321, 528

Length	3.2mm ID cat.#	4.6mm ID cat.#
<b>5µm Columns with Trident Integral Inlet Fittings</b>		
200mm	9159523-700	9159525-700



**Frank Dorman**  
 Director of Technical Development  
 12+ years of service!



# Allure® Biphenyl and Organic Acids Columns

Application	Page #
Antidepressants	551
Cannabinoids	526
Corticosteroids	553, 554
NSAIDs	557
Steroids	554, 555

## free literature

### Allure® Biphenyl HPLC Columns

Download your free copy from [www.restek.com](http://www.restek.com)!

Flyer  
lit. cat.# 580015A

## Allure® Biphenyl Columns (USP L11)

### Physical Characteristics:

particle size: 3µm or 5µm, spherical  
pore size: 60Å  
carbon load: 23%

endcap: yes  
pH range: 2.5 to 7.5  
temperature limit: 80°C

### Chromatographic Properties:

Highly retentive and selective for aromatic and unsaturated compounds. Increased retention and selectivity, compared to phenyl phases. Excellent selectivity for steroids, tetracyclines, explosives, and other unsaturated compounds.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#		cat.#		cat.#		cat.#	
<b>3µm Columns</b>								
30mm	9166331		9166332		9166333		9166335	
50mm	9166351		9166352		9166353		9166355	
100mm	9166311		9166312		9166313		9166315	
<b>5µm Columns</b>								
30mm	9166531		9166532		9166533		9166535	
50mm	9166551		9166552		9166553		9166555	
100mm	9166511		9166512		9166513		9166515	
150mm	9166561		9166562		9166563		9166565	
200mm	9166521		9166522		9166523		9166525	
250mm	9166571		9166572		9166573		9166575	

## Allure® Organic Acids Columns

### Physical Characteristics:

particle size: 5µm, spherical  
pore size: 60Å

endcap: no  
pH range: 2.5 to 7.5  
temperature limit: 80°C

### Chromatographic Properties:

Allure® Organic Acids columns provide enhanced retention and selectivity for polar organic acids, allowing the separation to be performed on a single 30cm column. An Allure® Organic Acids column effectively resolves key organic acids such as tartaric and quinic acids, using the chromatographic conditions specified in AOAC method 986.13. Retention is stable and reproducible, even with the 100% aqueous mobile phase specified in the AOAC method.

Application	Page #
Carboxylic Acids	540
Fruit Juice Acids	541

Length	3.2mm ID		4.6mm ID	
	cat.#		cat.#	
<b>5µm Column</b>				
150mm	9165563		9165565	
250mm			9165575	
300mm			9165585	

Note: Other dimensions available on request.

## ordering note

To order a 2.1mm, 3.2mm, or 4.6mm ID column with a Trident Integral Inlet Fitting, add "-700" to the catalog number for the column.

Nominal additional charge

Example: 100mm x 4.6mm ID Ultra C18 column with Trident Integral Inlet Fitting: 9174315-700

Also order an XG-XF fitting (cat.#25026 or 25062), see page 337.

For guard cartridges for these columns, see page 339.



## Allure® PFP Propyl Columns (USP L43) Excellent Columns for LC/MS and ELSD

### Physical Characteristics:

particle size: 3µm or 5µm, spherical      endcap: fully endcapped  
pore size: 60Å      pH range: 2.5 to 7.5  
carbon load: 17%      temperature limit: 80°C

### Chromatographic Properties:

A pentafluorophenyl phase with a propyl spacer. Highly retentive for basic analytes. An excellent phase for separating nucleosides, nucleotides, purines, pyrimidines, halogenated compounds, β-blockers, and tricyclic antidepressants.

	1.0mm ID	2.1mm ID	3.2mm ID	4.6mm ID
Length	cat.#	cat.#	cat.#	cat.#
<b>3µm Columns</b>				
30mm	9169331	9169332	9169333	9169335
50mm	9169351	9169352	9169353	9169355
100mm	9169311	9169312	9169313	9169315
<b>5µm Columns</b>				
30mm	9169531	9169532	9169533	9169535
50mm	9169551	9169552	9169553	9169555
100mm	9169511	9169512	9169513	9169515
150mm	9169561	9169562	9169563	9169565
200mm	9169521	9169522	9169523	9169525
250mm	9169571	9169572	9169573	9169575

## Allure® Aqueous C18 Columns (USP L1) Excellent Columns for LC/MS and ELSD

### Physical Characteristics:

particle size: 5µm spherical      endcap: no  
pore size: 60Å      pH range: 2.5 to 7.5  
temperature limit: 80°C

### Chromatographic Properties:

Highly retentive and selective phase for separating polar analytes, including polar acidic compounds. Compatible with highly aqueous (up to 100%) mobile phases. Highly base deactivated. An excellent choice when analyzing a wide range of compounds, as in LC/MS screening methods.

	1.0mm ID	2.1mm ID	3.2mm ID	4.6mm ID
Length	cat.#	cat.#	cat.#	cat.#
<b>5µm Columns</b>				
30mm	9168531	9168532	9168533	9168535
50mm	9168551	9168552	9168553	9168555
100mm	9168511	9168512	9168513	9168515
150mm	9168561	9168562	9168563	9168565
200mm	9168521	9168522	9168523	9168525
250mm	9168571	9168572	9168573	9168575

Application	Page #
Antibiotics .....	549
Catecholamines .....	525
Cocaine, Ecgonine Methyl Ester .....	527
Nucleic Acid Bases .....	523
Opiates .....	526
Pesticides .....	533



**Melissa Decker**  
Customer Service  
Representative  
4+ years of service!

## ordering note

To order a 2.1mm, 3.2mm, or 4.6mm ID column with a Trident Integral Inlet Fitting, add "-700" to the catalog number for the column.

Nominal additional charge

Example: 100mm x 4.6mm ID Ultra C18 column with Trident Integral Inlet Fitting: 9174315-700

Also order an XG-XF fitting (cat.#25026 or 25062), see page 337.

For guard cartridges for these columns, see page 339.

# Ultra Quat Columns

restek  
innovation!

An Ultra Quat column and Ultra Quat Reagent Solution eliminate the need for ion pairing reagents in paraquat/diquat analysis.

## Ultra Quat Columns

### Physical Characteristics:

particle size: 5 $\mu$ m, spherical  
pore size: 100Å

pH range: 2.5 to 7.5  
temperature limit: 80°C

restek  
exclusive!

### Chromatographic Properties:

A retentive, high-purity, base deactivated reversed phase packing. Ideal for the analysis of paraquat and diquat or other quaternary amines when used with Ultra Quat Reagent Solution mobile phase additive (cat.# 32441).

Length	4.6mm ID cat.#
5 $\mu$ m Column	
150mm	9181565

Application	Page #
Paraquat, Diquat . . .	328, 534

### Ultra Quat Reagent Solution

Use with Ultra Quat HPLC column. Dilute to 1 liter water, per instructions.

In water, 20mL/bottle

cat. # 32441 (ea.)

### Paraquat & Diquat Calibration Mix

diquat dibromide                      paraquat dichloride  
1,000 $\mu$ g/mL each in water, 1mL/ampul  
cat. # 32437 (ea.)

## for more info

For more information about the Ultra Quat system, review the *Restek Advantage* 2004 vol. 3 at [www.restek.com](http://www.restek.com).

lit. cat.# 59051

## free literature

### Simple, Sensitive HPLC/UV Analysis for Paraquat and Diquat

Download your free copy from [www.restek.com](http://www.restek.com).

Applications Note  
lit. cat.# 580006

## Consistent resolution, retention times, and peak symmetry for paraquat and diquat on an Ultra Quat column.

