



Pro ezGC software will save you time and money by greatly enhancing your productivity and increasing sample throughput.



For Fast GC, Windows® NT, 2000, XP, Vista, or Windows® 7 (compatibility mode).

Pro ezGC Methods Development Software

- Optimize temperature and flow programs with a single analysis.
- Reduce analysis time and improve sample resolution.
- Model retention gap and guard column applications, including Restek Integra-Guard® columns.
- Optimize dual-column run conditions, columns in parallel or in series.

Take the guesswork out of selecting the best column and conditions for your GC analysis. Pro ezGC software accurately predicts separations on any capillary column, and is useful for selecting a column and conditions from a single GC run. Using your retention data, or the extensive library, you can automatically evaluate thousands of combinations of column dimensions, oven temperature programs, and carrier gas pressure programs to determine the best separation with the fastest analysis time.

Pro ezGC includes a master set of retention index libraries at no extra charge! These libraries contain more than 3,000 compounds analyzed on the most commonly used stationary phases, in ten application areas, including pesticides, PCBs, dioxins/furans, flavor and fragrance compounds, drugs of abuse, FAMES, semivolatile and volatile pollutants, petroleum hydrocarbons, and solvents and chemicals. The libraries permit computer simulation without entering actual laboratory data.

Description	qty.	cat.#	price
Pro ezGC Method Development Software CD-ROM	ea.	21487	\$201

Table I lists the Kovats retention indices for the more common stationary phases. Assigning a retention index to each probe listed provides a basis for comparing several stationary phases and their relative retention to one another for a set of molecular probes. For example, when Kovats indices are identical on two column phases, then the resulting separations will be identical. If, however, a Kovats value of one probe varies significantly from the value on another phase for the same probe, then the resulting compound elution order will differ. Thus, the Kovats indices are useful for comparing column selectivity for different types of compounds among different phases.

Table I Retention indices for Restek phases

Phase	Benzene	Butanol	Pentanone	Nitropropane
Rtx-1	651	651	667	705
Rtx-5/Rtx-5MS	667	667	689	743
Rtx-20	711	704	740	820
Rtx-1301/Rtx-624	689	729	739	816
Rtx-35	746	733	773	867
Rtx-200	738	758	884	980
Rtx-50	778	769	813	921
Rtx-1701	721	778	784	881
Rtx-65TG	794	779	825	938
Rtx-225	847	937	958	958
Stabilwax	963	1158	998	1230

Retention, k

The capacity of the column relates to how much material can be injected onto a column without adversely affecting peak shape. If the amount of a compound (mass) exceeds the capacity of a column (WCOT), the peak will front, which sometimes can look like a “shark fin”. The goal is to select a column with sufficient capacity such that peak shape will not suffer. Peak symmetry is typically used to calculate the degree of sample overload. There are two primary column-related dimensions that affect capacity, assuming the proper column phase was selected: column internal diameter (ID) and phase film thickness (μ).

When selecting column ID, consideration should include the type of injection, the detector being used, and the concentration of sample (amount on-column). The injection technique is an important consideration because the ID of the column may need to be selected based on whether a split, splitless, cool on-column injection, or other sample transfer to the column is being used. The second consideration is the detector and how much flow it can optimally work under. For example, some MS detectors can only handle column flow rates of up to 1.5 mL/min.; therefore, a 0.53 mm ID column, which requires higher flows for proper chromatography, is not an option for this detector. The third consideration is sample capacity of the column. If the concentration of the sample exceeds the column capacity, loss of resolution, poor reproducibility, and peak distortion will result. Table II shows several typical column characteristics for various column IDs.

Table II Typical characteristics for columns with the same phase ratio, such as 0.10 mm ID x 0.10 μ m and 0.18 mm ID x 0.18 μ m, etc.

Characteristic	Column ID					
	0.10mm	0.15mm	0.18mm	0.25mm	0.32mm	0.53mm
Helium Flow (@ 20cm/sec.)	0.16mL/min.	0.3mL/min.	0.3mL/min.	0.7mL/min.	1.2mL/min.	2.6mL/min.
Hydrogen Flow (@ 40cm/sec.)	0.32mL/min.	0.6mL/min.	0.6mL/min.	1.4mL/min.	2.4mL/min.	5.2mL/min.
Sample Capacity (max load per component)	<10ng	<40ng	<50ng	50–100ng	400–500ng	1000–2000ng
Theoretical Plates/Meter	8000	4000	3500	3200	2500	1800