

Guard Columns and Retention Gaps

Guard columns and retention gaps are widely used in gas chromatography. The concept of the guard column is to trap nonvolatile material at the head of the column, not allowing the material to reach the analytical column. The concept of the retention gap is to help focus the compounds transferred from the inlet to a small band at the head of the analytical column in order to reduce chromatographic peak broadening. Both concepts (trapping nonvolatile material and refocusing the target analytes) may take place when a piece of deactivated tubing is connected to an analytical column as in Figure 1.

did you know?

We test our guard columns/ transfer lines with a comprehensive test mix to ensure high inertness.




please note

For superior inertness, try our **Siltek® guard columns!**
See page 33 for details.

Having trouble making a leak-free connection? Try our “built in” **Integra-Guard® columns!**
See page 35 for details.

Figure 1 A guard/retention gap column connected to an analytical column




Analyte Focusing

There are two injection techniques where the retention gap is used to help focus target analytes at the beginning of the analytical column, cool on-column injection and split-less injection.

For cool on-column injection, the purpose of a retention gap is to help focus the sample components when introducing a liquid sample directly into the retention gap. The cool on-column injection is performed by inserting the syringe needle into the retention gap (this can be accomplished with a 0.53mm ID retention gap and a 26s gauge syringe) and transferring the liquid sample directly into the retention gap. The injection is made with the injector and column oven set below the boiling point of the solvent. As the solvent is evaporated, the volatile target analytes migrate in the solvent towards the analytical column, and the heavier analytes will be distributed over the retention gap. As the oven temperature increases, the target analytes vaporize and move unretained down the retention gap column until the compounds reach the liquid stationary phase of the analytical column. At this juncture, the target analytes are trapped/focused by the liquid phase forming a narrow injection band.

The retention gap may also be useful in hot vaporization injections when the transfer of the compounds from the inlet to the column does not form a focused band. Typical applications include water injections or injections using small ID columns, where split or tailing peaks would indicate an unfocused band. In these applications, the target analytes are trapped in a nonuniform or longitudinally diffuse band at the head of the retention gap (Figure 2a). As the oven temperature is increased, the solvent and target compounds are vaporized and move unretained through the retention gap (Figure 2b). When the target compounds come in contact with the stationary phase, they are refocused in a narrow band (Figure 2c), improving the chromatography.

Figure 2 Retention gaps are used to focus components in a tight band at the beginning of the analytical column.



a) Sample introduction: a liquid film of solvent and sample is deposited in the first length of capillary.

b) As oven temperature increases, the solvent evaporates and the target compounds elute unretained through the retention gap until they contact the analytical column.

c) When target compounds come in contact with the stationary phase, they are refocused on the analytical column, resulting in a narrow initial band width.