

Siltek® Coated Vials

For a highly inert surface, we recommend Siltek® deactivation for your vials:

- Maximizes inertness, minimizes sample breakdown—ideal for difficult matrices and reactive compounds.
- Inert over a wide sample pH range.
- Low bleed.
- Thermally stable.

Siltek® deactivation produces a highly inert glass surface that features high temperature stability, extreme durability, and low bleed. Try Siltek® deactivated vials for better recovery of sample analytes.

For Siltek® vials, add the corresponding suffix number to the vial catalog number.

Qty.	Siltek Deactivation
100-pk.	-222 addl. cost
1000-pk.	-223 addl. cost

**Septum Selection Guide**

Materials	Compatibility	Incompatibility	Resealability	Max. Temp.
Red Rubber (synthetic)	acetone, alcohols, DMF, DMSO, ether	ACN, benzene, chloroform, heptane, hexane, pyridine, THF, toluene	very good	90 °C
PTFE/ Natural Rubber	PTFE: resistance until punctured Rubber: acetone, ACN, alcohols, diethylamine, DMF, DMSO, phenol	aromatics, carbon disulfide, chlorinated solvents, hydrocarbon solvents	very good	90 °C
PTFE/Silicone PTFE/Silicone/PTFE	PTFE: resistance until punctured Silicone: acetone, alcohols, DMF, DMSO, ether	ACN, benzene, chloroform, heptane, hexane, pyridine, THF, toluene	very good	205 °C
Polyethylene	Good resistance to solvents and weak acids or bases. Unreactive with most chemicals, but some solvents cause softening or swelling.	hydrocarbon solvents	one-time use	175 °C
Gray Chlorobutyl	acids or bases, water solutions, buffer solutions, oxygenated solvents, vegetable oils	aliphatic or aromatic hydrocarbons, halogenated solvents, mineral oils, strong acids	very good	100 °C

Abbreviations: ACN = acetonitrile, DMF = dimethylformamide, DMSO = dimethylsulfoxide, THF = tetrahydrofuran

NOTE: This chemical resistance chart is intended only as a guideline. It does not cover all compounds or all solvents. Tests were done at room temperature on pure, single solvents, and there is no data on solvent combinations. Always confirm the compatibility of your vial, closure, and chemical combination prior to sample preparations.

instrument reference**For 2.0 mL, 11 mm Crimp-Top Vials**

The following instruments are compatible with 2.0 mL, 11 mm crimp-top vials. For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #	Manufacturer	Instrument/Model #	Manufacturer	Instrument/Model #
A.I.M.	CPS -100, 200	Gilson	231 XL, 232 XL, 233 XL, Asted XL, Aspec XL	Sedere	
A.I.	42 vial tray	Gynkotech	Gina, others	Shimadzu	SIIL-10A, 10Ai, 10Axl, AOC 20i
Agilent	1042, 1050, 1080, 1082, 1084, 1090, 1100, 5890, 6850, 6890, 7670A, 7671A, 7672, 7673A/B, 7683, 7890, 8042	Hitachi	AS-2000, AS-6000	Siemens	AS 32, AS 200
Alcott	738	IBM		Spark Holland	SPH 125, Marathon, Promis, Triathlon, Midas
Alltech	570, 580 (standard tray)	ICI	other than LC 1600	Spectra-Physics	8875, 8880
Altex		Infochroma		Talbot	
AMS	42	Jasco	851-AS, AS-950	Thermo Scientific	TRACE GC 2000, AS2000
Antek	736 Unisampler, 738	Kipp		TosoHaas	TSK-6080, AS-8010, AS-8020
ASC		Kontron	360	Tosca 1	
Carlo Erba/	42 vial tray	LDC	Marathon, Promis, other than 713	Tracor	770, 771, 772
Fisons	AS 800	L.E.A.P. Technologies	CTC A200S, CTC A105S	TSP	8875, 8880, AS 100/1000, AS 300/3000
Carnegie	CMA-250/200	Magnus Scientific		Unicam	4710, LC-XP, 4247
Chrompack	CP 9000 GC Series	PerkinElmer	Autosystem GC, AI-1, Integral 4000, ISS 100, ISS 200, LC 600 42 vial tray, 420/B, 4900	United Technologies	
CTC	CTC A105S	Pharmacia LKB	2157-010	Varian/Rainin	Dynamax AI-IA, AI-200, Dynamax AI-3
Dani	ALS 39.80, 86.80	Precision Sampling	GC111, GC 311, LC 241-60	Varian	8100/8200, 9100/9090/9095, Marathon
Delsi		Phillips	LC-XP, 4247, 4710	Waters	Alliance 2690
Dynatech	42 vial tray	Polymer-Labs	GPC 110/210		
Fisons	42 vial tray, AS800	Pye	LCXP		
GBC		S.G.E.	M280D		
Gerstel					