


Polycyclic Aromatic Hydrocarbons
Rxi®-5Sil MS

Peak List	Retention Time
1. naphthalene	4.93
2. 2-methylnaphthalene	5.58
3. 1-methylnaphthalene	5.68
4. 2-fluorobiphenyl (SS)	5.93
5. acenaphthylene	6.45
6. acenaphthene	6.62
7. fluorene	7.12
8. phenanthrene	8.06
9. anthracene	8.11
10. fluoranthene	9.23
11. pyrene	9.45
12. p-terphenyl-d14 (IS)	9.61
13. benzo(a)anthracene	10.65
14. chrysene	10.69
15. benzo(b)fluoranthene	11.96
16. benzo(k)fluoranthene	12.00
17. benzo(a)pyrene	12.42
18. perylene-d12 (IS)	12.51
19. indeno(1,2,3-cd)pyrene	14.19
20. dibenzo(a,h)anthracene	14.23
21. benzo(ghi)perylene	14.65



Column: Rxi®-5Sil MS, 30m, 0.25mm ID, 0.25µm (cat.# 13623)
 Sample: PAH mix, 1µL of 0.005µg/mL (IS 2µg/mL)
 SV Calibration Mix #5 (cat.# 31011)
 1-methylnaphthalene (cat.# 31283)
 2-methylnaphthalene (cat.# 31285)
 2-fluorobiphenyl (cat.# 31091)
 Inj.: 1.0µL (5pg on-column concentration),
 4mm Drilled Uniliner® (hole near top) inlet liner w/wool (cat.# 21055-200.5),
 pulsed splitless: pulse 20psi @ 0.2 min., 60mL/min. @ 0.15 min.
 Inj. temp.: 300°C
 Carrier gas: helium, constant flow
 Flow rate: 1.4mL/min.
 Oven temp.: 50°C (hold 0.5 min.) to 290°C @ 25°C/min. to 320°C @ 5°C/min.
 Det.: MS
 Transfer line temp: 290°C
 Ionization: EI
 Mode: SIM

Single Ion Monitoring Program

Group	Time	Ion(s)	Dwell (ms)
1	4.00	128	100
2	5.25	142	100
3	5.80	172	100
4	6.25	152	100
5	6.90	166	100
6	7.60	178	100
7	8.75	202, 244	100
8	10.2	228	100
9	11.5	252, 264	100
10	13.5	276, 278	100



free literature

Accurately Quantify PAHs Down to 5 pg On-Column: GC/MS Sim Analysis with the New Rxi®-5Sil MS Column
 Polycyclic aromatic hydrocarbons (PAHs) are common environmental pollutants, affecting air, water, and soil quality. Although naturally occurring, human impact has created a steady increase in environmental levels of PAHs and their byproducts. PAHs are typically formed through the incomplete combustion of organic materials, such as wood, coal, and oil, but are also used in manufacturing of some medicines, plastics, and pesticides. Many PAHs are known or suspected carcinogens. The United States Environmental Protection Agency currently lists and mandates testing of the 16 PAHs they deem most hazardous.

Download your free copy from www.restek.com

Applications Note
 lit. cat.# EVAN1284