


TCDD Isomers
Rt®-2330



Column: Rt®-2330, 60m, 0.25mm ID, 0.20µm (cat.# 10726)
 Sample: 2.0µL splitless injection of TCDD isomers, approximately 0.5ng each isomer
 Splitless hold: 0.45 min.
 Oven temp.: 200°C (hold 1 min.) to 250°C @ 8°C/min. (hold 15 min.), to 275°C @ 15°C/min. (hold 5 min.)
 Inj./det. temp.: 275°C
 Carrier gas: hydrogen
 Linear velocity: 40cm/sec.
 ECD sensitivity: 21 kHz full scale

- 1. 1,4,7,8-TCDD
- 2. 2,3,7,8-TCDD
- 3. 1,2,3,4-TCDD
- 4. 1,2,3,7-TCDD
- 5. 1,2,3,8-TCDD
- 6. 1,2,7,8-TCDD
- 7. 1,2,6,7-TCDD

also **available**

For Rt®-Dioxin2 columns, designed specifically for dioxin/furan analysis, see **page 96**.

Get More!

Environmental
Solutions Online

www.restek.com/enviro

