



restek
innovation!

Basic Compounds Analysis

Rtx®-5 Amine Columns (fused silica)

(low polarity phase; Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

- Application-specific columns for amines and other basic compounds, including alkylamines, diamines, triamines, ethanolamines, and nitrogen-containing heterocyclics.
- Stable to 315 °C.

Active basic compounds that otherwise require derivatization, or an alternative analytical technique, can be analyzed on an Rtx®-5 Amine column. The tubing surface is chemically altered to reduce tailing of basic compounds, eliminating the need for column priming. An Rtx®-5 Amine column is ideal for analyzing a wide variety of basic compounds, but breakthrough technology also allows the analysis of neutral compounds, adsorptive compounds with oxygen groups susceptible to hydrogen bonding, or even weakly acidic compounds such as phenols. Every Rtx®-5 Amine column is tested to ensure that it exceeds the requirements for analyzing ppm levels of amines, without priming, and to ensure low bleed at maximum operating temperature.

similar phase

PTA-5, CP-Sil CB

also available


See **page 65** for Rtx®-35 Amine columns.

please note

We recommend using base-deactivated fused silica guard columns (**page 34**) and base-deactivated liners (**page 213**) with Rtx®-5 Amine columns.

ID	df	temp. limits	15-Meter	30-Meter
0.25mm	0.25µm	-60 to 300/315°C	12320	12323
	0.50µm	-60 to 300/315°C	12335	12338
	1.00µm	-60 to 300/315°C	12350	12353
0.32mm	1.00µm	-60 to 300/315°C	12351	12354
	1.50µm	-60 to 290/305°C	12366	12369
0.53mm	1.00µm	-60 to 290/305°C	12352	12355
	3.00µm	-60 to 280/295°C	12382	12385

Ethanolamines on an Rtx®-5 Amine column.



Column: Rtx®-5 Amine, 15m, 0.25mm ID, 0.50µm (cat.# 12335)
 Sample: 1.0µL split injection of ethanolamine mix in methanol
 On-column conc.: 34ng
 Oven temp.: 50°C (hold 2 min.) to 180°C @ 10°C/min. (hold 2 min.)
 Inj./det. temp.: 280°C/300°C
 Carrier gas: hydrogen
 Linear velocity: 43cm/sec. set @ 50°C
 FID sensitivity: 6.4 x 10⁻¹¹ AFS
 Split ratio: 58:1

Chromatogram Search Tool

Search by compound name, synonym, CAS # or keyword

www.restek.com/chromatograms

