


NEW!

Cocaine and Metabolites (TMS Derivatives)

Rxi®-5Sil MS (100 ng/mL)

Column Sample Rxi®-5Sil MS, 30 m, 0.25 mm ID, 0.25 μ m (cat.# 13623)Diluent: Butyl chloride
Conc.: 100 ng/mL

Injection

Inj. Vol.: 1 μ L splitless (hold 1 min.)
Liner: Single Gooseneck w/Wool (cat.# 22286-200.1)
Inj. Temp.: 250 °C
Purge Flow: 20 mL/min.

Oven

Oven Temp: 100 °C to 200 °C at 30 °C/min. to 300 °C at 15 °C/min.

Carrier Gas He, constant linear velocity

Linear Velocity: 40 cm/sec., 12.5 psi, 86.2kPa @ 100 °C

Detector

Mode: MS
SIM

Transfer Line

Temp.: 310 °C

Source Temp.: 250 °C

Solvent Delay

Time: 4 min.

Tune Type: PFTBA

Ionization Mode: EI

Instrument

Notes Shimadzu 2010 GC & QP2010+ MS


Samples were prepared as follows:

Standards brought to dryness under nitrogen, then 50 μ L BSTFA + 1%TMCS (cat.# 35606) added. 50 μ L pyridine was then added, and samples were incubated at 70°C for 30 min. After incubation, samples were diluted with butyl chloride.

* Used as derivatization check


Cocaine & Metabolites (TMS Derivatives)

Rtx®-5

Column: Rtx®-5, 30m, 0.25mm ID, 0.25 μ m (cat.# 10223)
Sample: 1.0 μ L split injection of cocaine and cocaine metabolites
Oven temp.: 150°C to 320°C @ 10°C/min.
Inj./det. temp.: 250°C/300°C
Carrier gas: helium
Linear velocity: 30cm/sec. set @ 50°C
FID sensitivity: 2.56 x 10¹⁰ AFS
Split ratio: 30:1

Opiates (TMS Derivatives)

Rtx®-5

Column: Rtx®-5, 30m, 0.25mm ID, 0.25 μ m (cat.# 10223)
Sample: 2.0 μ L split injection of opiates
Conc: 2,000ng/ μ L
Oven temp.: 200°C to 325°C @ 7°C/min.
Inj./det. temp.: 250°C/300°C
Det. type: MS
Ionization: EI
Carrier gas: helium
Mode: full scan
Linear velocity: 30cm/sec. set @ 200°C
Split ratio: 50:1

RESTEK

HROMalytic +61(0)3 9762 2034
ECHnology Pty LtdAustralian Distributors
Importers & Manufacturers
www.chromtech.net.au

11/12

www.restek.com 673