

PRIMARY APPLICATIONS

Aldehydes
CFCs
Explosives
Ketones
PAHs
Silanes
Unsaturated compounds

VB-Fluoro Capillary Columns *NEW*

- 100% bonded Fluorosilicone
- High thermal stability
- Unique selectivity

VB-Fluoro capillary columns feature unique selectivity created by high fluorine affinity to analyte lone pair electrons. This is coupled with thermal stability similar to low polarity phases such VB-1 and VB-5.

Low bleed characteristics make VB-Fluoro columns well suited for MS and ECD applications, and the high thermal stability allows their use as a complementary column for most high temperature applications

which commonly utilize low polarity stationary phases.

Primary applications include ketones, aldehydes, explosives, PAHs, silanes, CFCs, and unsaturated compounds.


VB-Fluoro columns are a good replacement for Rtx-200, DB-200, DB-210, and VF-200 columns.

VB-Fluoro

0.25 mm ID	df	Prod No
30 meters	0.25	CFS-N03025-025
0.53 mm ID		
30 meters	1.00	CFS-N03053-100

VB-Fluoro**ACRYLATE IMPURITIES**

Column: VB-Fluoro
30 m x .53 mm x 1.00 μ m
Column temp: 40°C for 5 min,
to 200°C @ 10°C/min
Carrier: Hydrogen at 35 cm/sec
Injector: Split 100:1
Detector: FID

VB-Fluoro**EXPLOSIVES**

Column: VB-Fluoro
30 m x .53 mm x 1.00 μ m
Column temp: Isothermal 200°C
Carrier: Hydrogen at 35 cm/sec
Injector: Split 100:1
Detector: FID