

Fast Analysis of Semi-Volatile Compounds: US EPA Method 8270

Gary Stidsen and Frank Dorman
Jarl Snider

Restek Corporation
www.restekcorp.com



HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Product Definition

Design a GC analysis for 8270 that increases laboratory sample throughput by:

- Decreasing analysis time
- Resolution of key analytes
- Utilizing run parameters that can be used for current mass specs including HP59XX.
- 106 compounds including internal standards and surrogates plus 30 chlorinated and organophosphorus pesticides

Classes of Compounds

8270 Calibration Mix #1

benzoic acid
4-chloro-3-methylphenol
2-chlorophenol
2,4-dichlorophenol
2,6-dichlorophenol
2,4-dimethylphenol
4,6-dinitro-2-methylphenol
2,4-dinitrophenol
dinoseb
2-methylphenol

3-methylphenol
4-methylphenol
2-nitrophenol
4-nitrophenol
pentachlorophenol
phenol
2,3,4,6-tetrachlorophenol
2,4,5-trichlorophenol
2,4,6-trichlorophenol

8270 Calibration Mix #2

aniline
benzidine
4-chloroaniline
3,3'-dichlorobenzidine
diphenylamine
2-nitroaniline

3-nitroaniline
4-nitroaniline
N-nitrosodimethylamine
N-nitrosodi-n-propylamine
pyridine

Classes of Compounds

8270 Calibration Mix #3

aramite

bis (2-chloroethyl) ether

bis (2-chloroethoxy) methane

bis (2-chloroisopropyl) ether

4-bromophenyl phenyl ether

chlorobenzilate

2-chloronaphthalene

4-chlorophenyl phenyl ether

1,2-dichlorobenzene

1,3-dichlorobenzene

1,4-dichlorobenzene

1,3-dinitrobenzene

hexachlorobenzene

hexachlorobutadiene

hexachlorocyclopentadiene

hexachloroethane

hexachloropropene

isodrin

kepone

pentachlorobenzene

pentachloronitrobenzene

1,2,4,5-tetrachlorobenzene

1,2,4-trichlorobenzene

Classes of Compounds

8270 Calibration Mix #4

acetophenone

azobenzene

benzyl alcohol

bis (2-ethylhexyl) phthalate

butyl benzyl phthalate

dibenzofuran

diethyl phthalate

dimethyl phthalate

di-n-butyl phthalate

di-n-octyl phthalate

2,4-dinitrotoluene

2,6-dinitrotoluene

ethyl methanesulfonate

isophorone

isosafrole (cis & trans)

methyl methanesulfonate

1,4-naphthoquinone

nitrobenzene

4-nitroquinoline-1-oxide

phenacetin

safrole

Classes of Compounds

8270 Calibration Mix #5

acenaphthene
acenaphthylene
anthracene
benzo(a)pyrene
benzo(ghi)perylene
benzo(a)anthracene
benzo(b)fluoranthene
benzo(k)fluoranthene
chrysene
dibenz(a,h)anthracene

fluoranthene
fluorene
ideno(1,2,3-cd)pyrene
1-methylnaphthalene
naphthalene
3-methylcholanthrene
2-methylnaphthalene
phenanthrene
pyrene

Classes of Compounds

8270 Calibration Mix #6

diallate (cis & trans)
dimethoate
disulfoton
famphur
methyl parathion


parathion
phorate
pronamide
thionazine
0,0,0-triethyl phosphorothioate

Organochlorine Pesticide Mix AB #1

aldrin
a-BHC
a-chlordane
b-BHC
4,4'-DDD
4,4'-DDE
4,4'-DDT
d-BHC
dieldrin
endosulfan I

endosulfan II
endosulfan sulfate
endrin
endrin aldehyde
endrin ketone
g-BHC (lindane)
g-chlordane
heptachlor
heptachlor epoxide (B)
methoxychlor

Current Analysis



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Analysis Conditions

30m, 0.25mm ID, 0.50 μ m Rtx-5MS (cat.# 12638)
2 μ L splitless injection. On-column concentration: 20ng.

Oven temp.: 40°C (hold 2 min.) to 300°C @
10°C/min. (hold 4 min.), to 330°C
@ 10°C/min. (hold 10 min.).

Inj. / det. temp.: 280°C / 300°C

Det. type: MS

Carrier gas: helium

Linear velocity: 32 cm/sec. set @ 40°C

Ionization: EI

Splitless hold: 0.5 min.

Scan Range: 35-500AMU

Peak Identification

- | | | | | | |
|-----|-------------------------------|-----|----------------------------|-----|-----------------------------|
| 1. | 2-fluorophenol (Surr) | 27. | 4-chloroaniline | 53. | 2,4,6-tribromophenol (Surr) |
| 2. | phenol-d5 (Surr) | 28. | hexachlorobutadiene | 54. | 4-bromophenyl-phenylether |
| 3. | phenol | 29. | 4-chloro-3-methylphenol | 55. | hexachlorobenzene |
| 4. | bis(2-chloroethyl)ether | 30. | 2-methylnaphthalene | 56. | pentachlorophenol |
| 5. | 2-chlorophenol-d4 (Surr) | 31. | hexachlorocyclopentadiene | 57. | phenanthrene-d10 (IS) |
| 6. | 2-chlorophenol | 32. | 2,4,6-trichlorophenol | 58. | phenanthrene |
| 7. | 1,3-dichlorobenzene | 33. | 2,4,5-trichlorophenol | 59. | anthracene |
| 8. | 1,4-dichlorobenzene-d4 (IS) | 34. | 2-fluorobiphenyl (Surr) | 60. | carbazole |
| 9. | 1,4-dichlorobenzene | 35. | 2-chloronaphthalene | 61. | di-n-butylphthalate |
| 10. | 1,2-dichlorobenzene-d4 (Surr) | 36. | 2-nitroaniline | 62. | fluoranthene |
| 11. | 1,2-dichlorobenzene | 37. | dimethylphthalate | 63. | pyrene |
| 12. | 2-methylphenol | 38. | acenaphthylene | 64. | terphenyl-d14 (Surr) |
| 13. | 2,2'-oxybis-(1-chloropropane) | 39. | 2,6-dinitrotoluene | 65. | butylbenzylphthalate |
| 14. | 4-methylphenol | 40. | 3-nitroaniline | 66. | 3,3-dichlorobenzidine |
| 15. | N-nitroso-di-n-propylamine | 41. | acenaphthene-d10 (IS) | 67. | benzo(a)anthracene |
| 16. | hexachloroethane | 42. | acenaphthene | 68. | bis(2-ethylhexyl)phthalate |
| 17. | nitrobenzene-d5 (Surr) | 43. | 2,4-dinitrophenol | 69. | chrysene-d12 (IS) |
| 18. | nitrobenzene | 44. | 4-nitrophenol | 70. | chrysene |
| 19. | isophorone | 45. | dibenzofuran | 71. | di-n-octylphthalate |
| 20. | 2-nitrophenol | 46. | 2,4-dinitrotoluene | 72. | benzo(b)fluoranthene |
| 21. | 2,4-dimethylphenol | 47. | diethylphthalate | 73. | benzo(k)fluoranthene |
| 22. | bis(2-chloroethoxy)methane | • | Fluorene | 74. | benzo(a)pyrene |
| 23. | 2,4-dichlorophenol | • | 4-chlorophenyl-phenylether | 75. | perylene-d12 (IS) |
| • | 1,2,4-trichlorobenzene | 50. | 4-nitroaniline | 76. | indeno(1,2,3-cd)pyrene |
| 25. | naphthalene-d8 (IS) | 51. | 4,6-dinitro-2-methylphenol | 77. | dibenzo(a,h)anthracene |
| 26. | naphthalene | | | | |

Fast Analysis Concerns

- Flow rate considerations for diffusion pumps (<1.3 mL/min)
- Closely eluting compounds w/same quantitation ions
 - phenol / aniline / bis(2-chloroethyl)ether
 - 1,3- & 1,4-dichlorobenzene
 - 2- & 1-methylnaphthalene
 - 2,4,6- & 2,4,5-trichlorobenzenes
 - phenanthrene / anthracene
 - benz(a)anthracene / chrysene
 - benzo(b)fluoranthrene / benzo(k)fluoranthrene

Reducing Analysis Time

- Important criteria to reducing run times
 - initial hold time
 - resolve early eluting compounds
 - eluting compounds on ramp rate vs isothermal
 - fast ramp rate through non critical areas


Initial Work

- Constant flow
- Maximize ramp rate
 - Start out with fast ramp rate and work backwards
 - Adjust oven temperature program so all compounds after initial hold elute on ramp rate
- Overload of high standard
 - Testing performed with 160ng due to potential overload causing coelution of closely eluting compounds

Resulting Run Conditions

- Constant flow rate @ 1.0 mL/min
- Temperature program:
 - 40°C (2 min)
 - 20°C/min
 - 290°C (0 min)
 - 2°C/min
 - 303°C (0 min)
 - 6°C/min
 - 330°C (1 min)

Rtx-5Sil MS (30m x 0.25mm ID, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA


Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA


Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)




Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.5 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13


Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Can We Go Faster?

- Thinner film columns
- Shorter columns
- Smaller ID

Rtx-5Si1 MS

(30m x 0.25mm ID, 0.25 μ m film)




HRMalytic +61(0)3 9762 2034 **Australian Distributors**
ECHnology Pty Ltd **Importers & Manufacturers** **12/13**
www.chromtech.net.au

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA


Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)




Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)



Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd


Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA


Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)




Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)




Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)



Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Separation of Critical Pairs

Rtx-5Sil MS (30m x 0.25mm, 0.25 μ m film)

indeno(1,2,3-cd)pyrene

19.82


19.91

benzo(g,h,i)perylene

20.45

dibenz(a,h)anthracene

19.91




HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Rtx-5 MS Chromatogram (30m x 0.25 mm ID, 0.25 μ m film)



HRMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Conclusion

- Faster analysis can be accomplished for 8270
- Keep compounds eluting on the temperature program ramp rate
- Faster elution of PAHs helps with column capacity due to lower k value