

The Analysis of Semi-Volatiles with Various Inlet Liner Deactivations

Gary Stidsen, Frank Dorman and
David Smith

Restek Corporation
www.restekcorp.com



GC Inlet Liner Deactivation is Critical for Demanding Analyses


- Chlorinated pesticides
- Amines (basic compounds)
- Environmental semi-volatiles

Criteria for Evaluating Inlet Liner Performance

- Breakdown of endrin and DDT
- Linear response for acids and bases
- Short conditioning times and low bleed
- Thermal stability
- Resistance to sample degradation


Pesticide Breakdown, On-Column

Rtx[®]-CLPesticides2
30m, 0.53mm ID, 0.42µm
On-column injection



1. TCX
2. endrin
3. 4,4'-DDT
4. methoxychlor
5. DCB

Raw Glass Inlet Liner



Types of Liner Deactivation Chemistry


- Pinpoint
- Polymeric
- Base or amine
- Siltek™ Deactivation

Siltek™ Deactivation

- Inertness for acids, bases, neutrals, pesticides (developed for pesticides analysis)
- Low bleed
- Thermal stability
- Durability (acids, bases, water)
- Regeneration (solvent sonication: explosives analysis and customer feedback)
- Developing Siltek II for SemiVols+


Siltek™ Inertness

1.5% Endrin Breakdown




Acid Response Splitless XTI Injection

Rtx[®]-5
30m, 0.25mm ID, 0.25 μ m
Splitless



1. 1,2-hexanediol
2. nitro-di-n-propylamine
3. benzoic acid
4. C14
5. 2,4-dinitrophenol
6. nitrophenol
7. nitroaniline
8. pentachlorophenol
9. carbazole

Thermal Stability at 330°C




XTI-5, 30m, 0.25mm ID, 0.5 μ m, splitless Siltek™ liner; Inj.: 330°C; Det.: 330°C;
Oven temp.: 40°C (hold 2 min.) to 190°C @ 6°C/min to 330°C @ 30°C/min. (hold 10 min.)

Siltek Sleeve Resistance to Acidic Matrix



Siltek Sleeve Resistance to Basic Matrix

Direct Injections on to Siltek Unlined



Chromatographic System

- Agilent 6890 GC
- Agilent 5973 GC/MS
- Restek Rtx-5Sil MS 30m x 0.28mm x 0.5um
- Constant flow rate @ 1.0 mL/min
- Temperature program:
 - 40°C (2 min); 20°C/min to 290°C (0 min); 2°C/min to 303°C (0 min); 6°C/min to 330°C (1 min)
- Splitless hold time 0.4min
- Single gooseneck 4mm ID liners

Classes of Compounds

8270 Calibration Mix #1

benzoic acid
4-chloro-3-methylphenol
2-chlorophenol
2,4-dichlorophenol
2,6-dichlorophenol
2,4-dimethylphenol
4,6-dinitro-2-methylphenol
2,4-dinitrophenol
dinoseb
2-methylphenol

3-methylphenol
4-methylphenol
2-nitrophenol
4-nitrophenol
pentachlorophenol
phenol
2,3,4,6-tetrachlorophenol
2,4,5-trichlorophenol
2,4,6-trichlorophenol

8270 Calibration Mix #2

aniline
benzidine
4-chloroaniline
3,3'-dichlorobenzidine
diphenylamine
2-nitroaniline

3-nitroaniline
4-nitroaniline
N-nitrosodimethylamine
N-nitrosodi-n-propylamine
pyridine

Classes of Compounds

8270 Calibration Mix #3

aramite

bis (2-chloroethyl) ether

bis (2-chloroethoxy) methane

bis (2-chloroisopropyl) ether

4-bromophenyl phenyl ether

chlorobenzilate

2-chloronaphthalene

4-chlorophenyl phenyl ether

1,2-dichlorobenzene

1,3-dichlorobenzene

1,4-dichlorobenzene

1,3-dinitrobenzene

hexachlorobenzene

hexachlorobutadiene

hexachlorocyclopentadiene

hexachloroethane

hexachloropropene

isodrin

kepone

pentachlorobenzene

pentachloronitrobenzene

1,2,4,5-tetrachlorobenzene

1,2,4-trichlorobenzene

Classes of Compounds

8270 Calibration Mix #4

acetophenone

azobenzene

benzyl alcohol

bis (2-ethylhexyl) phthalate

butyl benzyl phthalate

dibenzofuran

diethyl phthalate

dimethyl phthalate

di-n-butyl phthalate

di-n-octyl phthalate

2,4-dinitrotoluene

2,6-dinitrotoluene

ethyl methanesulfonate

isophorone

isosafrole (cis & trans)

methyl methanesulfonate

1,4-naphthoquinone

nitrobenzene

4-nitroquinoline-1-oxide

phenacetin

safrole

Classes of Compounds

8270 Calibration Mix #5

acenaphthene

acenaphthylene

anthracene

benzo(a)pyrene

benzo(ghi)perylene

benzo(a)anthracene

benzo(b)fluoranthene

benzo(k)fluoranthene

chrysene

dibenz(a,h)anthracene

fluoranthene

fluorene

ideno(1,2,3-cd)pyrene

1-methylnaphthalene

naphthalene

3-methylcholanthrene

2-methylnaphthalene

phenanthrene

pyrene

Classes of Compounds

8270 Calibration Mix #6

diallate (cis & trans)
dimethoate
disulfoton
famphur
methyl parathion

parathion
phorate
pronamide
thionazine
0,0,0-triethyl phosphorothioate

Organochlorine Pesticide Mix AB #1

aldrin
a-BHC
a-chlordane
b-BHC
4,4'-DDD
4,4'-DDE
4,4'-DDT
d-BHC
dieldrin
endosulfan I

endosulfan II
endosulfan sulfate
endrin
endrin aldehyde
endrin ketone
g-BHC (lindane)
g-chlordane
heptachlor
heptachlor epoxide (B)
methoxychlor


Criteria for Trace Analysis

EPA Method 8270C

- Linearity of response relative to ISTD 20 - 160ng/ μ L
- Minimum relative response factor (RRF) and % relative standard deviation (RSD)


	Min. RRF	Min. %RSD
2,4-dinitrophenol	0.050	15
4-nitrophenol	0.050	15
N-nitroso-di-n-propylamine	0.050	15
hexachlorocyclopentadiene	0.050	15

Sample Chromatogram




Standard Liner Response

Response Factors for SPCCs




Standard Siltek Liner Response

Response Factors for SPCCs




Experimental Siltek-A Response

Response Factors for SPCCs



Experimental Siltek-B Response

Response Factors for SPCCs




Relative Standard Deviation Comparison

	Standard	Siltek	Siltek A	Siltek B
DNP	13.6	35.1	24.7	36.1
4-NP	9.0	16.2	12.4	12.1
NdnPA	7.0	3.2	7.6	7.6
HCCP	14.8	24.5	19.1	7.4


Results of Other Test Probes - Standard Liner and Siltek Liner

Response Factors



Standard


Response Factors



Siltek


Results of Other Test Probes - Siltek A & B

Response Factors



Siltek A

Response Factors



Siltek B

Other Test Probes: RSDs

	Standard	Siltek	Siltek A	Siltek B
PCP	17.8	22.3	18.0	17.9
Acenaphthylene	4.0	9.7	4.3	2.8

Conclusions

- Current Siltek™ products allow a wide variety of analyses, except at acidic and basic extremes
 - standard liners for acidic (i.e., semivolts)
 - base-deactivated for amines (ethanolamines)
 - Siltek for pesticides, explosives, sulfurs, volatiles, petrochemical, high temp., caustic, etc.
- Progressing on Siltek II which will also allow highly efficient Semi-Volatile analyses
- Future - Siltek III to include amines (i.e., one surface for all analyses)