

# The Design of High Temperature Capillary Gas Chromatography Columns Based on Polydimethylsiloxane

Jarl Snider, D. J. Hotnisky, Kristi Sellers,  
Dinesh V. Patwardhan Ph. D.

Restek Corporation  
[www.restekcorp.com](http://www.restekcorp.com)



**HROM**alytic +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

**12/13**

Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

# Outline

Background



Definition of Terminology

Unique Attributes of High Temperature Column


Applications

Conclusions

# Polydimethylsiloxane




# Thermal Stability of PDMS





Adapted from Siloxane Polymers, ed. Clarson & Semlyen, 1993.

# Bleed

- Chromatography of “loose” silicones from the column



# Typical Chemical Structure of Bleed



# Challenge

- Since the decomposition temperature for PDMS is 343°C, it is a challenge to make a column that is stable and has low bleed at 380°C.


# Origin of Bleed

- Polymer Synthesis
- Oligomers that are “created” in a column’s lifetime



# Origin of Bleed...

- Polymer Synthesis



Adapted from Inorganic Polymers, Mark, Allcock, & West 1992.

# Polymer Synthesis


- Silicone Polymers can be synthesized using well known catalysts (KOH, HCl) under equilibration polymerization conditions.
- We use newer, more advanced catalysts, better synthetic techniques.

# Origin of Bleed...


- Polymer Synthesis
- Oligomers that are “created” in a column’s lifetime

# Origin of Bleed...


- Oligomers that are “created” in a column’s lifetime



# Enhancement of Thermal Stability by Using “Additives”




# Additives...



- A column with true “One” selectivity cannot be obtained by having silphenylene in the backbone.


# Stx™-1 HT

- The polymer or the phase is PDMS.



# Origin of Bleed...

deactivation



fused silica

**HROM**alytic +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au


**12/13**

Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA



# Origin of Bleed...

deactivation



fused silica

**HROM**alytic +61(0)3 9762 2034  
**ECH**nology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au

**12/13**


Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

# Deactivation of the Fused Silica Surface


- We are using an advanced deactivation called Siltek™. It is a deposition process, unlike silazane or silicone deactivation.




# Aroclor<sup>®</sup> 1242



# Aroclor<sup>®</sup> 1254



# Aroclor<sup>®</sup> 1260



**HROM**alytic +61(0)3 9762 2034 **Australian Distributors**  
Importers & Manufacturers **12/13**  
**ECH**nology Pty Ltd [www.chromtech.net.au](http://www.chromtech.net.au)  
Website NEW : [www.chromalytic.com.au](http://www.chromalytic.com.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

# Aroclor<sup>®</sup> Run Conditions

Column:	Stx-HT1	Splitless Hold Time:	1.0min.
Serial #:	215041	Split Vent Flow:	40cc/min.
Description:	30m, 0.32mmID, 0.10 $\mu$ m	Septa Purge:	5cc/min.
Sample:	Aroclor Solutions	Carrier Gas:	Helium
	Part #32009, 32011, 32012	Head Pressure:	8.5psi
		Column Flow Rate:	1.9cc/min.
Concentration:	400ppb	Linear Velocity:	31cm/sec.
		Detector:	ECD/310C
Solvent:	Hexane	Make up Gas Flow:	40cc/min.
Sample Size:	1.0 $\mu$ L		
		Temp. Program:	75°C(1.0min.) 15°C/min. to 150°C(0) 5°C/min. 300°C(10min.)
Instrument:	HP5890		
Injector:	Splitless/275°C		

# PCB Congeners

Column: Stx-HT1  
Serial #: 215041  
Description: 30m, 0.32mmID, 0.10µm


Sample: PCB Congener Standard #2  
Part #32294  
Concentration: 100ng  
Solvent: Iso-Octane  
Sample Size: 1.0µL

Instrument: HP5890  
Injector: Splitless/275°C  
Splitless Hold Time: 0.80min.  
Split Vent Flow: 40cc/min.  
Septa Purge: 5cc/min.  
Carrier Gas: Helium  
Head Pressure: 8.5psi  
Column Flow Rate: 1.9cc/min.  
Linear Velocity: 31cm/sec.  
Detector: ECD/310°C  
Make up Gas Flow: 40cc/min.


Temp. Program: 100°C(0.80min.) 15°C/min.  
to 325°C(10min.)

Elution Order:

1. BZ#28
2. BZ#52
3. BZ#101
4. BZ#118
5. BZ#153
6. BZ#138
7. BZ#180



# 10ppm Aliphatics



1. C24
2. C26
3. C28
4. C30
5. C36

30m, 3.32mmID, 0.01 $\mu$ m Stx<sup>TM</sup>-1HT (cat.# 11709)

Solvent: Hexane

Oven temp.: 30°C (hold 0.5 min.) to 240°C at  
15°C/min. to 380°C at 10°C/min.  
(hold 10 min.)

Det.: FID at 340°C

Inj.: splitless at 250°C


Carrier: hydrogen

Flow rate: 4.0cc/min. at 30°C


Make-up gas flow: 40cc/min.



# Bleed Comparison



# 50ppm Hydrocarbons



## Aliphatic hydrocarbons:

- |        |         |
|--------|---------|
| 1. C9  | 6. C18  |
| 2. C10 | 7. C19  |
| 3. C12 | 8. C20  |
| 4. C14 | 9. C22  |
| 5. C16 | 10. C24 |
|        | 11. C26 |
|        | 12. C28 |
|        | 13. C30 |

## Aromatic Hydrocarbons:

1. naphthalene
2. 2-methylnaphthalene
3. acenaphthylene
4. acenaphthene
5. fluorene
6. phenanthrene
7. anthracene
8. fluoranthene
9. pyrene
10. benzo(a)anthracene
11. chrysene
12. benzo(b)fluoranthene
13. benzo(k)fluoranthene
14. benzo(a)pyrene
15. indeno(1,2,3-cd)pyrene
16. dibenzo(a,h)anthracene
17. benzo(g,h)perylene

# Analysis Conditions

30m, 3.32mmID, 0.01 $\mu$ m Stx<sup>TM</sup>-1HT (cat.# 11709)

Oven temp.: 30°C (hold 0.5 min.) to 240°C at  
15°C/min. to 380°C at 10°C/min.  
(hold 10 min.)

Det.: FID at 340°C

Inj.: splitless at 250°C


Carrier: hydrogen

Flow rate: 4.0cc/min. at 30°C

Make-up gas flow: 40cc/min.

1.0 $\mu$ L injections of MA EPH Aliphatic Hydrocarbon  
standard (cat.# 31459, solvent: hexane) & MA EPH  
Aromatic Hydrocarbon standard (cat.# 31458, solvent:  
methylene chloride) @ 50ppm concentration.

# ASTM D2887



30m, 3.32mmID, 0.01 $\mu$ m Stx™-1HT (cat.# 11709)

Oven temp.: 30°C to 380°C at 10°C/min.

Det.: FID at 380°C

Inj.: cool on-column

Carrier: helium

Flow rate: 1.3cc/min. at 30°C

Linear velocity: 22cm/sec. at 30°C

Make-up gas flow: 45cc/min.

Head pressure: 6.0psi

1.	C7	13.	C32
2.	C8	14.	C36
3.	C9	15.	C40
4.	C10	16.	C44
5.	C11	17.	C46
6.	C12	18.	C48
7.	C14	19.	C50
8.	C16	20.	C52
9.	C18	21.	C54
10.	C20	22.	C56
11.	C24	23.	C58
12.	C28		


**HROM**alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

Australian Distributors  
Importers & Manufacturers  
www.chromtech.net.au


12/13

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

# XTI Mix on Rtx-1 & Stx-1HT



1. 1,2-hexanediol
2. nitro-di-N-propylamine
3. benzoic acid
4. C14
5. 2,4-dinitrophenol
6. nitrophenol
7. nitroaniline
8. pentachlorophenol
9. carbazole
10. C20
11. C21
12. C22
13. benzo(b)fluoranthene
14. benzo(k)fluoranthene



# Summary

- Definitions
- The design of low bleed & stable column system was discussed
- Several applications were shown