

The Analysis of Trace Level Sulfurs in Beverage Grade CO₂

Barry Burger, David Smith
Restek Corp.
110 Benner Circle
Bellefonte, Pa. 16823

Project Objective

- Develop a robust micro-packed column for the analysis of trace level sulfurs in CO₂ with the following critical characteristics:
 - Thermal stability (310°C)
 - Inert (sub ppb levels)
 - Low bleed (< 20pa @ 310°C)
 - High sample capacity (15,000 ng)
 - Rapid analysis
 - Column longevity

Porous Polymer Column Optimization

- Porous polymer surface modification for the maximum degree of inertness
- Particle classification resulting in a tighter sample bandwidth
- Must resolve SO_2 from all other target sulfur compounds (TSC)

Column tubing considerations

- Low ID RMS value resulting in increased theoretical plates/efficiency
- Sulfinert deactivation resulting in state-of-the-art inertness
- Flexibility, easily installs in any GC
- Unreactive to sulfurs @ 20 ppbv levels


Advantages of the Packed Sulfinert™ column

- No cryogenic oven cooling required
- Highly breakable glass column eliminated
- Critical resolution of SO₂, H₂S & COS
- Column has the flexibility to be installed in any model GC
- Cost effective

GC system requirements

- Detection limit of < 0.02 ppm for SO_2 plus the other TSC
- Method of detection: SCD, HECD, FPD or PFPD detector
- All wetted sample pathways Sulfinert™ passivated for maximum inertness


Sulfinert™-Treated Inlet System



Chromatographic Parameters for the Rt-XL Sulfur Column

- 1 meter x 0.75mm ID micro-packed column with Sulfinert™ deactivation
- 10mL/min. helium
- 60°C > 260°C @ 15°C/min. > hold 10 min.
- GC/SCD (Sievers) Detection


Block diagram of Analytical System



Sampling system integrity for CO₂ liquid sample acquisition


- Passivated vessel to collect the liquid CO₂, conforming to the International Society of Beverage Technologists (ISBT) procedure 1.0
- Connections between the bulk CO₂ and the sample vessel must be passivated to prevent adsorption of the sulfur compounds
- What are the options, and how are the ISBT criteria adhered to?

Sample Acquisition System



Sulfinert™ Vessel Inertness

Sulfur Stability @ 1.5ppbv




Response Factors @ 1.5ppbv

		Day 0	Day 1	Day 2	Day 3	Day 6		
		RRF/DMS	RRF/DMS	RRF/DMS	RRF/DMS	RRF/DMS		% RSD
Hydrogen Sulfide		0.62	0.9	0.85	0.77	0.7		14.6
Carbonyl Sulfide		0.42	0.53	0.6	0.5	0.62		15.07
Methyl Mercaptan		1	1	1.1	1.2	1		8.43
Ethyl Mercaptan		0.42	0.39	0.77	0.77	0.35		39.15
Dimethyldisulfide		1.2	1.3	1.3	1.2	0.8		17.87


Beverage Sampling Data

- Beverage grade CO₂ blank
- 20ppb sulfur standard in beverage grade CO₂
- 20 ppb SO₂ in beverage grade CO₂ standard
- Headspace analysis of beer
- Analysis of a hard lemon alcoholic beverage
- Analysis of a top brand of cola


Beverage Grade CO₂ blank




20 ppb Sulfur Standard in Beverage Grade CO₂




20 ppb SO₂ in CO₂ Standard




Brand "F" Beer Headspace Sample



Hard Lemon Beverage Headspace Sample



Top Brand Cola Headspace Sample



Conclusion

- The Rt-XL Sulfur micro-packed column is a robust, low cost, rapid turn-around analytical tool for ppb level analysis of sulfur compounds.
- Sulfinert™ treatment for steel surfaces is unsurpassed for low-level ppb containment and transfer of highly reactive sulfur compounds.

Acknowledgements

- Seivers Instrument Inc. for their cooperation.
- Lou D'agostaro of DCG Partnership 1 Ltd. of Pearland Texas for supplying us with the sulfur standards.

For More Information:

Author e-mail: Barry@restekcorp.com

Restek Corporation

110 Benner Circle

Bellefonte, PA 16823

(800) 356-1688