


Analysis of Volatile Organic Compounds by Purge and Trap using GC/MS and a unique Stationary Phase

Restek Corporation
www.restekcorp.com



GC/MS Volatiles Analysis EPA Methodology

- ◆ EPA Method 8260B GC/MS
- ◆ Over 100 Organic Compounds
- ◆ Dynamic Headspace Sampling
- ◆ Extracted Ion Quantitation



Dynamic Headspace Sampling – Purge & Trap

- ◆ Purge & Trap increases sensitivity compared to other introduction methods
- ◆ However, purging analytes from the sample to the trap, then desorbing the trap, can result in broad peaks for the most volatile compounds.
- ◆ Column design needs to take into account the early eluting peak shapes and isomeric pairs (ex. Chlorotoluenes, dichlorobenzenes).

Column Design Criteria

- ◆ Analysis Time
- ◆ Low Bleed
- ◆ Critical Resolution
- ◆ No Cryofocusing
- ◆ GC/MS Compatible
- ◆ Computer Aided Stationary Phase Design (CASPD) can help achieve these criteria




Computer Aided Stationary Phase Design

- ◆ Phase modeling involves comparing analyte retention times collected on different stationary phases.
- ◆ This data can be used to predict the optimum resolution of compounds with different phase combinations.
- ◆ For example, Rtx-1 and Rtx-35 can be modeled to predict a composition for the Rtx-502.2 stationary phase.


Modeling

	Rtx-1	Rtx-502	Rtx-35
Chlorobenzene	17.79	18.57	19.27
1112te ethane	17.78	18.58	19.29
E benzene	18.26	18.78	19.25
m/p-xylene	18.48	18.90	19.27
o-ClF benzene	18.16	18.93	19.61
o-xylene	19.07	19.60	20.07
Styrene	18.93	19.63	20.25
Bromoform	18.58	19.96	21.19




Stationary Phases Used for Modeling

trifluoropropylmethyl
polysiloxane




dimethyl
polysiloxane




Stationary Phases Used for Modeling


diphenyl
polysiloxane



cyanopropylphenyl
polysiloxane





Predicted vs. Actual Using a 4 Dimensional Phase



Rtx-VMS


- ◆ Stationary Phase created using CASPD and optimized for GC/MS.
- ◆ Resolution and bleed levels were evaluated across a range of column dimensions and run conditions for maximum versatility.





Rtx[®]-VMS


Critical Pairs with Common Ions



Rtx[®]-VMS

Volatiles GC/MS Column


ACETATES ion 43



Time-->

Rtx[®]-VMS

Volatiles GC/MS Column




Column Variations for GC/MS – 0.18mm ID

- ◆ 0.18mm ID columns are attractive due to increased efficiency and resolution enhancements
- ◆ An Rtx-VMS 20m x 0.18mm x 1.0um column gives excellent resolution of EPA 8260B components in 10 minutes without cryogenic cooling
- ◆ Desorb flow rate set at 40ml/min for 1 min
- ◆ Column flow rate optimized at 1 ml/min

Volatile Organics
EPA Method 8260B
Rtx®-VMS

Innovations at Work

- 1st choice for use with dual purge & traps¹
- EPA recommended surrogate used.



20m, 0.18 mm ID, 1.00µm Rtx®-VMS (cat.# 49914)
 Compounds in at 10ppb in 5mL of RO water
 unless otherwise noted, ketones in at 2.5X
 Concentrator: Tekmar LSC-3100 Purge and Trap
 Trap: Vocarb 3000 (type K)
 Purge: 11 min. @ 40mL/min. @ ambient temperature
 Dry purge: 1 min. @ 40mL/min.
 Desorb preheat: 245°C
 Desorb: 250°C for 2 min., Flow 40mL/min.
 Bake: 260°C for 8 min.
 Interface: transfer line 0.53mm ID Silcosteel® tubing
 1:40 split at injection port. 1mm ID sleeve.
 Oven temp.: 50°C (hold 4 min.) to 100°C @ 18°C/min. (hold 0 min.)
 to 230°C @ 40°C/min. (hold 3 min.)
 Carrier gas: helium @ ~1.0mL/min. constant flow
 Adjust dichlorodifluoromethane to a retention time of 1.03 min. @ 50°C.
 Detector: HP 5973 MSD
 Scan range: 35-300amu

- | | | | |
|---|-------------------------------------|---------------------------------------|---------------------------------|
| 1. dichlorodifluoromethane | 26. ethyl acetate | 51. toluene | 76. <i>n</i> -propylbenzene |
| 2. chloromethane | 27. carbon tetrachloride | 52. pyridine (250ppb) | 77. 1,1,2,2-tetrachloroethane |
| 3. vinyl chloride | 28. methyl acrylate | 53. tetrachloroethene | 78. 2-chlorotoluene |
| 4. bromomethane | 29. propargyl alcohol (500ppb) | 54. 4-methyl-2-pentanone | 79. 1,3,5-trimethylbenzene |
| 5. chloroethane | 30. dibromofluoromethane (SMC) | 55. <i>trans</i> -1,3-dichloropropene | 80. 1,2,3-trichloropropane |
| 6. trichlorofluoromethane | 31. tetrahydrofuran | 56. 1,1,2-trichloroethane | 81. 4-chlorotoluene |
| 7. ethanol (2500ppb) | 32. 1,1,1-trichloroethane | 57. ethyl methacrylate | 82. <i>tert</i> -butylbenzene |
| 8. 1,1-dichloroethene | 33. 2-butanone | 58. dibromochloromethane | 83. pentachloroethane |
| 9. carbon disulfide (40ppb) | 34. 1,1-dichloropropene | 59. 1,3-dichloropropane | 84. 1,2,4-trimethylbenzene |
| 10. allyl chloride | 35. benzene | 60. 1,2-dibromoethane | 85. <i>sec</i> -butylbenzene |
| 11. methylene chloride | 36. pentafluorobenzene (IS) | 61. <i>n</i> -butyl acetate | 86. <i>p</i> -isopropyltoluene |
| 12. acetone | 37. <i>tert</i> -amyl-methyl ether | 62. 2-hexanone | 87. 1,3-dichlorobenzene |
| 13. <i>trans</i> -1,2-dichloroethene | 38. 1,2-dichloroethane | 63. 2-picoline (250ppb) | 88. 1,4-dichlorobenzene-d4 (IS) |
| 14. methyl <i>tert</i> -butyl ether | 39. isobutyl alcohol (500ppb) | 64. chlorobenzene-D5 (IS) | 89. 1,4-dichlorobenzene |
| 15. <i>tert</i> -butyl alcohol (100ppb) | 40. isopropyl acetate | 65. chlorobenzene | 90. <i>n</i> -butylbenzene |
| 16. diisopropyl ether | 41. trichloroethene | 66. ethylbenzene | 91. 1,2-dichlorobenzene |
| 17. 1,1-dichloroethane | 42. 1,4-difluorobenzene (SMC) | 67. 1,1,1,2-tetrachloroethane | 92. 1,2-dibromo-3-chloropropane |
| 18. acrylonitrile | 43. dibromomethane | 68. <i>m</i> -xylene | 93. nitrobenzene (250ppb) |
| 19. vinyl acetate | 44. 1,2-dichloropropane | 69. <i>p</i> -xylene | 94. hexachlorobutadiene |
| 20. allyl alcohol (250ppb) | 45. dibromodichloromethane | 70. <i>o</i> -xylene | 95. 1,2,4-trichlorobenzene |
| 21. ethyl- <i>tert</i> -butyl ether | 46. methyl methacrylate | 71. styrene | 96. naphthalene |
| 22. <i>cis</i> -1,2-dichloroethene | 47. <i>n</i> -propyl acetate | 72. bromoform | 97. 1,2,3-trichlorobenzene |
| 23. 2,2-dichloropropane | 48. 2-chloroethanol (2500ppb) | 73. isopropylbenzene | |
| 24. bromochloromethane | 49. <i>cis</i> -1,3-dichloropropene | 74. 4-bromo-1-fluorobenzene (SMC) | |
| 25. chloroform | 50. toluene-d8(SMC) | 75. bromobenzene | |

¹A.L. Hilling and G. Smith, *Environmental Testing & Analysis*, 10(3), 15-19, 2001.


Column Variations – 0.25mm ID

- ◆ Rtx-VMS 60m x 0.25mm x 1.4um column
- ◆ 60°C starting temperature
- ◆ 20:1 split injection
- ◆ Constant Flow @ 1.3 ml/min.
- ◆ 97 components in EPA 8260B in 18 minutes

Volatile Organics
EPA Method 8260B
Rtx®-VMS

60m, 0.25 mm ID, 1.40µm Rtx®-VMS (cat.# 19916)
 Compounds in at 10ppb in 5mL of RO water (unless noted)
 ketones 2.5X
 Concentrator: Tekmar LSC-3100 Purge and Trap
 Trap: Vocarb 3000 (type K)
 Purge: 11 min. @ 40mL/min. @ ambient temperature
 Dry purge: 1 min. @ 40mL/min.
 Desorb preheat: 245°C
 Desorb: 250°C for 2 min., Flow 26mL/min.
 Bake: 260°C for 8 min.
 Interface: transfer line 0.53mm ID Silcosteel® tubing
 1:20 split at injection port. 1mm ID sleeve.
 Oven temp.: 60°C (hold 2 min.) to 180°C @ 12°C/min. (hold 0 min.) to
 225°C @ 45°C/min. (hold 6 min.)
 Carrier gas: helium @ ~1.3mL/min. constant flow
 Adjust dichlorodifluoromethane to a retention time of 4.03 min. @ 60°C.
 Detector: HP 5973 MSD
 Scan range: 35-300amu

Innovations at Work
 60°C starting temperature



- | | | | |
|---------------------------------|--------------------------------|-----------------------------------|---------------------------------|
| 1. dichlorodifluoromethane | 26. ethyl acetate | 51. toluene | 76. 1,1,2,2-tetrachloroethane |
| 2. chloromethane | 27. methyl acrylate | 52. 4-methyl-2-pentanone | 77. bromobenzene |
| 3. vinyl chloride | 28. propargyl alcohol (500ppb) | 53. pyridine (250ppb) | 78. 1,3,5-trimethylbenzene |
| 4. bromomethane | 29. dibromofluoromethane (SMC) | 54. trans-1,3-dichloropropene | 79. 2-chlorotoluene |
| 5. chloroethane | 30. tetrahydrofuran | 55. ethyl methacrylate | 80. 1,2,3-trichloropropane |
| 6. trichlorofluoromethane | 31. carbon tetrachloride | 56. tetrachloroethene | 81. 4-chlorotoluene |
| 7. ethanol (2500ppb) | 32. 2-butanone | 57. 1,1,2-trichloroethane | 82. tert-butylbenzene |
| 8. 1,1-dichloroethene | 33. 1,1,1-trichloroethane | 58. dibromochloromethane | 83. 1,2,4-trimethylbenzene |
| 9. carbon disulfide (40ppb) | 34. 1,1-dichloropropene | 59. 1,3-dichloropropane | 84. pentachloroethane |
| 10. allyl chloride | 35. pentafluorobenzene (IS) | 60. n-butyl acetate | 85. sec-butylbenzene |
| 11. methylene chloride | 36. tert-amyl methyl ether | 61. 1,2-dibromoethane | 86. p-isopropyltoluene |
| 12. acetone | 37. benzene | 62. 2-hexanone | 87. 1,3-dichlorobenzene |
| 13. trans-1,2-dichloroethene | 38. isobutyl alcohol (500ppb) | 63. 2-picoline (250ppb) | 88. 1,4-dichlorobenzene-d4 (IS) |
| 14. tert-butyl alcohol (100ppb) | 39. 1,2-dichloroethane | 64. ethylbenzene | 89. 1,4-dichlorobutadiene |
| 15. methyl tert-butyl ether | 40. isopropyl acetate | 65. chlorobenzene-D5 (IS) | 90. n-butylbenzene |
| 16. diisopropyl ether | 41. 1,4-difluorobenzene (SMC) | 66. chlorobenzene | 91. 1,2-dichlorobenzene |
| 17. 1,1-dichloroethane | 42. trichloroethene | 67. 1,1,1,2-tetrachloroethane | 92. 1,2-dibromo-3-chloropropane |
| 18. acrylonitrile | 43. dibromomethane | 68. m-xylene | 93. nitrobenzene (250ppb) |
| 19. vinyl acetate* | 44. bromodichloromethane | 69. p-xylene | 94. hexachlorobutadiene |
| 20. allyl alcohol (250ppb) | 45. 1,2-dichloropropane | 70. o-xylene | 95. 1,2,4-trichlorobenzene |
| 21. ethyl-tert-butyl ether* | 46. methyl methacrylate | 71. styrene | 96. naphthalene |
| 22. cis-1,2-dichloroethene | 47. n-propyl acetate | 72. bromoform | 97. 1,2,3-trichlorobenzene |
| 23. 2,2-dichloropropane | 48. 2-chloroethanol (2500ppb) | 73. isopropylbenzene | |
| 24. bromochloromethane | 49. cis-1,3-dichloropropene | 74. 4-bromo-1-fluorobenzene (SMC) | |
| 25. chloroform | 50. toluene-d8 (SMC) | 75. n-propylbenzene | |

*These compounds can be resolved using a lower starting temperature.

EPA Method 8240

- ◆ Developed to monitor 73 compounds in hazardous waste
- ◆ Rtx-VMS 30m x 0.25mm x 1.4um is a good choice for this shorter list
- ◆ Analysis times of 14 minutes

Volatile Organics

EPA Method 8240 (8260 Short List)

Rtx®-VMS


Rtx®-VMS

30m, 0.25mm ID, 1.40µm
(cat.# 19915)

Instrumentation & Conditions Common to BOTH RUNS...

Carrier gas: 1.3mL/min. @ constant flow
 Concentrator: Tekmar LSC-3000 Purge and Trap
 Trap: Vocabr 3000 (type K) — see concentrations in 5mL/RO
 Purge: 11 min. @ 40mL/min. @ ambient temp.
 Dry purge: 1 min. @ 40mL/min. (MCS bypass)
 Desorb preheat: 245°C
 Desorb: 250°C for 2 min., flow 26mL/min.
 Bake: 260°C for 8 min.
 Interface to GC: transfer line 0.32mm ID Siltek fused silica
 1:20 split at injection port w/ 1mm ID sleeve.
 GC: HP 6890
 Detector: HP 5973 MSD

Innovations at Work




Top chromatogram:

Oven temp.: 40°C (hold 4 min.) to 90°C @ 16°C/min. (no hold)
 to 210°C @ 32°C/min. (hold 5 min.)
 Adjust dichlorodifluoromethane to a retention time of 2.27 min. @ 40°C.
 MS Scan Range: 35-300amu
 compound concentrations, by mix: (in 5mL of RO water)
 Compounds at 100ppb (cat.# 30213, 30004, 30006, 30011, 30042)
 Alcohols at 1ppm (cat.# 30214) except 2Cl ethanol at 10ppm.
 vinyl acetate at 500ppb (cat.#30216)
 8240 Nitrile Mix at 200ppb (cat.# 30215)
 8240 Mix 1A at 300ppb (cat.# 30217)
 8240 Mix 2A at 500ppb (cat.# 30218)

Bottom chromatogram:

Oven temp.: 45°C (hold 4 min.) to 110°C @ 19°C/min. (hold 5 min.)
 to 220°C @ 32°C/min. (hold 5 min.)
 Adjust dichlorodifluoromethane to 2.23 min. @ 45°C.
 MS Scan Range: 29-260amu, for 2Cl ethanol response
 compound concentrations, by mix: (in 5mL of RO water)
 Compounds at 100ppb (cat.# 30213, 30004, 30006, 30011, 30042)
 Alcohols at 1ppm (cat.# 30214) (see MS scan)
 vinyl acetate at 100ppb (cat.# 30216)
 8240 Nitrile Mix at 400ppb (cat.# 30215)
 8240 Mix 1A at 300ppb (cat.# 30217)
 8240 Mix 2A at 500ppb (cat.# 30218)



Slower Run
OPTIMIZED
Separations

1. dichlorodifluoromethane
2. chloromethane
3. vinyl chloride
4. bromomethane
5. chloroethane
6. trichlorofluoromethane
7. ethanol
8. 1,1-dichloroethene
9. carbon disulfide
10. iodomethane
11. allyl chloride
12. methylene chloride
13. acetone
14. *trans*-1,2-dichloroethene
15. 1,1-dichloroethane

16. acrylonitrile
17. allyl alcohol
18. vinyl acetate
19. bromochloromethane
20. chloroform
21. carbon tetrachloride
22. propargyl alcohol
23. 1,1,1-trichloroethane
24. 2-butanone
25. benzene
26. propionitrile
27. methacrylonitrile
28. 1,2-dichloroethane-d4
29. isobutyl alcohol
30. 1,2-dichloroethane

31. trichloroethene
32. 1,4-difluorobenzene
33. dibromomethane
34. 1,2-dichloropropane
35. bromodichloromethane
36. methyl methacrylate
37. 1,4-dioxane
38. 2-chloroethanol
39. 2-chloroethyl vinyl ether
40. *cis*-1,3-dichloropropene
41. toluene-d8
42. toluene
43. pyridine
44. 4-methyl-2-pentanone
45. tetrachloroethene

46. *trans*-1,3-dichloropropene
47. ethyl methacrylate
48. 1,1,2-trichloroethane
49. dibromochloromethane
50. 1,2-dibromoethane
51. 2-hexanone
52. 2-picoline
53. chlorobenzene-D5
54. ethylbenzene
55. chlorobenzene
56. 1,1,1,2-tetrachloroethane
57. *m*-xylene
58. *p*-xylene
59. *o*-xylene
60. styrene

61. bromoform
62. 4-bromo-1-fluorobenzene
63. *cis*-1,4-dichloro-2-butene
64. 1,1,2,2-tetrachloroethane
65. 1,2,3-trichloropropane
66. *trans*-1,4-dichloro-2-butene
67. pentachloroethane
68. 1,3-dichlorobenzene
69. 1,4-dichlorobenzene
70. benzyl chloride
71. malononitrile
72. 1,2-dichlorobenzene
73. 1,2-dibromo-3-chloropropane

EPA Method 624


- ◆ 37 aromatic and halogenated compounds in wastewater
- ◆ Previously used packed columns, but can be performed with many different capillary column dimensions
- ◆ Rtx-VMS 30m x 0.25mm x 1.4um gives separation in less than 15 minutes, starting at 40°C

EPA Method 624 Rtx®-VMS

1. chloromethane
2. vinyl chloride
3. bromomethane
4. chloroethane
5. trichlorofluoromethane
6. 1,1-dichloroethene
7. methylene chloride
8. *trans*-1,2-dichloroethene
9. 1,1-dichloroethane
10. bromochloromethane
11. chloroform
12. carbon tetrachloride
13. 1,1,1-trichloroethane
14. benzene
15. pentafluorobenzene
16. 1,2-dichloroethane
17. fluorobenzene
18. trichloroethene
19. 1,2-dichloropropane
20. bromodichloromethane
21. 2-chloroethyl vinyl ether
22. *cis*-1,3-dichloropropene
23. toluene
24. tetrachloroethene
25. *trans*-1,3-dichloropropene
26. 2-bromo-1-chloropropane
27. 1,1,2-trichloroethane
28. dibromochloromethane
29. ethylbenzene
30. chlorobenzene
31. bromoform
32. 4-bromofluorobenzene
33. 1,4-dichlorobutane
34. 1,1,2,2-tetrachloroethane
35. 1,3-dichlorobenzene
36. 1,4-dichlorobenzene
37. 1,2-dichlorobenzene

30m, 0.25mm ID, 1.40 μ m Rtx®-VMS (cat#19915)
Conc.: 20 ppb in 5mL of RO water
Concentrator: Tekmar LSC-3000 Purge and Trap
Trap: Vocabr 3000 (type K)
Purge: 11 min. @ 40mL/min. @ ambient temperature.
Dry purge: 1 min. @ 40mL/min. (MCS bypassed using Silcosteel® tubing)

Desorb preheat: 245°C
Desorb: 250°C for 2 min., Flow 10mL/min.
Bake: 260°C for 8 min.
GC Interface: 1:10 split at injection port. 1mm ID sleeve.
GC: HP 6890
Oven temp.: 40°C (hold 4 min.) to 95°C @ 24°C/min. (hold 3 min.), to 210°C @ 40°C/min. (hold 6 min.)
Carrier gas: helium @ ~1mL/min. constant flow
Adjust dichlorodifluoromethane to a retention time of 2.54 min. @ 40°C
Detector: HP 5973 MSD
Scan range: 25-300amu




EPA Method 524.2 revision IV

- ◆ Monitors 84 compounds in drinking water
- ◆ Rtx-VMS overcomes closely eluting pairs found on other phases (ex. “624”) such as 1,1-dichloro-2-propanone/4-methyl-2-pentanone and methyl acrylate/propionitrile
- ◆ Rtx-VMS 30m x 0.25mm x 1.4um performs the analysis in less than 11 minutes
- ◆ Only difficult pair on the Rtx-VMS is 2-nitropropane/ 1,1-dichloro-2-propanone which share ion 43

EPA Method 524.2, Revision 4 Rtx®-VMS

Innovations at Work
Excellent Gas Resolution!



30m, 0.25mm ID, 1.4µm Rtx®-VMS (cat.# 19915)
 Carrier gas: helium @ ~1.3mL/min. constant flow
 Adjust dichlorodifluoromethane to a retention time of 2.29 min. @ 45°
 Concentrator: Tekmar LSC-3000 Purge and Trap
 Oven temp.: 45°C (hold 2 min.) to 85°C @ 14°C/min. to 210°C @ 40°C/min. (hold 4 min.)
 GC: HP 6890 Series II
 Trap: Vocabr 3000
 Purge: 11 min. @ 40mL/min.
 Dry purge: 1 min. @ 40mL/min. (MCS bypassed)
 Desorb preheat: 245°C
 Desorb: 250°C for 2 min.
 Bake: 260°C for 8 min.
 Interface: 1:10 split in port
 Transfer line: 5m, 0.32mm ID Siltek™ tubing (cat.# 10027)
 Detector: HP 5973 MSD
 Scan range: 35-300amu

Standards:
 20ppb in 5mL of RO water (unless otherwise noted); ketones at 40ppb.
 502.2 Cal Mix #1 (cat.# 30042)
 502.2 Cal2000 MegaMix™ (cat.# 30431)
 524 Cal Mix 7A & 7B (cat.# 30202)
 524 Cal Mix #8 (cat.# 30203)
 524 IS/SS Mix (cat.# 30201)

1. dichlorodifluoromethane	23. carbon tetrachloride	45. <i>trans</i> -1,3-dichloropropene	67. 1,2,3-trichloropropane
2. chloromethane	24. tetrahydrofuran (40ppb)	46. ethyl methacrylate	68. <i>trans</i> -1,4-dichloro-2-butene
3. vinyl chloride	25. 1,1,1-trichloroethane	47. 1,1,2-trichloroethane	69. 4-chlorotoluene
4. bromomethane	26. 2-butanone	48. dibromochloromethane	70. <i>tert</i> -butylbenzene
5. chloroethane	27. 1,1-dichloropropene	49. 1,3-dichloropropane	71. 1,2,4-trimethylbenzene
6. trichlorofluoromethane	28. 1-chlorobutane	50. 1,2-dibromoethane	72. pentachloroethane
7. diethyl ether	29. benzene	51. 2-hexanone	73. <i>sec</i> -butylbenzene
8. 1,1-dichloroethane	30. propionitrile	52. ethylbenzene	74. <i>p</i> -isopropyltoluene
9. carbon disulfide (40ppb)	31. 1,2-dichloroethane	53. chlorobenzene	75. 1,3-dichlorobenzene
10. iodomethane (40ppb)	32. fluorobenzene	54. 1,1,1,2-tetrachloroethane	76. 1,4-dichlorobenzene
11. allyl chloride	33. trichloroethene	55. <i>m</i> -xylene	77. <i>n</i> -butylbenzene
12. methylene chloride	34. dibromomethane	56. <i>p</i> -xylene	78. hexachloroethane
13. acetone	35. 1,2-dichloropropane	57. <i>o</i> -xylene	79. 1,2-dichlorobenzene-d4
14. <i>trans</i> -1,2-dichloroethane	36. bromodichloromethane	58. styrene	80. 1,2-dichlorobenzene
15. methyl <i>tert</i> -butyl ether	37. methyl methacrylate	59. bromoform	81. 1,2-dibromo-3-chloropropane
16. 1,1-dichloroethane	38. <i>cis</i> -1,3-dichloropropene	60. isopropylbenzene	82. nitrobenzene
17. acrylonitrile	39. toluene	61. 4-bromofluorobenzene	83. hexachlorobutadiene
18. <i>cis</i> -1,2-dichloroethane	40. chloroacetonitrile	62. <i>n</i> -propylbenzene	84. 1,2,4-trichlorobenzene
19. 2,2-dichloropropane	41. 2-nitropropane*	63. bromobenzene	85. naphthalene
20. bromochloromethane	42. 1,1-dichloropropanone*	64. 1,1,2,2-tetrachloroethane	86. 1,2,3-trichlorobenzene
21. chloroform	43. 4-methyl-2-pentanone	65. 1,3,5-trimethylbenzene	
22. methyl acrylate	44. tetrachloroethene	66. 2-chlorotoluene	

*These peaks (41 and 42) share a quantitation ion (43)

Conclusions

- ◆ Computer modeling can be used to optimize phase development, column dimensions and run conditions
- ◆ Column dimensions to interface with the GC/MS is a matter of customer preference
- ◆ Analysis of 100 + compounds typically use longer columns for better resolution and higher starting temperatures
- ◆ Rtx-VMS has excellent selectivity and rapid cycle time for all commonly used EPA GC/MS volatile methods

For More Information...

Author e-mail: cenglish@restekcorp.com

Restek Corporation
110 Benner Circle
Bellefonte, PA 16823
(800) 356-1688