

DESIGN OF NEW LOW-BLEED GAS CHROMATOGRAPHIC STATIONARY PHASES AND CAPILLARY COLUMNS

Frank L. Dorman, Gary B. Stidsen, Chris M.
English, Rick Morehead, Jack Cochran, Eric
J. Reiner, and Karen MacPherson

Restek Corporation
www.restekcorp.com



When is a Low-Bleed Column Important?

- Low concentration GC-MS analyses
 - Full scan
 - SIM
 - High-resolution
- When detector is especially sensitive to bleed
 - ECD
- When bleed directly interferes with quantitation
 - Bleed and analyte share the same ions (rare)

How Do You Make a Low-Bleed Column ?

- Thinner stationary phase thickness
 - Not really purchasing what you think...
- Higher conditioning temperatures
 - Can lead to higher reactivity and phase loss
- Improve deactivation chemistry
 - Phase has more “affinity” for tubing
- Move polymer functionality from side chain to backbone

Considerations for Backbone-modified Polysiloxane Stationary Phases:

- By using a “low-bleed” column, you do not want to give up selectivity.
 - Many “low bleed” columns are dimethyl or phenyl/methyl-type
- Computer Assisted Stationary Phase Design (CASPD) utilized to optimize selectivity of backbone-stabilized polymers for specific applications
 - *Anal. Chem.* **74**(9), 2133-2138 2002.
 - *LC*GC* **18**(9), 928, 2000.
 - *American Laboratory*, **31**(6), 20-26, 1999.

Applications

Rtx-PCB Capillary GC Column

- Proprietary polysiloxane designed for PCB congener separation for GC-MS analysis
- 380 maximum operating temperature in standard high-temperature fused silica tubing
- Standard column dimensions, standard film thickness
- Generally increased retention for aromatics over “5”, “5Si1 MS” and “XLB” type phases

Analysis of PCB Congeners – GC-MS

- European congener list:
 - BZ # 28,52,101,118,153,138 and 180 resolved by CI level or chromatographically from all others
- 1668
 - All 13 congeners resolved
- McFarland and Clarke congeners:
 - 33 of 36 resolved (1 more than best previous column – XLB)
- Aroclor congeners out of 139:
 - 127 resolved (3 more than XLB)
 - 9 pairs are unresolved
- 209 Congeners:
 - 168 resolved (same number but different congeners as XLB)

GC-ToF Analysis

- Pyrethroids analysis
 - Late-eluting pyrethroids can be difficult to keep well resolved on standard columns due to limited maximum operating temperature.
 - Peak width is narrower if the compounds are eluted during the temperature program, not on an isothermal hold
 - Preference is for column with high maximum operating temperature
 - Selectivity must be appropriate to achieve separation between the target compounds which can have similar mass spectra


GC Conditions

- Splitless injection
 - One microliter at 250°C
 - 60 sec valve time
 - 30 m x 0.25 mm x 0.25 μm Restek Rtx-PCB column
 - Constant flow helium, 1 mL/minute
- GC oven program
 - 60°C (1 min), 30°/min to 120°, 15°/min to 360° (1 min)
 - Total run time: 20 min

MS Conditions


- Source temperature: 225°C
- Electron ionization: 70 eV
- Stored mass range: 45 to 550 u
- Acquisition rate: 10 spectra/sec

Pyrethroids on Rtx-PCB




Re

Pyrethroids on Rtx-PCB



Separation of Cyfluthrins and Cypermethrins on Rtx-PCB




Pyrethroids on Rtx-PCB											
	RT	RT		Elution					Ions		
Pesticide	sec	min	RRT	Temp °C	MW	Formula	CAS#	1	2	3	
Pentachloronitrobenzene	680.5	11.34	1.000	245	293	C6Cl5NO2	82-68-8	237	249	295	
Bifenthrin	916.0	15.27	1.346	304	410	C22H22O2ClF3	82657-04-3	165	166	181	
Fenpropathrin	933.5	15.56	1.372	308	349	C22H23NO3	39515-41-8	97	181	265	
lambda-Cyhalothrin	961.0	16.02	1.412	315			91465-08-6	181	197	208	
cis-Permethrin	1002.2	16.70	1.473	326	390	C21H20Cl2O3	54774-45-7	127	163	183	
trans-Permethrin	1007.9	16.80	1.481	327	390	C21H20Cl2O3	61949-77-7	127	163	183	
Cyfluthrin	1025.5	17.09	1.507	331	433	C22H18Cl2FNO3	68359-37-5	163	206	226	
Cyfluthrin	1029.4	17.16	1.513	332	433	C22H18Cl2FNO3	68359-37-5	163	206	226	
Cyfluthrin	1032.5	17.21	1.517	333	433	C22H18Cl2FNO3	68359-37-5	163	206	226	
Cyfluthrin	1033.9	17.23	1.519	333	433	C22H18Cl2FNO3	68359-37-5	163	206	226	
Cypermethrin	1040.5	17.34	1.529	335	415	C22H19Cl2NO3	52315-07-8	163	181	209	
Cypermethrin	1044.8	17.41	1.535	336	415	C22H19Cl2NO3	52315-07-8	163	181	209	
Cypermethrin	1047.5	17.46	1.539	337	415	C22H19Cl2NO3	52315-07-8	163	181	209	
Cypermethrin	1048.9	17.48	1.541	337	415	C22H19Cl2NO3	52315-07-8	163	181	209	
Esfenvalerate	1085.6	18.09	1.595	346	419	C25H22ClNO3	66230-04-4	167	225	419	
Fenvalerate	1092.5	18.21	1.605	348	419	C25H22ClNO3	51630-58-1	167	225	419	

Rtx-500 Capillary GC Column


- Carborane-stabilized stationary phase
 - Maximum temperature of 380 C in “standard high-temperature” tubing
 - Maximum temperature of 440 C in passivated metal columns (Mxt-500)
- Extremely low bleed levels
 - Surpasses phenyl/methyl phases, and silphenylene stationary phases
- Common dimensions available

Carborane-polysiloxane Stationary Phase

Dimethyl Polysiloxane



Carborane Dimethyl Polysiloxane




Analysis of PCB Congeners by GC-MS

- Toxic PCB's by USEPA method 1668 commonly run on “octyl” phases
 - Very high background from bleed decreases sensitivity
 - Phase loss can cause retention order changes
 - “5” phases can have coelution issues
 - 118/123
 - 156/157
 - Rtx-500 Column can improve upon these...

Toxic PCB Congeners


dlp-cs5 on htc
Feb06_02samp1



Analysis of Siloxanes in the Environment

- Rare case, but becoming important
- Since most all GC phases are siloxanes, bleed ions directly interfere with quantitation ions for GC-MS analyses
- Rtx-500 Carborane-based capillary column does not have any siloxane character in its bleed spectrum
- Rtx-500 also useful for PBDE flame retardant analyses (Pittcon 2003 Workshop)

Polybrominated Diphenylethers



Br 1-10


- Named similarly to PCB congeners (BDE 1-209)
J. High Resolut Chromatogr **15**:260
- Human exposure via the food chain is 0.2 – 0.7 mg/day
Organohalogen Compounds **35**:411
- Listed as Endocrine Disruptors
Environ Health Perspect **101**:378

Wellington Laboratories BDE Mix-C

4-bromodiphenyl ether (3)	250	141,115
2,4-Dibromodiphenyl ether (7)	168	328,139
4,4'-dibromodiphenyl ether (15)	328	168,221
2,2',4-tribromodiphenyl ether (17)	248	408,406
2,4,4'-tribromodiphenyl ether (28)	406	248,246
2,2',4,5'-tetrabromodiphenyl ether (49)	326	486,328,324
2,3',4',6-tetrabromodiphenyl ether (71)	326	486,324,328
2,2',4,4'-tetrabromodiphenyl ether (47)	326	486,328,324
2,3',4,4'-tetrabromodiphenyl ether (66)	326	486,328,324
3,3',4,4'-tetrabromodiphenyl ether (77)	326	486,328,324
2,2',4,4',6-pentabromodiphenyl ether (100)	406	564,566
2,3',4,4',6-pentabromodiphenyl ether (119)	404	406,564
2,2',4,4',5-pentabromodiphenyl ether (99)	406	564,566
2,2',3,4,4'-pentabromodiphenyl ether (85)	406	564,566
3,3',4,4',5-pentabromodiphenyl ether (126)	566	564,568,406
2,2',4,4',5,6'-hexabromodiphenyl ether (154)	484	644,486
2,2',4,4',5,5'-hexabromodiphenyl ether (153)	644	484,486,482
2,2',3,4,4',5'-hexabromodiphenyl ether (138)	642	484
2,2',3,4,4',5',6-heptabromodiphenyl ether (183)	722	564
decabromodiphenyl ether (209)	956	

Baseline separation of Tri, Tetra, Penta, Hexa, Hepta, and DecaBDE's congeners in 44 minutes!

bde209 13C
Feb04_02bde2



Feb05_02bde2



Feb05_02bde2




Feb05_02bde2




Feb05_02bde2



Feb05_02bde2




Feb05_02bde2



Time

bde209 13C
Feb04_02bde2

Voltage SIR 7 Channels EI+
969.2102
1.18e6




13c-Deca-BDE in 43.4 minutes


GC-ECD Analysis

- Electron capture detector is common to many laboratories
- Compounds have excellent response by ECD
- Higher flow rates may allow for more rapid separation using larger-diameter columns
- Instrumentation less expensive than HRMS
- Instrumentation is also field portable


Wellington Laboratories BDE Mix-C




0.53 mm i.d. column partially resolves BDE 154 and 126 under these conditions



Commercially Available PBDE Mixes



Commercially Available PBDE Mixes



Conditions

Column: Rtx-500 30m x 0.53mm x 0.15 Ser# 233548
Flows: Hydrogen 7.69ml/min (66.7cm/sec) @ 100°C
Dead Time: MeCl₂ Headspace 0.76min @ 100°C
Oven: 100(1) 15/260(5) 15/380 (15) ~ 40 min runtime
Standards: Wellington Laboratory BDE-Mix C
1ul injection 30pg on column.
Injection: Drilled Press-Tight Uniliner.

Rtx-Dioxin Capillary Column

- Proprietary polysiloxane designed for replacement of “5”-type columns as primary, or high-cyano secondary columns for toxic dioxin and furan analysis by GC-HRMS
- 380 maximum operating temperature in standard high-temperature fused silica tubing

Dioxin and Furan Analysis


- Dual column method
 - Usually 5% diphenyl column and a high-cyano column (eg Rtx-225)
 - Cyano columns have poorer lifetimes and lower maximum operating temperatures
 - 5% diphenyl phases do not have the selectivity to accurately quantitate most samples
 - USEPA requires 2,3,7,8-tcdf to be confirmed on a X-225
- Desirable to have both columns in the same oven, and to improve the separation of the “5”

Feb06_0

TCDD Resolution Mix

Voltage SIR 14 Channels EI+

321.8936
3.53e5




JUL26_QC-S003

TCDD Resolution Mix

2: SIR of 20 Channels EI+


321.894
2.35e5



OCT24_QCS001

SIR of 12 Channels EI+
303.902
2.88e6


db225 TCDF RESOLUTION MIX
30M x .25mm x .25um



feb18_37

2: Voltage SIR 20 Channels EI+
303.9016
8.24e6

RTX-DIOXIN TCDF RESOLUTION MIX
40M x .18mm x .1um



Time

Fly Ash Round Robin Samples

>110 laboratories participating

	DB-5	DB-225	RTX-DIOXIN	MEDIAN	MEAN
Ash A	250	21	30	28	32
Ash B	2100	300	378	390	390
Ash C	170	19	28	27	32
All results reported as pg/g					

Median and Mean agreement gives good confidence in “true” value

Rtx-Dioxin Conditions


Micromass Altima High Resolution GC-MS			
Rtx-Dioxin 40m x .18 x .1			
Initial Temp 130 C			
	Time	Rate C/min	Temp
	0	52	200
	10.2	2.9	235
	10	6.9	300
	24		
Constant Pressure of 1.2 mL/min			
Injector Temp = 270C			

JUN18 Flyash extract
DB5
60M x .25mm x .25um

tcdf ion window


Up to 6 congeners coeluting

2: SIR of 20 Channels EI+
303.902
2.28e6




feb18 Flyash extract
RTX-DIOXIN
40M x .18mm x .10um

2: Voltage SIR 20 Channels EI+
303.9016
7.01e5




OCT24 Flyash extract
DB225
30M x .25mm x .25um

SIR of 12 Channels EI+
303.902
2.31e4



feb18_33

Flyash extract
RTX-DIOXIN
40M x .18mm x .10um



2: Voltage SIR 20 Channels EI+
303.9016
7.01e5


feb18_33



2: Voltage SIR 20 Channels EI+
317.9389
1.33e6


OCT24_CONFIRM2S019 Sm (Mn, 1x1)

Flyash extract
DB225
30M x .25mm x .25um



SIR of 12 Channels EI+
303.902
2.31e4

OCT24_CONFIRM2S019



SIR of 12 Channels EI+
317.939
4.54e5

Time

Rtx-Dioxin Column

- Replaces a 5% diphenyl column for high-resolution dioxin and furan analyses
 - Improves sensitivity due to lower bleed levels
 - Improves separation of many congeners
- May replace high-cyano columns as confirmation column to the 5% diphenyl
 - All but 2 congeners were within 10% of the fly ash median values
 - These 2 can be quantitated using the 5% diphenyl column
 - May eliminate the need for –2330, -2331, -225 confirmation columns
 - Final characterization currently underway

Low-Bleed Columns Summary:

- **Rtx-XLB**
 - Low bleed column for pesticides and PCB Congeners, similar to DB-XLB
- **Rtx-PCB**
 - Low bleed column for pesticides and PCB Congeners, resolves a few more congeners than XLB phases
 - Unbiased GC-MS results for European PCB congener method
- **Rtx-500**
 - Carborane-based column, exceptionally high thermal stability
- **Rtx-Dioxin**
 - Improvement over the 5% diphenyl columns, and possible replacement for high-cyano confirmation phases

Acknowledgements

- Reference materials courtesy of Wellington Laboratories – Guelph, Ontario
- Dr. Eric Reiner and Karen MacPherson of the Ontario Ministry of the Environment – Etobicoke, Ontario provided the HRMS analyses, and guidance.
- Jack Cochran of LECO Corporation – Pesticide and PCB analyses