

# Improving Method Performance through Fast LC

C. Vernon Bartlett, B.A.; Terrence S. Reid,  
M.S.; Rebecca E. Wittrig, Ph.D

Restek Corporation  
[www.restekcorp.com](http://www.restekcorp.com)



## Fast LC Methods

- Use of columns that can operate at high flow rates with reduced pressures (increased  $k'$ )
- Use of a reduced particle size substrate ( $<3\mu\text{m}$ ) to force higher theoretical plate counts (increased  $N$ )
- Use of shortened standard packed columns with optimized and possibly unique stationary phases (increased  $\alpha$ )

# Fast LC Technique

- Highly selective stationary phase is desired to maximize alpha values.
- Elution of components is typically accomplished through the use of gradients to reduce retention of highly retained components.
- Simple resolution of methylene substitutions /additions may be accomplished isocratically.
- Good screening technique for unknowns.

## Fast LC Technique – Advantages


- Faster re-equilibration (when using gradients)
- Sensitivity improvements.
- Older qualitative techniques can be adapted to a highly automated quantitative technique.
- Allows potentially high increases in sample throughput.
- Great technique when performed by LC-MS
- Shorter analysis times reduce solvent consumption and waste.

## Fast LC Technique – Disadvantages

- Critical separations are more sensitive to extra-column volume (as post column reactors.)
- Extremely selective stationary phase must be used to maximize selectivity – especially for structural isomers.
- May not be well suited to normal phase or ion pairing separations (with gradients).

# Principles and Theory of HPLC

## General Resolution Equation

$$R = \frac{1}{4} \left( \frac{a-1}{a} \right) \sqrt{N} \left( \frac{k'}{k'+1} \right)$$
The equation is displayed with three colored underlines: a green underline under the fraction (a-1/a), a yellow underline under the square root of N, and a pink underline under the fraction (k' / (k'+1)). Three arrows originate from these underlines: a green arrow points to the Selectivity box, a yellow arrow points to the Efficiency box, and a pink arrow points to the Retention box.

### Selectivity

- stationary phase
- mobile phase composition
- additives


### Efficiency

- particle size
- column length


### Retention

- chain length
- mobile phase strength

# Column Length vs Theoretical Plates



# Plate Count vs Particle Size and Flow Rate



4.6 mm Column




# Principles and Theory of HPLC

$$\frac{K'}{1+K'}$$


## Resolution Depends on $K'$

$K'$ Value	$K'$ Term	$K'$ Resolution ?
0	0	0
1	1/2	.50
2	2/3	.67
3	3/4	.75
10	10/11	.91
20	20/21	.95

$K'$  Ideal Range ?






# Resolution Increase vs Alpha



# Fast LC improvement of USP TLC and HPLC Method


## Digitalis Extracts and Derivatives




Digitoxin:  $X_1=X_2=H$

Digoxin:  $X_1=H; X_2=OH$

Gitoxin:  $X_1=OH; X_2=H$



# Fast LC Separation of Digitalis Derivatives Experimental Fluorophase (3 minutes)



# Advantages of Digitalis Fast LC over Current Methods

- Improved automation and analysis throughput.
- Reduction of analysis time – previously a 30 cm length C18 column was required for resolution.
- Ability to analyze all materials by HPLC.
- Ability to precisely quantitate materials vs TLC.
- Highly selective stationary phase.
- Technique can also be applied to purification and analysis of digoxin labeled materials for biological activity.
- Perfect technique for use in high speed cleaning validations.

# Fast LC Analysis of Carbamate Insecticides




Figure 3-Structures of Commonly Analyzed Carbamates

# Carbamate Analysis using Standard HPLC Methodology (About 40 minutes)

## Carbamate Pesticides on Pinnacle Carbamate

Applications Note: LC\_0192


### Peak List:

1. aldicarb sulfoxide
2. aldicarb sulfone
3. oxamyl
4. methomyl
5. 3-hydroxycarbofuran
6. aldicarb
7. propoxur
8. carbofuran
9. carbaryl
10. 1-naphthol
11. methiocarb

### Column: Pinnacle Carbamate

Catalog#: 9173575  
Dimensions: 250 x4.6mm  
Particle size: 5µm  
Pore size: 120Å  
Solvent A: Water  
Solvent B: Methanol  
Gradient Program:  
20-68% B; 0-20 min.  
68-79% B; 0-34 min.  
Flow Rate: 1.0mL/min  
Detection: Fluorescence  
Ex. 339 Em. 445nm

*Courtesy of Analytical Laboratories, Inc.  
Boise, Idaho*



# Fast LC Separation of Carbamates (About 13 minutes)

## Fast LC Separation of 11 Carbamates on Ultra Carbamate

### Peak List:

1. aldicarb sulfone
2. aldicarb sulfoxide
3. oxamyl
4. methomyl
5. 3-hydroxycarbofuran
6. aldicarb
7. propoxur
8. carbofuran
9. carbaryl
10. methiocarb
11. 4-bromo-3,5-dimethylcarbamate

### Sample:

Inj.: 5 $\mu$ L  
Conc.: 50 $\mu$ g/mL  
Solvent: methanol

### Restek standards:

Catalog# 32274 and 32273 mixed 50:50


### Column: Ultra Carbamate

Catalog #: 9177355  
Dimensions: 50 x 4.6mm  
Particle size: 3 $\mu$ m  
Pore size: 100 $\text{\AA}$

### Conditions:

Mobile Phase: A: 90:10 water:methanol  
B: 90:10 methanol:acetonitrile  
Time (min): %B  
0 10  
10 90

Flow: 1.5mL/min  
Temp.: 27 $^{\circ}$ C  
Det.: UV @ 220nm





# Fast LC Analysis of Carbamates with MS Detection

Table II. Experimental conditions for the LC/MS analysis of carbamate compounds.

## HPLC Conditions

Column: Ultra Carbamate, 100 mm x 4.6 mm, 3  $\mu$ m  
 Mobile Phase A: 90% water:10% methanol with 10mM ammonium formate  
 Mobile Phase B: 10% acetonitrile:90% methanol with 10mM ammonium formate  
 Gradient: 90%A:10%B to 10%A:90%B from 0-15 minutes  
 Inj. Volume: 10  $\mu$ L  
 Flow Rate: 0.75 mL/min to UV detector, 0.25 mL/min to MSD

## MSD Conditions

	Compound	Ion	Cone V
Detector:	1	223.3	25V
Mode:	2	207.3	18V
Capillary V:	3	237.2*	10V
Extractor:	4	163.2	15V
Ion Energy:	5	238.3	15V
Multiplier:	6	191.2	8V
Source Temp:	7	210.2	18V
Desolv. Temp:	8	222.3	22V
Gas Flow:	9	202.2	18V
	10	226.3	19V

\*Ammonium adduct (all other are [M+H]<sup>+</sup> ions)




Figure 6

# Conventional Vanilla Flavoring Analysis (About 20 minutes with gradient)

## Vanillin and Ethyl Vanillin on Ultra C8

Applications Note: LC\_0148

Peak List:	Conc. (mg/mL)
1. vanillin	0.12
2. ethyl vanillin	0.04


Sample:  
Inj.: 10mL  
Solvent: 40% ethanol

Column: Ultra C8  
Catalog#: 9103565  
Dimensions: 150x4.6mm  
Particle size: 5µm  
Pore size: 100Å

Conditions:  
Mobile phases: A: 1.2% acetic acid  
B: methanol  
Gradient:

Minutes	%B
0.0	20.0
5.0	20.0
15.0	40.0
20.0	40.0
25.0	20.0

Flow: 1.0mL/min.  
Temp.: 28°C  
Det.: UV @ 254nm



## Vanillin on Ultra C8

Application Note: LC\_0149

Column: Ultra C8  
Catalog#: 9103565  
Dimensions: 150x4.6mm  
Particle size: 5µm  
Pore size: 100Å


Peak List:  
1. vanillin

Sample:  
Inj.: 10mL  
Conc.: 5% solution of vanilla extract  
Solvent: 40% ethanol

Conditions:  
Mobile phases: A: 1.2% acetic acid  
B: methanol  
Gradient:

Minutes	%B
0.0	20.0
5.0	20.0
15.0	40.0
20.0	40.0
25.0	20.0

Flow: 1.0mL/min.  
Temp.: 28°C  
Det.: UV @ 254nm



# Highly Selective Fast LC Vanillin Analysis (Less than 5 minutes and isocratic)





Figure 9

## Vanillin Fast LC Analysis Advantages

- Takes advantage of improved methylene selectivity by the use of C18 phase versus a C8
- Analysis can be conducted in 3 to 5 minutes versus older 15 minute methods.
- High but not excessive selectivity of C18 phases toward simple methylene substitutions allow the use of an isocratic mobile phase.
- Reduction of column length also reduces the time needed for more hydrophobic analytes to elute from system.

# EPA Method 610 – PAHs on Pinnacle II PAH



# PAHs on Pinnacle II PAH- Fast LC


Column: Pinnacle II PAH  
Catalog #: 5219555  
Dimensions: 50mm x 4.6mm  
Particle Size: 5  $\mu\text{m}$   
Pore Size: 110  $\text{\AA}$

## Conditions:


Mobile Phase A: purified water  
Mobile Phase B: acetonitrile

Time	%B
0	50
1.5	60
4.0	100
5.5	100
5.6	50

Flow: 3.0 mL/min  
Temp.: ambient  
Det.: UV @ 254nm



# Fast LC Analysis of Carboxylic Acids on Ultra Aqueous C18



# Fast LC Analysis of Aromatic Amino Acids on Ultra Aqueous C18

## Column:

Catalog #: 9178535  
Dimensions: 30 mm x 4.0 mm  
Particle Size: 5  $\mu$ m  
Pore Size: 100 Å

## Conditions:

Mobile Phase: A: 50mM Potassium Phosphate, pH 2.5

B: Acetonitrile

Time:	%B
0.0	0%
1.0	40%
1.2	0
5.0	0

Flow: 1.0 mL/min  
Temp.: 50 °C  
Det.: UV @ 254 nm


3

Peak List:	Peak Retention Time	Conc. (mg/ml)
1. Tyrosine	1.286	1.25
2. Phenylalanine	2.197	0.55
3. Tryptophan	2.475	0.035

## Sample:

Inj.: 10  $\mu$ L

Sample Diluent: 50mM Potassium Phosphate, pH 2.5





## Conclusions

- Highly selective and sometimes unique stationary phases make Fast LC a reality using conventional hardware and techniques.
- Fast LC is viable, precise quantitative alternative for analysis previously performed by Thin Layer Chromatography.
- Fast LC can improve method sensitivity, reduce solvent waste, and enhance laboratory throughput.

# Conclusions

- The proper stationary phase selection can change a gradient system to a faster isocratic system.
- Sometimes only a simple phase change is needed to convert a method to a Fast LC separation.
- Selectivity is still the greatest factor in separation optimization.

# Acknowledgements

- The authors would like to thank Randy Romesberg, Larry Peters, and Rahul Patil of Restek Corporation for their participation in creating columns and hardware.