

# Food & Flavor Analysis: Sample Extraction and Introduction

**Becky Wittrig, Ph.D.**  
RESTEK CORPORATION



# Food & Flavor Analysis: Sample Extraction and Introduction

- I. Extraction Techniques
  - A. Solvent Extraction
  - B. Solid Phase Extraction (SPE)
  - C. SPME
  - D. SBSE
- II. Headspace Sampling
  - A. Static Headspace
  - B. Dynamic Headspace
  - C. Thermal Desorption

# Sampling Challenges

- Very Low Concentrations
  - ppm to ppt
- Matrix Issues
  - Volatiles can be intracellular
  - Disruption might be necessary
- Very Complex Volatile Profiles
  - e.g. coffee has ~800 components
- Wide Range of Volatilities
- Air + Heat Instability

# Solvent Extraction Techniques

- Solid-Liquid
  - Soxhlet extraction
- Liquid-Liquid
  - Separatory funnel
  - Liquid-liquid extractors
- Solvents
  - Methylene chloride – good all-purpose
  - Diethyl ether, pentane, hydrocarbons, Freons
- Concentration Techniques
  - Dry over sodium or magnesium sulfate
  - Concentrate on steam bath


# Distillation Techniques

- Prior to liquid-liquid extraction
  - Distillation can remove nonvolatile compounds
  - Can be done under vacuum – less heat required
- *Supercritical CO<sub>2</sub> extractions are gaining in popularity – “cleaner” final sample*

# Solid Phase Extraction (SPE)

- “Trap” specific compounds
  - Compound of interest
  - Interferences
- Wide range of materials available
  - Non-polar – C18
  - Polar – carbon materials, fluorisil
- Standardized procedures
  - i.e. EPA test methods
- Validation is very important
  - Recoveries, breakthrough levels

# HMF in Grape Juice: No SPE Clean-up




# SPE Clean-up Procedure

1. Conditioning
  - a. Apply 3 mL methanol
  - b. Apply 3 mL deionized water
2. Sample Application
  - a. Apply 4 mL sample to moist SPE tube, gravity feed
3. Wash
  - a. Pull remaining sample through tube, using vacuum
  - b. Apply 3 mL water
  - c. Remove excess water from bed under vacuum
4. Elution
  - a. Apply 2 mL elution solvent, gravity feed, dilute to volume




# HMF in Grape Juice

C18 SPE, HMF Eluted with 20% Methanol



# HMF in Grape Juice

C18 SPE, HMF Eluted with 50% Methanol



# Solid Phase Microextraction (SPME)


## ■ Solvent-less Extraction Technique

- ◆ Partitioning of organic components between an aqueous or a vapor phase and a thin polymeric film
- ◆ Independent of the matrix
- ◆ Applicable to liquids, solids, and gases

## ■ Factors Affecting the Extraction

- ◆ Organic/water partition coefficient
- ◆ Polarity & volatility
- ◆ Volume of sample/headspace, pH, temperature, etc.

# Principle of SPME



## Fiber types:

PDMS

Carboxen

Carbowax

Divinylbenzene

Combinations


*-expose fiber to liquid, headspace*  
*-thermally desorb or solvent extract*

# SBSE: The Gerstel Twister™

- 1.5cm long magnetic stir bar sealed in glass
- High capacity PDMS phase on glass
- Extremely rugged
- Preconditioned for low background
- Stirs and extracts in one step
- Splitless desorption of stir bar gives low detection limits



# Stir Bar Sorptive Extraction (SBSE)




- Construction:
  - Magnetic core material, sealed in glass (1.5 cm)
  - Glass deactivated and persilylated
  - PDMS phase (0.5 mm thick)
  - Solvent washed and thermally conditioned
- Amount of PDMS on stir bar:
  - Minimum, 1 cm bar: 24  $\mu$ l
  - Maximum, 2 cm bar: 126  $\mu$ l


# SBSE:

Lower Detection Limits due to Higher Analyte Recovery

- Typical SPME phase:  $<0.5\mu\text{l}$
- Typical Twister phase:  $24\text{-}126\mu\text{l}$
- Greater SBSE capacity gives higher recovery of analytes relative to SPME as polarity increases



# SBSE: How to use



- Add stir bar to vial
- Stir 1hr-overnight
- Remove stir bar with forceps
- Rinse briefly in distilled water
- Dry with lint-free tissue
- Transfer bar into thermal desorption tube
- Thermally extract
- Refocus analytes in cold inlet




# Stir Bar Sorptive Extraction

## Features/Benefits

- ▽ Immobilized liquid PDMS phase retains negligible water
- ▽ Linear sorption isotherms with minimal displacement effects
- ▽ PDMS phase selectivity enhances extraction of nonpolar analytes
- ▽ Selectivity eliminates polar matrix interference
- ▽ Extractions from lipophilic matrices (e.g. dairy, soaps) possible
- ▽ High recovery and splitless sample introduction provide extremely low detection limits
- ▽ Predictable analyte recovery
- ▽ Sample extraction done in parallel with minimal labor


# Stir Bar Sorptive Extraction

## Eliminates Matrix Interferences




- Polar matrix components (water, acetic acid, ethanol, PG, glycerol) do not partition into PDMS.
- Volatile components (ethanol, methanol, acetone, etc.) can be further reduced by allowing brief evaporation time.


# Overlays of 1 yr old Bourbon, Triplicate SBSE



# Overlays of 1yr, 2yr, 4yr, 6yr, 8yr and 10yr Bourbon



# SBSE of Baking Soda Impurities



# Flavor Profiling

## Headspace Sampling Techniques

- Advantages
  - “Cleaner” samples
  - Minimal sample preparation
- Types
  - Static headspace
  - Dynamic headspace
  - Thermal desorption

# Flavor Profiling


## Static Gas Extraction

- **Advantages**


- Excellent screening tool
- Inexpensive
- Minimal sample carryover
- Easy to perform

- **Disadvantages**

- Less sensitive
- Involves preparing calibration in sample matrix



## Flavor Profiling: Static Gas Extraction



- Constant ratio between liquid & gas phases at equilibrium



Where  $K = (C_{\text{gas}} / C_{\text{liquid}})$

$K$  = distribution coefficient

- Large  $K$  values favor the gas state
- $C_{\text{gas}}$  is directly proportional to peak area




# Flavor Profiling: Static Gas Extraction



Volatiles from 2 different batches of chewing gum. The headspace was sampled after heating to 60°C.

# Flavor Profiling: Static Gas Extraction

## Residual solvents in decaffeinated lemon tea




# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Inlet System

- Concentrates volatiles
- Dynamic extraction of solids and liquids
- Adsorbent trap
- Desorb (10-80 mL/min)
- Narrow bore column (split flow)
- Part of GC system



# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Sampling

## 1. Wet Purge




- Carrier gas bubbled through the matrix
- Volatiles in matrix diffuse into carrier gas & are carried away

Typical flows: 40-50 mL/min. (can heat)

Typical purge time: 10-12 min.

# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Sequence


## 2. Dry Purge



- Trap is dried by purging with gas only
- Typical time: 1-4 min.

# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Sequence

## 3. Desorb Preheat



- Trap is heated without flow
  - Typical temp: 5° below desorb temperature
- Minimizes retention on the trap

# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Sequence


## 4. Desorb

- Trap is backflushed into column
- Typical time: 2 - 4 minutes
- Typical flow: 10-80 mL/min.

Typical temp: 180° - 250°C

Restek


[www.restekcorp.com](http://www.restekcorp.com)



# Flavor Profiling: Dynamic Gas Extraction Purge & Trap Sequence

## 5. Trap Bake


- Trap is baked clean with flow
- Typical time: 8+ minutes
- Typical temp: higher than desorb temperature
- Avoid overheating adsorbents





# Flavor Profiling: Dynamic Gas Extraction

## Typical Adsorbents for Purge & Trap



# Flavor Profiling: Dynamic Gas Extraction

## Requirements of a Trap

- Retention of polar & non-polar compounds
- Hydrophobic characteristics
- Reproducible desorption
  - Increasing levels of adsorbency
- Able to withstand a broad temperature range

# Flavor Profiling: Dynamic Gas Extraction


## Connecting a Purge & Trap to the GC Column

- Via the injection port
- Directly to the column
  - 1/16" union
  - Silica transfer line from 6 port valve directly to the column
- Low volume injector

# Interfacing to a GC/MS System


- Narrow bore column flows (0.5-1.3 mL/min.) permit a direct interface
- 0.53mm ID columns require:

Open Split Interface




Excess Flow

or Jet Separator



Vacuum

# Purge & Trap GC/MS of Volatiles



# Purge & Trap GC/MS of Volatiles

## GC Parameters

- Column: Rtx-5MS, 30m x 0.25mm x 1.0um
- Injector: 250°C, 20:1 split
- Carrier gas: Helium at 1 mL/min, constant flow
- Oven: 50 °C to 92 °C at 3 °C/min, to 220°C/min. at 20°C/min. (1 min. hold)

## MSD Parameters

- Temperature: 280°C
- Scan Range: 35-260, 1 min. solvent delay
- Ionization: EI @ 70eV


## Purge & Trap Parameters

- Concentrator: Tekmar LSC-3100 with Vocarb 3000 (type K) trap
- Purge: 10 min. at 40 mL/min, 60°C
- Dry purge: 3 min. at 40 mL/min.
- Desorb: 2 min. at 40 mL/min, 245°C

# Thermal Desorption Techniques

- Heat sample to drive volatiles into headspace
- Can trap and concentrate volatiles
  - Cryofocusing step
- Minimal sample preparation
- Commercially available units

# Short-Path Thermal Desorption System





# Summary of Extraction & Introduction Techniques

- Sampling techniques for solid and liquid samples
- Food & flavor samples often have extremely challenging matrices
  - Need to sample trace level components
  - Matrix can vary significantly from sample to sample
- Newer technologies focus on extracting flavor compounds from sample matrix without significant dilution