

Selecting an LC/MS Interface

Becky Wittrig, Ph.D.
RESTEK CORPORATION



LC/MS Interfaces

I. Background of LC/MS

I. Historical Perspective

II. Reasons for use

II. Interfaces

I. Transport devices

II. Particle Beam

III. Thermospray

IV. Atmospheric Pressure Interfaces



Historical Perspective

- Goldstein – 1886
 - Existence of positively charged particles
- Wein – 1898
 - Positively charged ions can be deflected in electrical and magnetic fields
- J.J. Thomson – 1913
 - Demonstrated isotopes of Neon
 - “Father of mass spectrometry”
- First GC-MS interface - 1960's

Historical Perspective

- First LC-MS interface developed - 1969
 - 1 μ L/min flow into an EI source
- Transport devices applied to LC/MS - 1970's
 - Loss of volatile components
 - Thermally-reactive compound losses
- Thermospray (TSI) gains popularity - 1983
 - 1.0 – 1.5 mL/min
 - Mobile phases consist primarily of an aqueous buffer

Historical Perspective

Atmospheric Pressure Interfaces (API)

- Early 1990's (commercialization)
- Now most common interface
- Electrospray (ESI)
 - Initial interfaces required lower flows (1-5 $\mu\text{L}/\text{min}$)
 - Able to produce multiply-charged molecules
- Atmospheric Pressure Chemical Ionization
 - Similar to Thermospray
 - "Solvent-mediated" ionization

Why Use Mass Spectrometry?


- Spectral resolution is possible
 - Chromatographic coelutions
- Compound identification from spectral data
 - Mass spectrum is very dependant on the ionization
 - Limited availability of LC/MS libraries
- High degree of specificity

LC/MS Interfaces

- I. Background of LC/MS
 - I. Historical Perspective
 - II. Reasons for use
- II. Interfaces
 - I. Transport devices
 - II. Particle Beam
 - III. Thermospray
 - IV. Atmospheric Pressure Interfaces



Components of an MS Detector



Ions
Made

Ions
Selected

Ions
Detected

LC/MS Systems

Interfaces

- Transport Devices
- Particle Beam
- Thermospray
- Atmospheric Pressure
 - Electrospray
 - APCI
- Others

Analyzers

- Quadrupoles
- Ion Trap
- MSⁿ
- High Resolution
- Time-of-Flight (TOF)
- Others

Ion Sources

- Analytes must be charged (ions)
 - Needed to separate and detect
- Where ionization occurs
 - In the interface
 - In a separate ion source
- Types of ionization
 - “Hard” - considerable fragmentation
 - “Soft” - molecular ion is main product
- Ions can be positive or negative

Requirements of LC/MS Interfaces

- Elimination of the mobile phase
 - Most difficult step
 - Can use splitters
 - Volatilized solvent vapor removed under vacuum
- Often where ionization occurs
- Vacuum required by mass analyzers


Challenges in Interfacing LC to MS

- Flow Rate Concerns
 - Differential pumping can only handle 2 $\mu\text{L}/\text{min}$ of water
 - For maximum sensitivity, want to use all of the eluent
- Use of Buffers and Additives
 - Non-volatile buffers a concern
 - Some additives suppress ionization
- Wide Range of Analytes
 - Many are nonvolatile, thermally labile

Transport Devices

- One of the first commercial interfaces
- Sample deposited onto a moving belt or wire
- Sample passes through multiple vacuum zones
 - Solvent elimination
- Sample is desorbed into source using heat
 - Electron impact ionization
- Belt/wire cycles back

Moving Belt Interface




From: Niessen

Particle Beam Interface


- Column effluent is nebulized
 - Pneumatic or thermospray nebulization
- Desolvation chamber is under a moderate vacuum
- A momentum separator is used for analyte enrichment
 - High MW compounds favored
- Analytes into the EI or CI source as small particles
 - Evaporative collisions with the walls

Particle Beam Interface



Analyte Enrichment in PB Interfaces

- Analyte Ion
- Solvent




Analyte enrichment with a molecular beam approach

Heavier molecules are in the core of the vapor jet and are sampled through the skimmer

Electron Impact Ionization


- LC interfaces with the ion source
- Electrons are “boiled” from a hot wire (filament), and accelerated (70eV)
- As electrons pass neutral molecules, they may remove outer shell electrons



Electron Impact Ionization

- Produces positively charged ions
- Fragmentation is generally significant for most molecules – “hard ionization”
 - Masses of these fragments is the information used in interpretation
- Efficiency of ionization is $1/10^5$
 - Bulk of molecules are removed by vacuum pump – Use Traps

Electron Impact Ionization (EI)




Thermospray Interface

- Nebulization of the eluent from a heated transfer tube
- Uses a “reagent gas”
 - Mobile phase buffer
 - Added buffer solution
 - Similar spectra to GC/MS CI
- Reagent gas is ionized
 - Volatilization
 - EI with high energy electrons
- Charge transfer to the analyte(s)

Thermospray Interface

- Effectively replaced transport systems
- Ionization in a medium pressure environment
 - Approx. 1000 Pa or 0.01 atm
- Inlet flows of 1 to 2 mL/min
 - Can use standard LC column flows
- Positive and negative ions are possible
- Temperature optimization is critical
 - Maintain gas phase


Thermospray LC/MS System




From: McMaster and McMaster

Thermospray Ionization

- Positive Ion Mode



- Negative Ion Mode



$gpa = \text{gas phase acidity}$

Typical Proton Affinities

| Compound | Proton Affinity |
|-------------------------|-----------------|
| Water | 723 |
| Methanol | 773 |
| Acetonitrile | 797 |
| Ethers, esters, ketones | 630-670 |
| PAHs | 710-800 |
| Carboxylic Acids | <800 |
| Alcohols | 750-840 |
| Peptides | 880-1000 |

From: Niessen

Thermospray Modes of Operation

- Solvent-Mediated CI
 - Typical mobile phase: MeOH or MeCN + 0.1M NH₄OAc
 - Ionization by ion evaporation (pre-formed ions) or ion-molecule reactions
- Discharge or Filament On Mode
 - High energy electrons (0.5-1 keV) ionize the reagent gas
 - Ion-molecule reactions


Electrospray Interface

- High electric potential applied to eluent from transfer capillary
 - Atmospheric pressure
 - Droplet formation
 - Ionization in the solution phase
- Orthogonal sampling of ions
 - Reduces contamination of the sampling orifice
 - Z-Spray devices

Electrospray Interface


- Higher flows now possible
 - Pneumatic + thermal nebulizers
- Applicable to volatile and nonvolatile analytes
 - Need to be ionizable
- Can create multiply-charged ions
 - Allows for analysis of large molecules

Electrospray LC/MS System




From: McMaster and McMaster

Electrospray Interface



Electrospray Needle Design




Example of an electro-spray needle design (coaxial flow)

The nebulizer gas disrupts the liquid surface so that small droplets are formed and then dispersed by the gas


Theory of API Electrospray

Charged Droplets → Analyte Ions



Courtesy of Agilent Technologies

Agilent 1100 LC/MSD Electrospray



Courtesy of Agilent Technologies

Electrospray Interface


- Best vaporization with higher % organic and lower flow rates
- Cluster ion formation possible
 - Solvent clusters
 - Analyte/salt clusters
- Salts and sample impurities can affect the response
 - TFA causes signal suppression
 - TFA anion masks the analyte ion

APCI Interface


- Atmospheric Pressure Chemical Ionization
- Initially investigated in 1974
 - Popular in late 1980's
- Uses an API (ESI) source
- Column effluent nebulized into heated vaporizer tube
- Solvent vapor acts as a reagent gas
 - Charge transfer to the analytes
- Can be very sensitive

APCI Interface

- Results in a chemical ionization spectrum
[M+H]⁺ or [M-H]⁻
- Products depend on equilibrium
(concentration) conditions
- Analytes must have sufficient proton affinities
- May be simplest interface to operate
- Liquid flow rates of 0.2-2.0 mL/min




APCI Interface



From: Scott


APCI Probe (Heated)



The vaporizer tube temperature is optimized for complete transfer to the vapor state


Theory of APCI

Analyte containing aerosol



Courtesy of Agilent Technologies


Agilent 1100 LC/MSD - APCI



Courtesy of Agilent Technologies



Agilent 1100 LC/MSD APCI Ion Source





Corona Needle

Vaporize into the gas phase and ionize the gas with a discharge

Courtesy of Agilent Technologies

API-ES vs. APCI for Triamterene



Comparison of the mass spectra typical of
Electrospray vs. APCI

Optimizing Mobile Phases for API-MS

- Ion-Pair Agent Alternatives
 - Use highly bases-deactivated silica columns
 - Use low pH (3-4) to reduce tailing
 - Use columns that retain based on polar interactions (e.g. CN, IBD)
- Using Ion-Pair Reagents
 - Use low amounts (<0.02%)
 - Use post-column addition to negate the effect of the ion-pair agent


Sample Considerations for LC/MS

- The Analyte Must Have Ionizable Groups
 - Amines
 - Carboxylic Acids
 - Ketones, Aldehydes
- For Best Sensitivity, Work at a pH Where the Analyte is Ionized
 - Neutral to basic pH (7-9) for acids
 - Acidic pH (3-4) for bases


Sample Considerations

- Positive Ion Mode
 - Analyte = $(M+H)^+$
- Negative Ion Mode
 - Analyte = $(M-H)^-$


Basic Compound – Sensitive in Positive Ion Mode



Acidic Compound – Sensitive in Negative Ion Mode



Selecting an Interface



“Advances in LC/MS”, Waters Corporation, Milford, MA.

Electron Ionization (EI)

- Analytes Suitable for EI
 - Small molecules with rings and double bonds
 - Compounds that would need derivatization for GC/MS
 - Pesticides, PAHs, natural products
- Compound Identifications
 - Fragmentation is possible
- Poor Detection Limits
- Will Not Tolerate Non-Volatile Buffers

Comparison of API vs. EI

- Atmospheric Pressure Ionization (API)

- MW confirmation
- Good for fragile compounds
- Able to fragment in the source
- Low (ppb) LODs in SIM mode


- Electron Ionization (EI)

- Ionization occurs in a vacuum
- Standard libraries are available
- Classical EI spectra (similar to GC/MS)
- Higher LODs (ppm to high ppb)

Atmospheric Pressure Ionization (API)

- Electrospray Ionization (ESI) Uses Solution Phase Ionization
- Atmospheric Pressure Ionization (APCI) Uses Gas Phase Ionization
- Products are $[M+H]^+$ and $[M-H]^-$, adducts
- Suitable for Analyzing Drugs, Small Molecules, Dyes, Peptides
- Good Sensitivity
- Thermal Degradation is Possible (APCI)

Adduct Formation with API



Adducts can form between polar molecules and sample or solvent components. For example, adducts with Na^+ , K^+ , NH_4^+ , MeOH, MeCN, and H_2O are common.

Electrospray Ionization

- + Good molecular weight information, including high MW compounds
- + Can be used for volatiles, nonvolatiles, ionic/polar compounds
- + Good sensitivity
- Need relatively low flow rates
- Need to be able to form ions in solution
- Limited structural information
- Problems with high aqueous and buffer solutions

Atmospheric Pressure Chemical Ionization

- + Gives molecular weight information
- + Easy to use, rugged
- + Can use higher LC flow rates (up to 2 mL/min)
- Thermal degradation can occur
- Limited structural information
- Not appropriate for higher MW (e.g., >1000 Da)

MALDI

- Matrix-assisted laser desorption/ionization
 - Sample is deposited on a target and co-crystallized with a solid matrix (dihydroxybenzoic acid)
 - Desorption/ionization occurs using a laser such as Nd-YAG (266nm, v^4)
 - Energy transferred to matrix, then analyte
 - Useful in excess of 200kDa (biomacromolecules)

Fast Atom Bombardment (FAB)

- Cf-FAB = Continuous Flow FAB
- Column eluent mixes with a matrix (glycerol)
- Eluent + matrix deposited on a target
- Analyte film hit with fast atoms or ions
- Only low flow rates used

Acknowledgements

Charlie Schmidt - Thermo Instruments

Agilent Technologies

Liquid Chromatography-Mass Spectrometry, Marcel Dekker, New York, NY, W.M.A. Niessen ed. 1999

Chromatographic Detectors, Marcel Dekker, New York, NY, Raymond P.W. Scott ed. 1996

GC/MS A Practical Users Guide, Marvin C McMaster and Christopher McMaster, Wiley-VCH, New York, NY 1988