

GC-MS Analysis of Polychlorinated Biphenyl Congeners Using a New Capillary GC Column

Frank L. Dorman, Gary B. Stidsen, Chris M.
English and Lydia Nolan

Restek Corporation, Bellefonte, PA

Jack Cochran

LECO Corporation, Las Vegas, NV

Abstract

- The analysis of polychlorinated biphenyls (PCB's) by GC can be a very challenging analysis. Numerous researchers have demonstrated various separations on the common commercially available GC stationary phases, but all of these have had limitations when it comes to the separation of all 209 possible PCB congeners. While it may be argued that the complete separation of all 209 possible congeners is not necessary, separation of the environmentally significant and the toxic congeners is important.
- The separation requirements for the analysis of PCB congeners differ considerably depending on the type of detector used. Separation needs when using a mass spectrometer are different than when using standard GC detectors.
- This presentation will discuss the separation needs for PCB congeners when using a mass spectrometer, and demonstrate an improved separation of PCB congeners, using a new capillary stationary phase developed for this purpose.

Experiment


- Determine retention times for all 209 PCBs on a new low-bleed capillary column
 - AccuStandard 9 mixtures
 - Hexachlorobenzene for relative retention times
- Analyze mix of Aroclors 1221:1242:1254:1262
 - 1:2:2:2 ratio
 - Show chromatogram with congeners labeled
- GC-MS results obtained using a time-of-flight mass spectrometer to allow for possible spectral deconvolution


GC-TOFMS Conditions


- Gas Chromatography
 - One microliter splitless 250°C, 60 sec valve
 - 60 m x 0.18 mm x 0.18 µm Restek Rtx-PCB column
 - Constant flow helium, 1.5 mL/minute
 - 70°C (1 min), 50°/min to 120°, 3°/min to 360° (1 min)
- Mass Spectrometry LECO Pegasus III
 - Source temperature: 225°C
 - Electron ionization: 70 eV
 - Stored mass range: 120 to 520 u
 - Acquisition rate: 5 spectra/sec


Aroclor Mix Chromatogram


- Overlaid ion chromatograms are plotted
 - For example, hexachlorobiphenyl is 360
 - Peaks are labeled with congener numbers
 - Text size and format is only for space and not meant to indicate anything about PCBs, except:
 - **European indicator congeners** are in red
- Chromatograms are not to one scale
 - Each panel is adjusted so most of the PCB peaks can be seen




























Important Congener Separations

- European indicator PCBs
- May coelute on 5% phenyl methylsilicone and other column phases





Hexachlorobiphenyls




Table Key

- Congener classification according to their significance in Aroclors 1242, 1254, or 1260
 - **Bold Underline**: Congener > 1.0 Wt%
 - **Bold**: Congener 0.05 to 1.0 Wt%
 - *Italics*: Trace or undetected congener
- Congener naming according to Guitart et al. by chlorine positions
 - 107: 234-35
 - 108: 2346-3
 - 109: 235-34
 - 199: 2345-2356
 - 200: 23456-236
 - 201: 2346-2356
- RRT is relative retention time
- Resolutions of significance where MS won't help are

PCB#	CI#	RT (sec)	RT (min)	RRT
1	1	1403.3	23.39	0.7206
2	1	1601.3	26.69	0.8222
3	1	1642.1	27.37	0.8432
4	2	1678.1	27.97	0.8617
10	2	1679.7	28.00	0.8625
9	2	1808.3	30.14	0.9285
7	2	1815.3	30.26	0.9321
6	2	1849.3	30.82	0.9496
5	2	1882.3	31.37	0.9665
8	2	1894.3	31.57	0.9727
19	3	1953.3	32.55	1.0030
14	2	1963.3	32.72	1.0079
30	3	1993.5	33.22	1.0236
18	3	2053.3	34.22	1.0543
11	2	2061.1	34.35	1.0582
17	3	2067.3	34.46	1.0615
12	2	2088.1	34.80	1.0721
27	3	2093.3	34.89	1.0748
13	2	2109.3	35.15	1.0829
24	3	2110.9	35.18	1.0839

PCB#	CI#	RT (sec)	RT (min)	RRT
<u>16</u>	3	2135.1	35.58	1.0963
<u>32</u>	3	2148.7	35.81	1.1033
<u>15</u>	2	2158.9	35.98	1.1084
<u>34</u>	3	2176.1	36.27	1.1173
23	3	2187.3	36.45	1.1234
54	4	2188.9	36.48	1.1237
<u>29</u>	3	2202.7	36.71	1.1309
50	4	2237.1	37.28	1.1486
<u>26</u>	3	2239.3	37.32	1.1498
<u>25</u>	3	2252.7	37.54	1.1567
<u>53</u>	4	2286.5	38.11	1.1738
<u>31</u>	3	2290.7	38.18	1.1762
<u>28</u>	3	2303.5	38.39	1.1828
<u>33</u>	3	2309.5	38.49	1.1858
<u>51</u>	4	2312.9	38.55	1.1874
21	3	2312.7	38.54	1.1874
<u>20</u>	3	2313.3	38.56	1.1877
<u>45</u>	4	2349.5	39.16	1.2064
<u>22</u>	3	2364.3	39.40	1.2140
73	4	2374.3	39.57	1.2189

PCB#	CI#	RT (sec)	RT (min)	RRT
46	4	2374.5	39.58	1.2193
69	4	2391.5	39.86	1.2279
43	4	2397.7	39.96	1.2312
36	3	2398.9	39.98	1.2317
52	4	2406.9	40.11	1.2359
48	4	2417.3	40.29	1.2412
49	4	2426.5	40.44	1.2458
104	5	2440.3	40.67	1.2529
47	4	2444.7	40.75	1.2552
65	4	2448.5	40.81	1.2571
75	4	2451.3	40.85	1.2584
62	4	2450.7	40.84	1.2586
39	3	2452.9	40.88	1.2598
38	3	2461.3	41.02	1.2637
44	4	2487.3	41.45	1.2772
59	4	2501.3	41.69	1.2842
42	4	2506.7	41.78	1.2870
71	4	2521.5	42.02	1.2947
35	3	2527.9	42.13	1.2978
41	4	2535.7	42.26	1.3020

PCB#	CI#	RT (sec)	RT (min)	RRT
72	4	2542.1	42.37	1.3052
96	5	2545.1	42.42	1.3067
103	5	2552.3	42.54	1.3105
68	4	2560.5	42.67	1.3150
64	4	2561.9	42.70	1.3153
40	4	2569.9	42.83	1.3195
37	3	2581.3	43.02	1.3254
100	5	2581.9	43.03	1.3255
94	5	2589.1	43.15	1.3297
57	4	2594.1	43.23	1.3319
67	4	2613.5	43.56	1.3420
58	4	2614.7	43.58	1.3426
102	5	2618.9	43.65	1.3446
61	4	2635.9	43.93	1.3533
98	5	2637.9	43.96	1.3544
93	5	2646.1	44.10	1.3586
76	4	2646.1	44.10	1.3587
63	4	2650.9	44.18	1.3610
121	5	2653.9	44.23	1.3630
95	5	2657.5	44.29	1.3646

PCB#	CI#	RT (sec)	RT (min)	RRT
88	5	2658.9	44.31	1.3656
<u>74</u>	4	2667.9	44.46	1.3699
155	6	2675.7	44.59	1.3742
<u>70</u>	4	2678.3	44.64	1.3752
<u>91</u>	5	2688.7	44.81	1.3805
<u>66</u>	4	2695.7	44.93	1.3842
80	4	2709.7	45.16	1.3917
55	4	2723.3	45.39	1.3984
<u>92</u>	5	2726.5	45.44	1.3999
<u>89</u>	5	2743.7	45.73	1.4087
<u>84</u>	5	2745.9	45.77	1.4100
<u>101</u>	5	2750.9	45.85	1.4124
90	5	2751.5	45.86	1.4125
<u>56</u>	4	2751.5	45.86	1.4128
113	5	2757.3	45.95	1.4157
<u>99</u>	5	2772.9	46.21	1.4238
150	5	2776.1	46.27	1.4258
<u>60</u>	4	2776.7	46.28	1.4258
152	6	2798.9	46.65	1.4370
<u>119</u>	5	2802.3	46.71	1.4388

PCB#	CI#	RT (sec)	RT (min)	RRT
83	5	2807.5	46.79	1.4416
125	5	2811.3	46.85	1.4438
112	5	2812.1	46.87	1.4440
86	5	2817.7	46.96	1.4467
108	5	2821.5	47.02	1.4488
145	6	2825.9	47.10	1.4509
<u>97</u>	5	2828.9	47.15	1.4524
148	6	2835.5	47.26	1.4563
79	4	2841.9	47.36	1.4591
116	5	2855.5	47.59	1.4665
154	6	2867.7	47.79	1.4722
<u>87</u>	5	2868.3	47.80	1.4728
78	4	2870.3	47.84	1.4737
111	5	2871.1	47.85	1.4746
<u>136</u>	6	2877.3	47.96	1.4773
<u>117</u>	5	2878.9	47.98	1.4780
<u>115</u>	5	2883.5	48.06	1.4806
<u>85</u>	5	2889.1	48.15	1.4833
120	5	2895.7	48.26	1.4869
<u>110</u>	5	2906.7	48.44	1.4925

PCB#	CI#	RT (sec)	RT (min)	RRT
81	4	2927.7	48.79	1.5030
151	6	2929.1	48.82	1.5039
135	6	2937.3	48.96	1.5082
82	5	2948.3	49.14	1.5139
144	6	2949.1	49.15	1.5142
147	6	2966.3	49.44	1.5231
149	6	2968.3	49.47	1.5241
77	4	2972.3	49.54	1.5261
140	6	2980.3	49.67	1.5306
139	6	2981.3	49.69	1.5307
124	5	2988.9	49.81	1.5344
143	6	2990.5	49.84	1.5359
107	5	3005.3	50.09	1.5430
123	5	3009.5	50.16	1.5452
109	5	3011.3	50.19	1.5462
188	7	3011.5	50.19	1.5467
134	6	3018.5	50.31	1.5498
106	5	3018.9	50.31	1.5500
133	6	3023.1	50.38	1.5521
142	6	3027.5	50.46	1.5549

PCB#	CI#	RT (sec)	RT (min)	RRT
<u>118</u>	5	3032.5	50.54	1.5570
<u>131</u>	6	3034.3	50.57	1.5580
<u>184</u>	7	3040.7	50.68	1.5612
<u>165</u>	6	3044.1	50.73	1.5628
<u>146</u>	6	3050.3	50.84	1.5663
<u>122</u>	5	3060.1	51.00	1.5710
<u>161</u>	6	3060.3	51.00	1.5712
<u>168</u>	6	3072.5	51.21	1.5780
<u>114</u>	5	3076.1	51.27	1.5794
<u>153</u>	6	3077.1	51.28	1.5800
<u>132</u>	6	3094.1	51.57	1.5887
<u>179</u>	7	3109.3	51.82	1.5966
<u>141</u>	6	3125.7	52.10	1.6050
<u>176</u>	7	3137.5	52.29	1.6109
<u>105</u>	5	3140.9	52.35	1.6126
<u>137</u>	6	3150.7	52.51	1.6177
<u>186</u>	7	3157.1	52.62	1.6211
<u>127</u>	5	3157.7	52.63	1.6212
<u>130</u>	6	3168.7	52.81	1.6267
<u>164</u>	6	3171.1	52.85	1.6283

PCB#	CI#	RT (sec)	RT (min)	RRT
<u>178</u>	7	3187.7	53.13	1.6368
<u>138</u>	6	3191.7	53.19	1.6389
<u>163</u>	6	3196.7	53.28	1.6411
<u>160</u>	6	3198.9	53.31	1.6429
<u>129</u>	6	3203.9	53.40	1.6450
<u>158</u>	6	3208.5	53.48	1.6473
<u>182</u>	7	3209.7	53.49	1.6480
<u>175</u>	7	3210.5	53.51	1.6482
<u>187</u>	7	3221.1	53.68	1.6540
<u>183</u>	7	3243.3	54.06	1.6654
<u>166</u>	6	3264.1	54.40	1.6759
<u>159</u>	6	3277.3	54.62	1.6828
<u>126</u>	5	3290.3	54.84	1.6893
<u>185</u>	7	3293.3	54.89	1.6909
<u>162</u>	6	3295.9	54.93	1.6927
<u>128</u>	6	3307.9	55.13	1.6984
<u>174</u>	7	3307.7	55.13	1.6984
<u>202</u>	8	3321.1	55.35	1.7050
<u>167</u>	6	3322.7	55.38	1.7060
<u>181</u>	7	3328.9	55.48	1.7091

PCB#	CI#	RT (sec)	RT (min)	RRT
<u>177</u>	7	3348.5	55.81	1.7194
201	8	3352.3	55.87	1.7210
<i>204</i>	8	3353.9	55.90	1.7220
<u>171</u>	7	3366.1	56.10	1.7284
173	7	3379.7	56.33	1.7354
197	8	3383.1	56.39	1.7371
172	7	3403.1	56.72	1.7474
<i>192</i>	7	3414.7	56.91	1.7534
<u>156</u>	6	3415.3	56.92	1.7535
157	6	3430.1	57.17	1.7611
<u>180</u>	7	3432.1	57.20	1.7623
193	7	3437.5	57.29	1.7651
200	8	3447.7	57.46	1.7700
191	7	3456.5	57.61	1.7747
198	8	3534.3	58.90	1.8148
<u>199</u>	8	3539.9	59.00	1.8177
<u>170</u>	7	3546.9	59.12	1.8212
190	7	3564.7	59.41	1.8302
196	8	3564.7	59.41	1.8304
<u>203</u>	8	3569.3	59.49	1.8328

PCB#	CI#	RT (sec)	RT (min)	RRT
169	6	3583.5	59.72	1.8399
208	9	3643.7	60.73	1.8708
207	9	3676.5	61.28	1.8876
189	7	3683.3	61.39	1.8911
195	8	3692.9	61.55	1.8962
<u>194</u>	8	3766.9	62.78	1.9342
205	8	3791.5	63.19	1.9468
206	9	3872.3	64.54	1.9884
209	10	3952.9	65.88	2.0295

Conclusions

- Greater than 140 Aroclor PCBs can be determined with the new column and mass spectrometry
- The seven European indicator PCBs can be determined individually
- Spectral deconvolution available with TOFMS allows qualitative identification for coeluting PCBs