

Sample VIALS GC AutoSamplers Liquid & HeadSpace Vials



Crimp-Top Vials, 11mm Crimp Top
 Short-Cap 9mm Screw Thread Vials

These fits ALL popular
 2.0ml, 12x32mm Crimp-Top
 Vial based Autosamplers

Dynatech, Finnigan, Gilson, Agilent, Hitachi, LDC, L.E.A.P,
 Shimadzu, Spark, Spectra Physics, Varian, Waters

Vial Options : Clear or amber
 Deactivated Silcote CL7 or untreated
 Limited volume vials in clear glass or polypropylene
 Limited Volume Inserts

Seals - Standard Silver Al crimp seal - colors optional
 With or without white graduated marking spot
 Optional (to order) color marking spots

Septa - PTFE/Butyl or PTFE/Silicone or to order PTFE/Silicone/
 PTFE

Available seperately or in kits vial/cap/septa
 100pk or 1000pk

Other Vial Types include . . .

Crimp-Top WISP 4.0ml WISP48 Style 15x45mm Snap Seal
 Screw-Thread 4.0ml 48 WISP Style 13mm/425, Step
 2.0ml, 8mm/425 Thread
 2.0ml, 10ml/425 Thread, Step Design
 Short-Cap, 9mm Screw Thread Vials

HeadSpace AutoSampler Vials

Clear glass
 6, 10, 20 ,27ml
 Round or flat-bottom
 PTFE/Silicone or PTFE/Butyl(gray) Septa
 Steel (colored caps) for HeadSpace
 Autosamplers with magnetic head pick-up.

Enquire about the
HeadSpace Autosampler
 (for ANY GC)

**Model
 HS200**

Snap-Rack Storage System

50 vial polypropylene rack
 for 2.0 ml 12x32 vials
 double-rack stackable

Aluminium Seal Crimper / Decapper

8mm, 11mm, 13mm 20mm versions

MinInert Valves for Sample Vials

On/off Valves for sample vials
 minimum exposure to septa rubber
 For screw cap and crimp top vials
 sizes: 13, 15, 18, 20 and 24mm caps

These products are available through Restek, Ansys and Chromos
 Please see our Website <http://www.chromtech.net.au>
 or ask for printed catalog / information

Snap Seal Vials
 11mm Crimp top
 Short Cap Screw Thread Vials
 8mm Crimp Top Vials
 HeadSpace Vials
 8mm and 10mm Screw Thread Vial
 Seals
 PolyCrimp 11mm
 11mm Snap-Top Cap
 Septa
 PTFE/Silicone
 PTFE/Butyl
 Crimper / Decapper
 MinInert Valves for Vials

Deactivated Vials from Restek
 Magnetic Head Pick-Up Vials Col-
 ored-steel Caps
 Limited Volume Vial Inserts

Syringe Filters

Maximum sample recovery
 Low extractables
 Leak-Proof design
 Low Cost

Sizes: 4mm, 13mm, 25mm
Material: Nylon, PTFE, PVDF, Regener-
 ated Cellulose
 25mm : also in Polyethersulfone
 0.25um, 0.45um
 Available in 100-pk

**Economical
 LOW Cost**



Applications

Nylon - hydrophilic, compatible with aque-
 ous and organic solvents, limited use with
 THF, chlorinated solvent; but NOT strong
 acids/bases

PTFE - hydrophobic for most organic
 solvents; acids/bases.

Regenerated Cellulose - hydrophilic low
 protein-binding for filtration of proteins and
 biological fluids

Sterile versions available to special order –
 for medical use individually wrapped/
 packaged.

