

Supplies & Accessories for ASE[®] Systems

Extraction Cell Parts, Collection Vials, PEEK[®] Washers, Filters



REST

CHROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

into Reality[™]

Replacement Parts for ASE® Systems

Accelerated solvent extraction has become a common technique for fast and reliable extraction of organic materials from solid matrices. Restek offers a wide range of replacement parts to keep your extraction system running smoothly. All parts are economically priced to save you money and are designed to meet or exceed the performance of the original manufacturer's parts.

Replacement Parts for ASE® 200 Systems

Extraction Cell Parts for ASE® 200 Systems

- Choose original equipment-equivalent stainless steel, or Siltek® deactivation for improved inertness.
- Inner surfaces polished, for easier cleaning.
- Caps include frit, washer, PTFE O-ring, and threaded insert.

save money
buy in Bulk!



Funnel for Filling ASE® 200 Cells

Description	Similar to		Stainless Steel		Siltek® Treated	
	Dionex part #	qty.	cat.#	qty.	cat.#	
Extraction Cell Body for ASE® 200, 1mL	054973	ea.	26110	ea.	26111	
Extraction Cell Body for ASE® 200, 5mL	054974	ea.	26112	ea.	26113	
Extraction Cell Body for ASE® 200, 11mL	048820	ea.	26114	ea.	26115	
Extraction Cell Body for ASE® 200, 22mL	048821	ea.	26098	ea.	26099	
Extraction Cell Body for ASE® 200, 33mL	048822	ea.	26116	ea.	26117	
Replacement Extraction Cell End Caps for ASE® 200	049450	2-pk.	26096	2-pk.	26097	
Cap Inserts for ASE® 200		2-pk.	26166			
Replacement Frits for ASE® 200	049453	10-pk.	26100	10-pk.	26101	

Description	Similar to		
	Dionex part #	qty.	cat.#
Snap Rings for Caps for ASE® 200	049456	10-pk.	26184
Funnel for ASE® 200	056958	ea.	26180
PTFE O-Rings for ASE® 200 & ASE® 300 Caps	049457	100-pk.	26187
Viton® O-Rings for ASE® 200 & ASE® 300 Caps	056325	50-pk.	26188



PEEK Washers for ASE® 200 Systems

- Meet original equipment manufacturer's performance.
- Available in 3 quantities.

Description	Similar to Dionex part #	12-pk.	48-pk.	250-pk.
PEEK® Washers for ASE® 200	049454	25256	25257	26120



20mm Filters for ASE® 200 Extraction Cells

- Cellulose or glass fiber construction.
- Cellulose filters available in economical 1000-packs.

Description	Similar to Dionex part #	qty.	cat.#
Cellulose Filters for ASE® 200	049458	100-pk.	26118
Cellulose Filters for ASE® 200	049458	1000-pk.	26190
Glass Fiber Filters for ASE® 200	047017	100-pk.	26119



60mL Sample Collection Vials for ASE® 200 Systems

- Cleaned, assembled, and ready to use.
- Clear or amber glass.
- Caps and Teflon®-lined septa included.

Description	Similar to Dionex part #	qty.	cat.#
60mL Collection Vials, Clear Glass, for ASE® 200	048784	72-pk.	26121
60mL Collection Vials, Amber Glass, for ASE® 200	048781	72-pk.	26122

Replacement Parts for ASE® 300/100 Systems

Extraction Cell Parts for ASE® 300/100 Systems

- Choose original equipment-equivalent stainless steel, or Siltek® deactivation for improved inertness.
- Inner surfaces polished, for easier cleaning.
- Caps include frit, washer, PTFE O-ring, and threaded insert.

Description	Similar to Dionex part #	Stainless Steel		Siltek® Treated	
		qty.	cat.#	qty.	cat.#
Extraction Cell Body for ASE® 300, 10mL		ea.	26172	ea.	26173
Extraction Cell Body for ASE® 300, 34mL	056646	ea.	26176	ea.	26177
Extraction Cell Body for ASE® 300, 66mL	056696	ea.	26178	ea.	26179
Extraction Cell Body for ASE® 300, 100mL	056693	ea.	26132	ea.	26133
Replacement Extraction Cell End Caps for ASE® 300	056921	2-pk.	26170	2-pk.	26171
Cap Inserts for ASE® 300		2-pk.	26167		
Replacement Frits for ASE® 300/100		6-pk.	26174	6-pk.	26175

Description	Similar to Dionex part #	qty.	cat.#
Funnel for ASE® 300	—	ea.	26169
PTFE O-Rings for ASE® 200 & ASE® 300 Caps	049457	100-pk.	26187
Viton® O-Rings for ASE® 200 & ASE® 300 Caps	056325	50-pk.	26188



PEEK® Washers for ASE® 300 Systems

- Meet original equipment manufacturer's performance.
- Available in 2 quantities.

Description	Similar to Dionex part #	12-pk.	48-pk.
PEEK® Washers for ASE® 300	061687	25393	25394



30mm Filters for ASE® 300 Extraction Cells

- Cellulose or glass fiber construction.

Description	Similar to Dionex part #	qty.	cat.#
Cellulose Filters for ASE® 300	056780	100-pk.	26168
Glass Filters for ASE® 300	056781	100-pk.	26189



250mL Sample Collection Bottles for ASE® 300 Systems

- Cleaned, assembled, and ready to use.
- Clear glass.
- Caps and Teflon®-lined septa included.

Description	Similar to Dionex part #	qty.	cat.#
250mL Collection Bottles, Clear Glass, for ASE® 300	056284	12-pk.	26191



Tools for ASE® Systems

- Use with either ASE® 200 or ASE® 300 systems.



Description	Similar to Dionex part #	qty.	cat.#
Internal Retaining Ring Pliers for ASE® Systems	—	ea.	26185
Filter/O-Ring Insertion Tool Kit for ASE® 100/200/300 (Includes plug adapter handle and insertion tools)	—	kit	26181
Replacement Plug Adapter Handle for ASE® 100/200/300	—	ea.	26182
Replacement Insertion Tools for ASE® 100/200/300	—	4-pk.	26183



Diatomaceous Earth

- Improves extraction efficiency.
- Adsorbs moisture from samples.

Description	Similar to Dionex part #	qty.	cat.#
Diatomaceous Earth, 30/40 mesh	062819	1kg	26033

Please visit us at www.restek.com, or see our general catalog, for additional sample preparation products, including solid phase extraction tubes, bulk absorbents, and syringe filters.

Restek Trademarks: Siltek, Turning Visions into Reality, Restek logo

Other Trademarks: ASE (Dionex Corporation); PEEK (Victrex plc); Teflon, Viton (E.I. du Pont de Nemours & Co., Inc.)



Lit. Cat.# 59921A-INT

© 2006 Restek Corporation.

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
Website NEW : www.chromalytic.com.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au

