

Adjustable Hand Crimper 13mm, 20mm

Adjustable Hand Crimpers



The **Adjustable Crimper** has an adjustment feature which will allow it to be used with a variety of septa, discs, stoppers or liners with a thickness range of 0.25 to 3.0mm. The unit looks like the standard hand crimper but it is supplied with a wrench and special instructions for making adjustments.

The user may adjust the stroke or movement for each application as needed.

Available for 13mm and 20mm Diameter Crimp Vials Standard Caps and Flip-off Cap types

Ordering ;

CT-13/20Crimp-SA and CT-13/20Crimp-FOA

Chromalytic Technology Pty Ltd (Australia).
Ph.03 9762 2034 sales@chromtech.net.au