

# Biobase Meihua Trading Co., Ltd

Product Description

## Biobase 5L/h 10L/h 20L/h Electric-heating Distilled Water Device Laboratory Water Distiller

Reference FOB Price:

Min. Order:1 Piece

Port:Qingdao, China

Production Capacity: > Unit/Units Per Year

Payment Terms:L/C, T/T

Certification:ISO3834, ISO9001, CE

Wall Thickness:18mm

Installation:Hidden Type

Material:Stainless Steel

Type:Spherical Water Tank

Water Output:≥5 L/H

Hot



### Features:

1. All these products are made from stainless steel.
2. Condenser is made of seamless stainless

-steel pipes provide very good cooling effect & high of water production capacity.

3. The cooled water return to the vaporizing drum can save energy and water.

4. Drain valve can empty the water in container, minimize spillage.

5. When lack water, still will auto stop heating, and when the return water reach the normal level, then continue heating.

6. Immerge electric heating tube made from copper,for better heating efficiency.

### TECHNICAL PARAMETERS

Model	WD-5	WD-10	WD-20
Water Output	≥5 L/H	≥10 L/H	≥20 L/H
Consumption	5 KW	7.5 KW	15 KW
PowerSupply	AC110/220V,50/60Hz	380V/50Hz	380V/50Hz
ExternalSize (Dia.*H)mm	310*700mm	355*820mm	405*677mm
Package size(W*D*H)	360*360*820mm	360*350*900mm	540*390*740mm
Gross Weight	7.5kg	8.5kg	11.5kg

### Company Information

BIOBASE, Founded in 1999 and been through high and stable development for 20 years.

**Main Products:** Biological Safety Cabinet, Fume Hood, Laminar Flow Clean Bench, Refrigerator, Incubator, Autoclave, Centrifuge etc.

**HROM**alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

**Australian Distributors**  
Importers & Manufacturers  
www.chromtech.net.au

Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA