

# **Bottle Feed Line Kit for Filtration . . .**

**CT-BFLK**

connect to . . . Vacuum Pump  
or Nalgene FastCap Filter



Syringe Filter  
SFNY025045N

2x FTLL210-6

CT-TUBPVC05  
5mmID PVC 1metre  
LDPE1031 1ft

LastDrop 5um Filter

**Order GL-45 Bottle/Blank Cap Separately**

CT-BOT-GL45-1000(shown)

also -2000 .-500 ..-250, -100ml