

NOW SHIPPING!
EPA Approved Colitag™

colitag

The First Definitive Water Test System

Unprecedented Accuracy & Reliability
All In A Simple, One Step Process

The First Definitive Water Test System for Coliforms & E. coli.

CPI International is shipping the breakthrough Colitag™ water test system – the first truly definitive system for total coliforms and *E. coli*. One simple step. One accurate reading. Colitag™ gives you maximum confidence in your test results.

- Exceptional accuracy – eliminates false negatives.
- Patented self-adjusting pH resuscitates chlorine-damaged bacteria.
- Simultaneous testing for MUG-positive and MUG-negative *E. coli*.
- One step testing – just combine with water sample!
- Detects 1 CFU of coliform or *E. coli* in a 100mL sample.

See page 6 for more information.



**If You're Impressed With
Colitag™ Now, Wait Till You
Try It For Yourself... FREE!**

**Call 800-878-7654 x302
for details of this risk-
free offer of a free
sample of Colitag™**