

Calidus microGC from FalconFAST

Dual Columns

Dual "Oven"

400 degC Max

FAST

to 300degC per min

300Watts

24VDC or 220VAC

microFID

Micro TCD

Compact

43x21x29cm

Capillary

PLOT

microPacked

Restek MXT Columns



Typical Analysis . . . in 90seconds C2-C60+



Process Control GC
NeSSi Cart

UNIQUE !

Targetted Sample
Enhancement
(TSE)



larger peaks
x15 to x50

water/ methanol/ CO & CO2

C2-C4s; benzene

Post Column Cryo

VORTEX Cooling



Fast Laboratory GC - 100Hz 24bit Data System

