

## Portable air Compressor DA5001/4C



Model No. : DA5001/4C

Brand Name : DYNAIR



### Product Description

#### CHARACTERISTICS:

- OIL FREE: Do not need any lubricated oil, and harmless to the human body.
- SUPER SILENT: Noise level lower than 56DB(A) to reduce noise pollution.
- MULTI-PHASE FILTERATION: Advanced multi-phase filtration to ensure extremely clean air, and lengthen service life.
- SAFETY USING: Air compressor motors equipped with thermal prevention device to avoid overheating and protect motors.
- CLEAN TANK: All air receiver tanks have internal epoxy coating applied to avoid corrosion.
- AUTOMATIC ADJUSTMENT: High quality pressure switch used to control the power of air compressor automatically by adjustment of pressure level.
- PRESSURE ADJUSTMENT: Working pressure could be adjusted to meet the demand of different supporting equipments.
- HIGH DURABILITY: Serviceable time exceeds 15,000 hours.
- SIMPLE OPERATION: Using directly when connecting with power.
- EASY MAINTANANCE: No need any lubricated oil.
- WARRANTY: 2 years.

Payment Terms :	TT in advance/ Paypal/ LC at sight
Packing :	CARTON/WOODEN BOX
Lead Time :	within one month after receive order confirmation
Standards Certificate :	CE/FDA/ISO9001/ISO13485
Specifications :	Power: 550W [3/4HP] Air Flow: 115L/M [4.03CFM] Noise Level: 48dB(A) Max. Pressure: 8Bar [0.8Mpa] Air Tank: 4L [0.88Gallon] Gross/Net Weight: 25.5/22KG [56.1/48.4Pound] Package Dimension: 580*250*550mm [184.8*82.5*181.5inch]