

Silent Air Compressor DA5001CS



Model No. : DA5001CS

Brand Name : DYN AIR

Oil-Less, 4CFM, 22L
78/55kg net
500x650x770mm

Product Description

The oil-less DynAir compressors are available in protective metal cabinets. The cabinets reduce the sound emission from the already quiet oil-less from the already quiet oil-less compressor to a level that is approximately one quarter of the level of basic compressors. The cabinets are designed individually for each model with an aesthetic look and with usability, hygiene as well as easy maintenance in mind. All cabinets are equipped with the universal wheel with foot-brake, enabling mobility and re-arrangement of installations and thus providing a high degree of flexibility.

The cabinets are treated with a matt textured powder coating which is hard wearing and enables easy cleaning. The color of the cabinets is making the compressors blend naturally with the equipment found in most laboratory, medical and other environments.

Lead Time : within one month after receive order confirmation

Standards Certificate : FDA/CE

Specifications :
Power (W):550
Air displacement (L/M):115
Noise db(A):45
Maximum pressure (Kpa):800
Gasholder capacity (L):22
Net weight (KG):54.5
Gross weight (KG):78
Dimension (L/W/H)mm: 500*650*770

Payment Details : paypal \ L/C

Min Order : 1PC

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd

Website NEW : www.chromalytic.net.au E-mail : info@chromtech.net.au Tel: 03 9762 2034 . . . in AUSTRALIA

Australian Distributors
Importers & Manufacturers
www.chromtech.net.au