

# Accessories

Standard External Control Module for PP BT50-1J, BQ50-1J, BT100-1J, BT300-1J, BT100-1I, BT300-1F



### Applicable Drive Types


- BT50 - 1J
- BQ50 - 1J
- BT100 - 1J
- BT300 - 1J
- LEAD - 2
- BT100 - 1L
- BT300 - 1F

### Four Kinds of Control Modes


- 0 - 5V analog input
- 0 - 10V analog input
- 4 - 20mA analog input
- 0-10KHz Pulse Input
- RS 485

ACCESSORY


### ■ Pins Definition Of Standard External Control Interface



### ■ RS 485




### ■ Pulse Input External Control Interface



TTL pulse speed control input  
Max. rpm corresponds to 10KHz

### ■ 0 to 5V, 0 to 10V, 4 to 20mA External Control Interface



## Dispensing Controller



Type	Part Number	Function/Capability Introduction
FK-1C	0301001	Dispensing controller can be equipped to many peristaltic pumps to realize dispensing automatically.

## Footswitch



Type	Part Number	Applicable Products
JK-1(DB-9)	0301003	Peristaltic Pump
JK-2	0301004	Dispensing Controller
JK-3(DB-15)	0301005	Peristaltic Pump