

Dispensers

Bottle Top Dispenser



• Microlit Bottle top dispenser is another genius of its strong R&D base. Responding to the needs of the present generation, Microlit is proud to offer its NEW BOTTLE TOP DISPENSERS which is an unique combination of affordable pricing and high performance.

Vial Filling

Simple design clubbed with variety of useful features makes Microlit Bottle top Dispensers an obvious choice for all routine dispensing needs.

High precision and accuracy is ensured through several stages of strict quality checks during manufacturing process. Each instrument is individually calibrated in accordance with ISO 8655 standards and comes with Individual Calibration certificate.



FLUID DISPENSING SYRINGE

Vial Filling



The Syringe fluid dispensing system is effectively used when repetitious filling of a syringe to the same volume is needed. It is a manual, sterile device that has a two-way valve that attaches to both the syringe and the solution to be drawn into the syringe. Each time the grip is released, the syringe fills with solution. Each time the grip is compressed, the syringe volume is expelled. The Cornwall Syringe System is a sterile, disposable system that includes spring-loaded syringe, valve system, spike and tubing set for repetitive fluid dispensing. 10cc capacity.

CAN BE USED IN CONJUNCTION WITH SYRINGE FILTERS.

Syringe Fluid Dispensing System
10mL, Sterile, Disposable

Item No.
Price Each:

CT-FDS-10
Aud 55.00

NE-300 "Just Infusion"™ Syringe Pump

(Vial Filling)

Syringe Filter / Oils
Syringe Heater available
for viscous Fluids



NE-9000 Programmable Peristaltic Pump

- ✓ Pump from a reservoir through replaceable continuous tubing.
- ✓ Pump never comes into contact with fluid.
- ✓ Continuous Pumping. Up to 900 ml/min.
- ✓ Measured dispenses for production filling of bottles.
- ✓ "Learn and Repeat"™ for quick setup of production dispenses.
- ✓ Computer and foot switch interface.

Vial Filling
semi-Automatic

