

Distilled Water



Voltage V/Hz	Power W	Distilled water L/h	Chamber size Ømm	Outside Size mm	Weight Kg	Packing Size mm
220/50	750	1.5	180x200	290x290x390	3.5	245x285x510

Autoclave must adopt distilled water as operation medium. Even Purified water has various minerals, and can produce furring scale after evaporation at high temperature, which, after a period of time, can lead to leaking of steam due to improper closure. It can also lead to the clog in the narrow pipes and pressure sensor, to the malfunction of temperature sensor and to other failures. Once the furring scale enters handpiece and other instrument which are hollow inside or has apertures, it will clog the narrow pipes and axis, lowering the turning speed of handpiece, thus reducing its life span. Therefore, using distilled water is necessary.

HROMalytic +61(0)3 9762 2034
ECHnology Pty Ltd
Australian Distributors; Importers & Manufacturers

