


# Chromtech End Caps - Fittings, Tubing

**Stock Items (Melbourne)**

**Other sizes to ORDER  
(100-pk Min)**

- 1/16", 3mm (1/8"tube), 6mm(1/4"tube)
- Ideal for storage of chromatography columns



**Price  
(Ex. GST)**

These Chromtech CT PVC series caps are extremely flexible and durable. They can be stretched up to 15% larger than manufactured diameter. Offering a high quality glossy finish that makes them ideal to use on finished products, thread protectors, dust and moisture seals, pipe and stud end colour coding and for masking during painting and plating operations. The specially formulated PVC is unaffected by weather, oxidation, oils and moisture. PTs and PTFs can withstand temperatures up to +135C. They can be glued to other PVC items. Range of six basic colours. Black, Red, White, Yellow, Blue and Clear.

Part Number	A mm	B mm	Part Number	A mm	B mm
CT002-10	2.4	10	CT003-25	3	25
CT002-12	2.4	10	CT038-12	3	12
CT002-15	2.4	15	CT004-08	4.8	8
CT002-20	2.4	20	CT004-10	4.8	10
CT002-22	2.4	22	CT004-12	4.8	12
CT003-03	3	3	CT004-15	4.8	15
CT003-08	3	8	CT004-19	4.8	19
CT003-10	3	10	CT004-22	4.8	22
CT003-12	3	12	CT004-25	4.8	25
CT003-15	3	15	CT004-38	4.8	38
CT003-16	3	16	CT00406-05	7 gauge	5
CT003-20	3	20	CT00406-08	7 gauge	8
			CT00406-12	7 gauge	12