

GC DETECTORS

VICI microTCD

Unit (Cell & Electronics); 220V
TCD2-NiFe, VICI Microvolume TCD Detector



VICI Pulsed Discharge Detectors

- Non-radioactive,
Multiple Mode; Electron Capture;
HePhotoionisation
ASK for options pricing



Plasma Emission Detector,

Unique to LDETEK
High Purity Gas Analysis
- Many options
inc Custom GC Designs



DET - Thermionic Detectors

- Transducers Replacement Beads
NPD, PTID, TID-1, TID-3, TID-5, FTID, FID,
Remote FID, Tandem TID)
TID-1, NPD electronics, Ceramic TID beads
> Selectivity to Oxygenates and Branch HCs
see main site for details
DET Beads \$495ea
TID-1, NPD(TID-2 & TID-4)
specify GC Type; mounting flange/cable are different
Constant Current Power Supply



SRI GCS: full range FID, TCD, ECD, TID, FPD, HID, NPD, DELCD,