



BIOCHROM/HARVARD BIOSCIENCE/GE/PHARMACIA			
LIBRA S21, S22	L-990	P/XE	
ULTROSPEC 1000 E, I	LH-975	D2/W	INQUIRE
BIO-RAD			
BIOLOGIC DUAL FLOW			
760-1300	LH-25	D2	
760-1300	LHA-23	W	Lamp Ass'y
750-0216		HG	
1706, 1806A	L-295	D2	
1790, 1801, BIO-DIMENSION	LH-40	D2	
BIOFOCUS	LH-41	D2	
MICROSAMPLER 100	LH-41	D2	
REMEDEI	LH-115	D2	
SMARTSPEC 3000	L-60	P/XE	
BIOTAGE/DYAX/ARGONAUT			
HORIZON	L-735	HG	
ISOLERA SYSTEMS			
KILOPREP(100G)	LH-40	D2	
SP1, SP4, SPX	LH-60	D2	
FLASHMASTER II	LH-22	D2	
ISO-1SF, 1EF, 4SF, 4EF	L-735	HG	
ISO-1SV, 1SW, 1EV, 1EW	LH-60	D2	
ISO-4SV, 4SW, 4EV, 4EW	LH-60	D2	
BIOTEK INSTRUMENTS, INC			
EL-311	INQUIRE	W	
ELX-800 UV	L-1035	W	
ELX-808	INQUIRE	W	
FLX-800, 800T, 800TBI, 800TBID		W	
FLX-800TBX, 800 TC, 808/IU		W	
SYNERGY 2/AB, 2/FI/FP		W	
SYNERGY H4, HI/FI		W	
BISCHOFF			
LAMBDA 1000, 1010			INQUIRE
SP6, VISI-6			INQUIRE
BUCHI			
28062, 28063, 28067, 28068	LH-22	D2	
C-640			INQUIRE
BUCK			
BLC-20	LH-94	D2	
CAMAG			
TLC II SCANNER	L-285	D2	
CARLO ERBA/THERMO SCIENTIFIC			
	L-280	D2	
CARY/VARIAN			
1, 3, 4, 5, 100, 300, 400, 500	L-230	D2	
1, 3, 100, 300	L-930	W	
500	L-1015	W	
400	L-110	W	
50	L-1045	P/XE	
219	L-205	D2	
219	L-925	W	

CECIL			
CE 1200, 1220	LH-20	D2	
CE 1020, 2020, 2040, 2041	L-115	D2	
CE 4000 thru 9000 SERIES	L-115	D2	
CE 2501, 2502, 7250	L-115	D2	
COLE -PALMER			
01953-00	LH-40	D2	
DIONEX/LC PACKING/GYNKOTEK			
AD20, PDA20	LH-5	D2	
AD25, PDA100, PDA2996	LH-3	D2	
PDA20		W	INQUIRE
PDA100		W	
CES I			INQUIRE
DSA-1	LH-40	D2	
FDM-2		P/XE	INQUIRE
RF1002/2000	L-800	XE	
UDM-170S, 340S			INQUIRE
UDM-37182			INQUIRE
UVD 160, 170S, 170U, 320, 340S, 340U	LH-91	D2	
VDM-2	L-235	D2	
DAD 3000, MWD 3000	LH-45	D2	
VWD 064377, MWD 064654	LH-45	D2	
VWD 3000	LH-45	D2	
D-STAR			
DDW-10, DLC-20, DVW-10	LH-94	D2	
ECOM			
LCD 2073A	LH-60	D2	
LCD 2071.2, LCD 2071.3	L-735	HG	
SAPPHIRE, OPAL, TOPAZ	LH-60	D2	
FLASH 06	LH-60	D2	
DYAX/SEE BIOTAGE			
EM SCIENCE/SEE HITACHI			
ERMA/ERC			
7511, 7512, 7515, 7522, 7530, 7535	L-610	W	
ESA			
520, 522	LH-40	D2	
FISHER /SEE THERMO SCIENTIFIC			
GBC			
AA	L-200	D2	
CINTRA 906, 908, 911	L-200	D2	
CINTRA 916		W	INQUIRE
CINTRA 1155, 1255	LH-910	P/XE	
CINTRA 1200	L-280	D2	
CINTRA 1205, 1206, 1210	LH-40	D2	
CINTRA 5100	LH-82	D2	

<b>GE HEALTHCARE/AMERSHAM/BIOSCIENCE/PHARMACIA/LKB</b>			
2138, 2140, 2141			INQUIRE
2151	LH-20	D2	
4050, 4051, 4053, 4054	L-210	D2	
4050, 4060	L-835	W	
<b>AKTA SERIES</b>			
BASIC, EXPLORER, PILOT, PURIFIER	LH-990	P/XE	
AKTA MICROSYSTEM	LH-990	P/XE	
AKTA PRINME PLUS, PROCESS	LH-990	P/XE	
FPLC, AKTA PURIFIER 10, 100	L-1010 ***	HG	
ETAN LC 18-5050-50	LH-990	P/XE	
GENEQUANT II	LH-975	D2/W	INQUIRE
GENQUANT PRO			INQUIRE
RNA/DNA CALCULATOR			INQUIRE
NOVA SPEC II	L-835	W	
<b>ULTROSPEC SERIES</b>			
ULTROSPEC II, III, PLUS	L-210	D2	
ULTROSPEC II, III, PLUS	L-835	W	
ULTROSPEC 500 PRO	LH-125	D2/W	
ULTROSPEC 500 PRO	L-935	W	
ULTROSPEC 1000 e, 1000 i	LH-975	D2/W	INQUIRE
ULTROSPEC 2000, 3000, 4000	LH-125	D2/W	
ULTROSPEC 2000, 3000, 4000	L-935	W	
ULTROSPEC 3000 PRO	LH-125	D2/W	
ULTROSPEC 3000 PRO	L-935	W	
ULTROSPEC 1100 PRO	LH-975	D2/W	INQUIRE
ULTROSPEC 1100 PRO	L-935	W	
ULTROSPEC 2100 PRO, 3100 PRO	L-640	P/XE	
ULTROSPEC 3300, 4300	LH-125	D2/W	
ULTROSPEC 3300, 4300	L-935	W	
ULTROSPEC 5300 PRO, 6300 PRO	LH-125	D2/W	
UV-1	L-1030	HG	INQUIRE
UPC 900	L-1010 ***	HG	
UV 900	LH-990	P/XE	\$
***MUST RETURN SPENT LAMP AND HOUSING			
<b>GENTECH/SEE BEIJING RALEIGH</b>			
<b>GILFORD</b>			
200, 220, 222, 230, 240, 250, 252, 260	L-205	D2	
200, 220, 222, 230, 240, 250, 252, 260	L-875	W	
2000, 2400, 2600	L-205	D2	
2000, 2400, 2600	L-875	W	
FLUORO IV	L-801	XE	
RESPONSE I	L-210	D2	
RESPONSE I	L-900	W	
STASAR I, II, III		W	
SYVA S-3		W	

<b>GILSON</b>		
115, 116, 117	LH-150	D2
118, 119, 151, 152, 155, 156	LH-151	D2
122	L-800	XE
<b>GRACE/ALLTECH</b>		
200, 201, 205	LH-40	D2
200, 201, 205	LH-760	W
305	LH-910	P/XE
REVELERIS UV1, UV2	LH-40	D2
<b>GROTON</b>		
FD100	L-890	P/XE
<b>GYNKOTEK/SEE DIONEX</b>		
<b>HACH</b>		
DR/4000U	LH-110	D2
DR/4000U, 4000V	L-730	W
DR/5000, DR/6000	LH-78	D2
DR/5000, DR/6000	L-780	W
48000, 48100	L-730	W
<b>HEWLETT-PACKARD/SEE AGILENT</b>		
<b>HITACHI</b>		
100-10, 20, 30, 40, 60, 80	L-280	D2
101, 102	L-430	D2
101, 102, 111, 124	L-805	W INQUIRE
124	L-210	D2
181, 191	L-805	W INQUIRE
150-20	LH-51	D2
150-20	L-805	W INQUIRE
200, 300, 330, 340	LH-51	D2
655A-22, 23	LH-51	D2
F-1000 THRU F-7000	L-800	XE
GENESPEC III	LH-58	D2
GENESPEC III		W INQUIRE
L-2400, 2420, 2450, 2455	LH-56	D2
L-2420, 2450, 2455	L-835	W
L-2480, 2485	L-800	XE
L-2500	LH-51	D2
L-3000, 3100, 3200, 4200, 4250	LH-51	D2
L-4000, 4500	LH-56	D2
L-4200, 4250	L-815	W
L-7000, 7300, 7400, 7420, 7450A, 7455	LH-56	D2
L-7420	L-815	W
L-7480, 7485, 7490	L-800	XE
SERIES 100	L-805	W INQUIRE
U-0080D	LH-56	D2
U-1100, 1500	LH-51	D2
U-1100	L-815	W

HITACHI (CONTINUED)		
U-1800, 1900, 2800, 2810	LH-57	D2
U-1800, 2800, 2810	L-935	W
U-2900, 2910	LH-57	D2
U-1900, 2900, 2910	L-815	W
U-2000, 2001, 2010	LH-51	D2
U-2000, 2001, 2010	L-815	W
U-3000, 3010, 3110, 3210, 3300	LH-51	D2
U-3000, 3010, 3110, 3300	L-815	W
U-3310, 3400, 3410, 3501, 4100	LH-51	D2
U-3310, 4100	L-815	W
U-3410, 3501	LH-55	W
U-3500	LH-56	D2
U-4001, 4100	LH-51	D2
HORIBA JOBIN YVON/SEE JY		
IBM		
9420	L-200	D2
9523, 9563, 9565	L-205	D2
ISCO/TELEDYNE		
205	LH-40	D2
205	LH-760	W
3850	L-205	D2
COMBIFLASH COMPANION, TS, 4X, XL	LH-61	D2
COMBIFLASH RF 75, RF 200	LH-62	D2
COMBIFLASH SYSTEM SQ16X	L-205	D2
V4/1840	L-205	D2
JASCO		
870, 875	L-295	D2
FP-750, 777, 820, 821, 920, 921	L-800	XE
FP-1520, 1570, 1575, 2020, 3120	L-800	XE
FP-6200, 6300, 6500, 6600	L-800	XE
J-810 CD, J-815 CD	L-960	XE
CD-995, 2095, OR-2090	L-790	XE-HG
UVIDEC 100 THRU UVIDEC 100-VI	L-200	D2
UVIDEC 320, 340, 420, 430, 505	L-200	D2
UVIDEC 510, 610, 650, 710	L-200	D2
V-530, V-560	L-265	D2
V-530, V-550, V-560	L-835	W
UV-970A, UV975A	L-295	D2
V-650, V-660, V-7100	L-265	D2
JENWAY		
6105, 6405, 6505	LH-650	D2
6105, 6300, 6400	L-835	W
6405, 6500, 6505	L-835	W

JY(JOBIN YVON/SPEX/HORIBA)			
FLUOROLOG-2/F111XI, 112XI, F1T11I	L-855	XE	
FLUOROLOG-2/F212I, F222I	L-855	XE	
FLUOROLOG-3, TAU-3	L-855	XE	
FLUOROMAX-2, 3, 4	L-955	XE	
SYSTEM 5000	L-60	P/XE	
KNAUER/SONNTEK			
AZURA SERIES			
UVD 2.1S, 2.1L	A5193	D2	INQUIRE
BLUE SHADOW SERIES			
40D, 50D	A5193	D2	INQUIRE
SMARTLINE SERIES			
200	A4142	HG	INQUIRE
2500, 2600	A4071	D2	INQUIRE
2500	A4073	W	INQUIRE
2520, 2550	A5193	D2	INQUIRE
2550	A5195	W	INQUIRE
2800	A4447	D2	INQUIRE
2800	A4448	W	INQUIRE
WELLCHROM SERIES			
200	A4142	HG	INQUIRE
2000, 2001, 2500, 2501, 2600, 2650	A4071	D2	INQUIRE
2001, 2500, 2501, 2600	A4073	W	INQUIRE
6700, 8700, 9700, D871	A0877	D2	INQUIRE
KONTRON/BIOTEK/NORTH STAR SCIENTIFIC			
440	LHA-34	W	
720, 810, 820, 860, 900, 930, 940	L-205	D2	
820, 860, 900, 930, 940	L-110	W	
UVIKON	L-205	D2	
UVIKON XL, XS			INQUIRE
ALL SPECTROPHOTOMETERS	L-205	D2	
KRATOS/SEE ABI, PERKIN-ELMER			
LABALLIANCE/SCIENTIFIC SYSTEMS (SSI)			
200,201, 500, 504, 506, V520, 525	LH-40	D2	
200,201, 500, 504, 506, V520, 525	LH-760	W	
1200 SERIES			
AGG1314B	LH-14	D2	
AGG1315D	LH-13	D2	
AGG1321A	L-870	P/XE	
LC PACKINGS/DIONEX			
ULT-UVD, UVD-3000	LH-25	D2	
ULT-UVD	LHA-23	W	
LDC/SEE THERMO SCIENTIFIC			
LKB/SEE GE HEALTHCARE			
MCPHERSON			
FL-750/A & B	LH-101	XE	
FL-750/A & B	LH-100	D2	
FL-750/A & B	LH-102	XE-HG	

MERCK/SEE HITACHI		
<b>MICHROM BIORESOURCES/BRUKER</b>		
UV-200	LH-40	D2
MD1-MD4	LH-25	D2
<b>MICROTECH SCIENTIFIC, INC</b>		
LINEAR MT200/200CE	LH-40	D2
LINEAR MT201, 205	LH-40	D2
<b>MILTON ROY/SEE THERMO SCIENTIFIC/BAUSCH &amp; LOMB</b>		
<b>MOLECULAR DEVICES</b>		
SPECTRA MAX 340	L-65	P/XE
SPECTRA MAX 340 PC, 190, PLUS, 384	L-60	P/XE
SPECTRAMAX M2/M2E, M3, M4, M5/M5E	L-775	P/XE
SPECTRAMAX ANALYST GT	L-760	P/XE
SPECTRAMAX GEMINI, XPS, XS, EM	L-760	P/XE
SPECTRAMAX 250	L-760	P/XE
FLEX STATION, I, II, II 384	L-760	P/XE
<b>NICOLET/SEE THERMO SCIENTIFIC/ATI</b>		
<b>OLIS</b>		
RSM-1000	L-1000	XE
<b>PERKIN-ELMER/EG&amp;G</b>		
<b>HPLC UV DETECTORS</b>		
785A	LH-98	D2
INTEGRAL 4000	L-210	D2
LC-55, 55B, 65, 75, 85B, 95, 135	L-210	D2
LC-200	LH-98	D2
LC-200		W INQUIRE
LC-235	L-210	D2
LC-480	LH-92	D2
SERIES 200/DAD, PHASE I & 2	LH-97	D2
SERIES 200/DAD	L-725	W
SERIES 200/FLUORO	L-800	XE
SERIES 200/UV-VIS/785A	LH-98	D2
<b>FLUOROMETERS</b>		
SERIES 200, 203, 204, 512, 650	L-800	XE
MPF-2A, 3, 4, 43, 43A, 44, 44A, 44B	L-800	XE
<b>SPECTROPHOTOMETERS</b>		
46, 55, 57	L-200	D2
220, 320, 330, 340	LH-51	D2
550, 551, 552, 555	L-240	D2
552A, 553, 559A	L-200	D2
570, 571, 572, 575, 576	L-200	D2
3840	L-200	D2
BIO	LH-92	D2
BIO	LH-845	W
EZ150, 201, 210	LH-51	D2
EZ150, 201, 210	L-840	W

PERKIN-ELMER (CONTINUED)		
LAMBDA 1, 3, 3B	L-280	D2
LAMBDA 1, 3B	LH-845	W
LAMBDA 2, 5, 7, 9, 10, 11, 12	LH-92	D2
LAMBDA 2, 5, 7, 9, 10, 11, 12	LH-845	W
LAMBDA 4, 6	LH-93	D2
LAMBDA 4, 6	LH-845	W
LAMBDA 14 THRU 19	LH-92	D2
LAMBDA 14 THRU 19	LH-845	W
LAMBDA 20, 25, 35, 40, 45	LH-92	D2
LAMBDA 20, 25, 35, 40, 45	LH-845	W
LAMBDA 650,750,800,850,900,950,1050	LH-92	D2
LAMBDA 650,750,800,850,900,950,1050	LH-845	W
<b>ATOMIC ABSORPTION SPECTROMETERS</b>		
100, 300 ANALYSIS	L-335	D2
200, 400 ANALYSIS		D2
272, 280, 360, 370, 372, 373, 380	L-310	D2
380, 460, 560	L-310	D2
410, 420, 430		D2
600, 700, 800 ANALYSIS	L-405	D2
1100, 1100B, 2100, 4100	L-405	D2
2280, 2380, 3030B	L-310	D2
3100, 3110, 3300 ANALYSIS	L-335	D2
4000 ANALYSIS	L-310	D2
5000, 5100	L-310	D2
<b>PERSEPTIVE BIOSYSTEMS/SEE ABI</b>		
<b>PHARMACIA/SEE GE HEALTHCARE</b>		
<b>PHOTON TECHNOLOGY INT'L (PTI)</b>		
QM-2/2005, QM-4/2005, QM-6/2005	L-970	XE
QM-7/2005, QM-8/2005, QM-14/2005	L-970	XE
<b>POSTNOVA</b>		
S3210, 3240	LH-13	D2
<b>PYE UNICAM/SEE THERMO SCIENTIFIC</b>		
<b>RAININ/SEE VARIAN</b>		
<b>ROCHE/APPLIED SCIENCE</b>		
<b>LIGHT CYCLE</b>	<b>L-1040</b>	<b>XE/R</b>
<b>ROSEMONT ANALYTICAL</b>		
890	L-330	D2
<b>SCHOEFFEL/ SEE ABI, PERKIN-ELMER</b>		
<b>SCIENTIFIC COMPUTING INTERNATIONAL</b>		
<b>FILMTEK</b>		
4000, 4000-EM-DUV		
<b>SCINCO</b>		
S-3100, S-4100		
<b>SEDEX</b>		
55, 65, 75, 75C		
L-935		
<b>SEQUOIA TURNER</b>		
340, 360, 690		
W		

SHIMADZU			
HPLC UV DETECTORS			
LC2010, LC2010AHT, LC2010CHT	LH-84	D2	
SP-4	L-305	D2	
SPD-2A, 2AV, 6A, 6AV, 7A	L-305	D2	
SPD-6AV		W	INQUIRE
SPD-10A, 10AV, 10AVP, 10AVVP	LH-82	D2	
SPD-10AV, 10AVP, 10AVVP	L-830	W	
SPD-20A, 20AV	LH-82	D2	
SPD-20AV, SPD-M20AV	L-830	W	
SPD-M6A, M10A, M10AV		D2	INQUIRE
SPD-M6A, M10A,	L-935	W	
SPD-M10AV, M10AVP	L-245	W	
SPD-M10AVP	LH-88	D2	
SPD-M20A, M20AV	LH-87	D2	
SPD-M20A	L-245	W	
FLUOROMETERS			
RF-10A, 10AXL	L-800	XE	
RF-20A, 20AXS	L-795	XE	
RF- 535, 540, 551	L-800	XE	
RF-1501, 5000, 5301	L-800	XE	
SPECTROPHOTOMETERS			
BIOSPEC 1601	L-265	D2	
BIOSPEC 1601	L-110	W	
CS-920, 930,,	L-265	D2	
MULTISPEC 1501	L-265	D2	
MULTISPEC 1501	L-935	W	
UV-120, 120-01, 120-02, 140-02	L-265	D2	
UV-120	L-110	W	
UV-120-01, 120-02,		W	INQUIRE
UV-150, 150-02, 160A, 160U, 163	L-265	D2	
UV-160A, 160U	L-110	W	
UV-240, 250, 260, 265,	L-265	D2	
UV-240, 260, 265	L-110	W	
UV-1200, 1201, 1201PC, 1203	L-265	D2	
UV-1200, 1201, 1203	L-245	W	
UV-1240, 1501, 1600, 1601, 1601PC	L-265	D2	
UV-1240, 1240V, 1700, 1800	L-245	W	
UV-1501, 1600, 1601, 1601PC	L-110	W	
UV-1603, 1650, 1700, 1800, 2100	L-265	D2	
UV-1603, 1650	L-110	W	
2101, 2401PC, 2450, 2600, 2700	L-265	D2	
UV-2101PC, 2401, 2450, 2600, 2700	L-110	W	
UV-2501PC, 2550, 3100, 3100PC	L-265	D2	
UV-2501PC, 2550,	L-110	W	
UV-3101PC, 3600, 3700	L-265	D2	
UV-3101PC, 3600	L-110	W	
ATOMIC ABSORPTION SPECTROMETERS			
AA-625, 630, 640, 660, 670, 680U	L-265	D2	
AA-6200	L-265	D2	

SHODEX			
RI 73, 74, 101		W	
SIEVERS/IONICS INSTR. BUS.GROUP/GE ANALYT .INST.			
400, 800	L-700	HG	
900, 5310C, PPT	L-605	HG	
SKALAR			
6750	LH-22		
SPECTRA-PHYSICS/SEE THERMO SCIENTIFIC			
SPECTRONICS/SEE THERMO SCIENTIFIC			
SPEX/SEE JY (JOBIN YVON)			
TECAN			
SPECTRO FLUOR PLUS	L-1025	P/XE	
GENIOS	L-1025	P/XE	
INFINITE F200, M200 NANQUANT	L-1025	P/XE	
INFINITE F200 PRO, M200 PRO	L-1025	P/XE	
THERMO SCIENTIFIC			
AMINCO			
AB-2 (SERIES 2)			INQUIRE
SLM 2000	L-280	D2	
SLM 3000	L-270	D2	
SLM 4800, 8000, 8100	L-855	XE	
ATI			
ATI/240, 250	LHA-34	Lamp Ass'y	
	LH-34	Cartridge Only	
BAUSCH & LOMB			
AQUAMATE, AQUAMATE UV	L-910	W	
AQUAMATE UV	LH-54	D2	
BIOMATE 3		P/XE	INQUIRE
BIOMATE 5	L-910	W	
DE 50	L-225	D2	
EPSILON	L-915	W	
GENESYS 2, 5		D2	INQUIRE
GENESYS 2, 5	L-915	W	
GENESYS 6, 60, 10 BIO, 10 UV		P/XE	INQUIRE
GENESYS 8	LH-54	D2	
GENESYS 8	L-910	W	
GENESYS 10, 20	L-915	W	
EVOLUTION 60		P/XE	INQUIRE
EVOLUTION 60 S			INQUIRE
SPECTRONIC 20			INQUIRE
SPECTRONIC 21, -DV, MV, UVD, UVM	L-945	W	
SPECTRONIC 21D, 21 DUV	L-225	D2	
SPECTRONIC 200, 601, 700, 701, 710	L-225	D2	
SPECTRONIC 70, 100, 700, 701, 710	L-920	W	
SPECTRONIC 301, 401	L-915	W	
SPECTRONIC 501, 601 ROUND TOP	L-615	W	
SPECTRONIC 501, 601 POINTED TOP	L-915	W	
SPECTRONIC 1001, 1201, 2000	L-225	D2	
SPECTRONIC 1001, 1201	L-810	W	
SPECTRONIC 2000	L-920	W	
ENVIRONMENTAL INSTRUMENTS/SO2 ANALYZER			
43A, 43CT, 43CH	L-40	P/XE	
FISHER			
20-0344	LH-40	D2	



THERMO SCIENTIFIC (CONTINUED)			
<i>LDC</i>			
SC100-234	LH-40	D2	
SC100-234	LH-760	W	
SM I, II, III, D	L-205	D2	
SM 3000, 3100, 4000	L-205	D2	
SM 3100X, 4100	LH-42	D2	
SM 5000	LH-47	D2	
UV100-411, UV150-511	LH-40	D2	
UV100-411, UV150-511	LH-760	W	
<i>LINEAR</i>			
100, 200, 203, 204, 205, 206PHD	LH-40	D2	
100, 200, 201, 203, 204, 205	LH-760	W	
304, 305	LH-910	P/XE	
<i>NICOLET</i>			
EVOLUTION 100, 500, 600	LH-54	D2	
EVOLUTION 100	L-910	W	
EVOLUTION 500, 600	L-940	W	
<i>SPECTRA-PHYSICS/FINNIGAN</i>			
CE500, 1000, 2000	LH-40	D2	
CE500, 1000, 2000	LH-760	W	
FL2000, 3000	LH-910	P/XE	
SC100, 200, SF100, 200	LH-40	D2	
SC100, 200, SF100, 200	LH-760	W	
SF770, 970	LHA-1	D2 Lamp Ass'y	
	LH-1	Cartridge Only	
SP8400, 8440, 8440XR, SP8450, 8490	L-215	D2	
SP8480XR	LHA-3	D2 Lamp D2	
	LH-1	Cartridge Only	
SP8773XR	LH-1	Cartridge Only	
UV100, 150, 1000, 2000, 3000	LH-40	D2	
UV100, 150, 1000, 2000, 3000	LH-760	W	
UV3000HR	LH-40	D2	
<i>SURVEYOR</i>			
UV6000LP, 6057-60030, 60040	LH-46	D2	
UV6000LP, 6057-60030, 60040	L-765	W	
PDA1P, 5P, UV1P, 5P	LH-46	D2	
PDA1P, 5P, UV1P, 5P	L-765	W	
FL PLUS DETECTOR	L-800	XE	
<i>TURNER/BARNSTEAD</i>			
110, 111, 112	L-185	HG	
SP-830 PLUS, 850, 870 PLUS	W		INQUIRE
<i>UNICAM</i>			
UV 1-100	LH-54	D2	
UV 1-100	L-910	W	
UV 2-100, 200, 300, 400	LH-54	D2	
UV 2-100, 200, 300, 400	L-940	W	
UV 3-100, 200, UV 4-100, 200	LH-54	D2	
UV 3-100, 200, UV 4-100, 200, 500	L-940	W	
UV 300, 500, 550	LH-54	D2	
UV 300, 500, 550	L-940	W	
4225	LH-40	D2	
HELIOS ALPHA, BETA, DELTA	LH-54	D2	
HELIOS ALPHA BETA, DELTA	L-910	W	

THERMO SCIENTIFIC (CONTINUED)			
HELIOS EPSILON	L-915	W	
HELIOS GAMMA	LH-54	D2	
HELIOS GAMMA	L-910	W	
<i>TRANSGENOMIC</i>			
FDX-99-2100	L-800	XE	
WAVE 3500, 3500A, 3500HT	LH-56	D2	
WAVE 4500, 4500 HT, HT-HS	LH-56	D2	
WAVE HS, 4500, HT-HS	L-800	XE	
WAVE OLIGO, DNS	LH-56	D2	
<b>TURNER/BARNSTEAD /UNICAM SEE THERMO SCIENTIFIC</b>			
<b>VARIAN/CARY/RAININ/AGILENT</b>			
1, 3, 5, 14, 15, 16, 17	L-230	D2	
1, 3, 14, 15, 16, 17	L-930	W	
50	L-1045	P/XE	
87/8700, 8700	LH-20	D2	
100, 210, 219, 220, 300, 400, 500	L-230	D2	
100, 210, 219, 220, 300	L-930	W	
400	L-110	W	
500	L-1015	W	
2050	L-290	D2	INQUIRE
2200, 2300, 2400,	L-205	D2	
2550	L-295	D2	
9070, 9075	LH-910	P/XE	
AA200, 240, 275, 280,300, 400	L-230	D2	
AA600, 640, 775, 800, 875	L-230	D2	
AA975, 1275, 1475, 3000	L-230	D2	
DMS 70 THRU DMS 300	L-230	D2	
FL-1, 2	LH-910	P/XE	
LC-5000, 5500	L-305	D2	
PDA-1		D2	INQUIRE
PDA-2	LH-82	D2	
PRO-STAR, PRO-STAR 320	LH-950	P/XE	
PRO STAR 310	L-305	D2	
PRO STAR 325, 335	LH-325	D2	
PRO STAR 330	LH-56	D2	
PRO STAR 340, 345	LH-40	D2	
PRO STAR 340, 345	LH-760	W	
PRO STAR 360	LH-910	P/XE	
PRO STAR 363	L800	XE	
SPECTRA A10, A20, A30, A40	L-230	D2	
SPECTRA A640, A880	L-230	D2	
STAR 9050	L-305	D2	
UV-1, UV-1A	LH-950	P/XE	
UV 50, 100, 200	L-305	D2	
UV- C, DII, M	LH-40	D2	
VISTA 5000, 5500	L-305	D2	
<b>WATERS</b>			
420	L-770	HG	
440, 441	L-820	HG	
450	LHA-1	D2 Lamp Ass'y	
	LH-1	Cartridge Only	
470, 474	L-800	XE	
2420, 2424 ELSD	L-120	W	

WATERS (CONTINUED)		
480, 481	L-205	D2
484	LH-30	D2
486	LH-31	D2
490E	P/XE	INQUIRE
990, 991, 994	LHA-34	D2 Lamp Ass'y
	LH-34	D2 Replacement
996, 2996	LH-33	D2
2475	LH-925	XE
2487, 2488	LH-36X2	D2
LC MODULE 1/34551 SYSTEM	LH-31	D2



**HROMalytic** +61(0)3 9762 2034  
**ECHnology** Pty Ltd

**Australian Distributors**  
 Importers & Manufacturers  
[www.chromtech.net.au](http://www.chromtech.net.au)

Website NEW : [www.chromalytic.net.au](http://www.chromalytic.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) Tel: 03 9762 2034 . . . in AUSTRALIA

**Ask Jim about our Deal of the Month !**

LEGEND			
<b>CD</b>	CADIUM	<b>NA</b>	SODIUM
<b>D2</b>	DEUTERIUM	<b>P/XE</b>	PULSED XENON
<b>HG</b>	MERCURY	<b>W</b>	TUNGSTEN
<b>L</b>	LAMP ONLY	<b>XE</b>	XENON
<b>LH</b>	LAMP & HOLDER	<b>X2</b>	SUFFIX DESIGNATES LONG LIFE
<b>LHA</b>	LAMP HOLDER ASSEMBLY	<b>ZN</b>	ZINC

Lamp types and prices subject to change without prior notice. Chromtech assumes no liability for representations made herein. All information is accurate to the best of our knowledge. Majority of light sources listed are genuine OEM parts made to meet or exceed original equipment manufacturer's specifications. Chromtech is not affiliated with any companies mentioned; company or product names may be trademarks of their respective owners. 01/01/2014.

"Inquire" . . . May be in short supply or even discontinued ASK  
 ALL other Prices are available; but still need to check for availability  
 NIL Stock Melbourne. Allow 2 weeks import time ex-USA OEM  
 Prices are ex GST but include import freight  
 © All Rights Reserved